


Aug 12, 2022

URA3 PCR

Brian Teague¹¹University of Wisconsin - Stout1 *Works for me* Share


This protocol is published without a DOI.

Yeast ORFans CURE

Brian Teague
University of Wisconsin - Stout

ABSTRACT

The polymerase chain reaction (PCR) amplifies linear DNA using a DNA polymerase enzyme and a pair of short single-stranded DNA "primers." This protocol is for amplifying the URA3 gene from the YTK74 plasmid and adding the homology sequences needed for homology-directed repair.



PROTOCOL CITATION

Brian Teague 2022. URA3 PCR. **protocols.io**
<https://protocols.io/view/ura3-pcr-ce6tthen>

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 11, 2022

LAST MODIFIED

Aug 12, 2022

PROTOCOL INTEGER ID

68531

MATERIALS TEXT

- YTK74 PCR template, 1 ng/ul
- YFG-URA3 F forward primer, [M]10 micromolar (μM) concentration
- YFG-URA3 R reverse primer, [M]10 micromolar (μM) concentration
-

[Q5 Hot Start High-Fidelity 2X Master Mix - 100 rxns](#) **New England**

Biolabs Catalog #M0494S Step 4

- [Nuclease free water](#) **Contributed by users** Step 4
- A [200 μL](#) PCR tube

SAFETY WARNINGS

None of the materials used in this lab are hazardous.

HOWEVER, we are shedding nucleases -- enzymes that degrade DNA -- all the time. Wear lab coats and gloves to keep your samples nuclease-free.

- 1 Check that the thermocycler is programmed and holding at 98°C. The thermocycler program we're using is the following:

A	B
98°C for 30 seconds	
Repeat 35 times:	
	98°C for 5 seconds
	60°C for 15 seconds
	72°C for 30 seconds
72°C for 2 minutes	
Hold at 8°C	

Do not skip this step -- you don't want to wait for a thermocycler to warm up!

- 2 Grab an ice bucket and fill it with ice. If you don't have an ice bucket, a beaker will work in a pinch.

- 3 If necessary, dilute the primers to a concentration of **10 micromolar (μM)** in **TE Buffer Contributed by users** . (Make 100 μl of each dilution.)

Remember, by convention the blue-capped tubes from IDT have a concentration of **100 micromolar (μM)** .

- 4 Mix the following in a PCR tube **on ice, in this order**:

- **7 μL** **Nuclease free water Contributed by users**
- **1 μL** forward primer
- **1 μL** reverse primer
- **1 μL** template DNA
- **10 μL**

Q5 Hot Start High-Fidelity 2X Master Mix - 100 rxns New England Biolabs Catalog #M0494S

- 5 Mix the reaction by gently flick the tube several times, then spin down in a microfuge.

Do this quickly and return the tube to the ice bucket ASAP. There aren't enough thermocycler blocks for every group -- you may need to wait to share a thermocycler with other groups.

- 6 Transfer the tube from ice to a pre-heated thermocycler holding at 98°C. Start the PCR program.

- 7 After the PCR program has run, store the tube at -20°C.