



JAN 24, 2024

Freely moving imaging

In 1 collection

Mai-Anh Vu¹, mwhowe¹

¹Department of Psychological and Brain Sciences, Boston University, Boston, MA, United States

Team Cragg



Mai-Anh Vu

OPEN ACCESS



ABSTRACT

We have developed a new micro-fiber array approach capable of chronically measuring and optogenetically manipulating local dynamics across over 100 targeted locations simultaneously in head-fixed and freely moving mice, enabling investigation of cell-type and neurotransmitter-specific signals over arbitrary 3-D volumes . This protocol includes the steps for imaging in freely-moving mice. Please contact us (mwhowe@bu.edu) if you are interested in using this technique.

DOI:

dx.doi.org/10.17504/protocols.io.261ged1bjv47/v1

Protocol Citation: Mai-Anh Vu, mwhowe 2024. Freely moving imaging. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.261ged1bjv47/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: Jan 09, 2024

Last Modified: Jan 24, 2024

PROTOCOL integer ID: 93206

Keywords: ASAPCRN

Funders Acknowledgement:

Aligning Science Across
Parkinson's (ASAP)
Grant ID: ASAP-020370

- 1 For visualizing fibers in behaving mice custom-built miniscopes (Liberti et al., 2016; Liberti et al., 2017) were positioned above the fiber bundle, focused with an XYZ manipulator, and cemented in place with Metabond (Parkell) under brief isoflurane anesthesia.
 - 1.1 Note: Permanent implantation was used for initial experiments in this study, but baseplate attachment is equally viable.
- 2 Mice were allowed to recover for 24-48 hrs after attachment, after which imaging was performed in a custom built acrylic arena as mice foraged for Froot Loops (Kellogg's).
- 3 The miniscope GRIN objective had an N.A. of 0.45 for high-efficiency light collection and the imaging area was approximately 800 x 600 μm , allowing visualization of a substantial fraction of the fiber bundle.
- 4 Both behavioral video and fiber signals were acquired at 30Hz using custom MATLAB acquisition software.
- 5 Neural imaging data were acquired from DAQ and CMOS PCBs from the UCLA Miniscope project (Pnevmatikakis and Giovannucci, 2017).

- 6 Illumination intensity was delivered by a 470 nm LED (Luxeon Rebel, Lumileds) at an intensity of 0.1-0.25 mW/mm², controlled by a D-A interface (National Instruments USB-6009) and BuckPuck LED controller (LEDdynamics, Inc).
- 7 Imaging data were motion-corrected with NormCorre73, corrected for minor background light leak using spatial high-pass filtering, and converted to $\Delta F/F$.
- 8 Intensity timeseries for each fiber were extracted using manually selected circular ROIs.
- 9 Locomotor activity was quantified using manual tracking in ImageJ.