



DEC 06, 2022

WORKS FOR ME

1

Fluorescently labeled polyamine uptake (via Flow Cytometry)

DOI

dx.doi.org/10.17504/protocols.io.n92ldp8qxl5b/v1Marine Houdou¹, Peter Vangheluwe¹¹KU Leuven

Marine Houdou

COMMENTS 0

ABSTRACT

Assess polyamine uptake capacities of a specific cell line after incubation with fluorescently labeled polyamines and mean fluorescence intensity acquisition *via* flow cytometry.

DOI

dx.doi.org/10.17504/protocols.io.n92ldp8qxl5b/v1

PROTOCOL CITATION

Marine Houdou, Peter Vangheluwe 2022. Fluorescently labeled polyamine uptake (via Flow Cytometry). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n92ldp8qxl5b/v1>



KEYWORDS

Fluorescently labeled polyamines

LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Dec 02, 2022

LAST MODIFIED

Dec 06, 2022

PROTOCOL INTEGER ID

73483

GUIDELINES

Make sure to stay protect from the light when manipulating fluorescent materials.

MATERIALS TEXT

0.25% Trypsin-EDTA: Gibco, 25200056

Albumin Fraction V: Carl Roth, 8076.4

Dulbecco's Phosphate Buffered Saline modified without calcium chloride and magnesium chloride (DPBS (-/-)): Gibco, D8537

TrypLE Express Enzyme: Gibco, 12604021


Versene Solution: Gibco, 15040

BEFORE STARTING

- Prepare appropriate dilution of fluorescently labeled polyamine in cell culture medium.
- Prepare Flow Cytometry (FC) buffer made of 1% Albumin Fraction V in DPBS (-/-).
- Prepare Eppendorf tubes and FC tubes by labeling them and keeping them on ice.

- 1 Seed cells in 12 well-plate that they reach 70%-80% confluency the day of the experiment.

Note

- 2 The day of the experiment, remove cell culture medium and add fluorescently labeled polyamines to the cells in a final volume of  500 µL / well .

Note

Note



3 Incubate the desired time at 37°C, 5% CO₂, protected from the light.


4 Harvest cells and prepare samples for Flow Cytometry (FC).

4.1 Discard medium.


4.2 Wash with  500 µL /well Versene or DPBS (-/-).

4.3 Discard Versene or DPBS (-/-).

4.4 Add  200 µL /well 0.25% Trypsin or TrypLE. Incubate  00:05:00 at RT protect from the light.

4.5 Add  500 µL /well 1%-Albumin Fraction V in DPBS (-/-), collect the cells and transfer in Eppendorf tubes on ice.

4.6 Centrifuge for  00:05:00 at 400 g and 4°C.

- 4.7 Discard supernatants and resuspend cell pellets in  250-500 µL 1%-Albumin Fraction V in DPBS (-/-), depending on the size of the cell pellet.
- 4.8 Filter cell suspension through Nylon filter into FC tubes.
- 5 Keep on ice until acquisition at the flow cytometer. Record 10,000 live events per sample.