

FEB 20, 2023



DOI:

dx.doi.org/10.17504/protocol s.io.5qpvorrndv4o/v1

Protocol Citation: David A Priestman, Reuben Bush, Maria Leondaraki, Danielle te Vruchte, Kerri-Lee Wallom, María E Fernández-Suárez, Frances M Platt 2023. Analysis of glycosphingolipids from human cerebrospinal fluid. protocols.io https://dx.doi.org/10.17504/p rotocols.io.5qpvorrndv4o/v1

MANUSCRIPT CITATION:

Analysis of fluorescently labeled glycosphingolipid-derived oligosaccharides following ceramide glycanase digestion and anthranilic acid labeling DCA Neville, V Coquard, DA Priestman, DJM te Vruchte, DJ Sillence, Raymond A Dwek, Frances M Platt, Terry D Butters Analytical biochemistry 331 (2), 275-282, (2004). doi: 10.1016/j.ab.2004.03.051.

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Analysis of glycosphingolipids from human cerebrospinal fluid

Forked from Analysis of glycosphingolipids from human plasma

David A Priestman^{1,2}, Reuben Bush¹, Maria Leondaraki^{1,3}, Danielle te Vruchte¹, Kerri-Lee Wallom¹, María E Fernández-Suárez¹, Frances M Platt^{1,2,4,5}

¹Department of Pharmacology, University of Oxford, Mansfield Road, Oxford, OX1 3QT, United Kingdom;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815;

³Faculty of Medicine, Goethe University Frankfurt am Main, Theodor-Stern-Kai, 60590 Frankfurt am Main, Germany;

⁴Royal Society Wolfson Research Merit Award Holder;

⁵Wellcome Trust Investigator in Science.



Reuben Bush

ABSTRACT

Interest in the role of cellular glycosphingolipids (GSLs) in health and disease led to us developing a sensitive method to analyse the full complement of GSL structures present in mammalian cells, fluids and tissues. The original qualitative method we developed was published in 2004 and measured the oligosaccharides selectively released from glycosphingolipids using a ceramide glycanase enzyme derived from the medicinal leech. We have now updated and refined this protocol with the focus on achieving sensitive and reproducible quantitation of GSLs in control and patient cerebrospinal fluid samples. The method uses the fluorescent compound anthranilic acid (2-AA) to label oligosaccharides prior to analysis using normal-phase high-performance liquid chromatography. The labelling procedure is rapid, selective, and easy to perform. With the inclusion of a 2AA-labelled chitotriose calibration standard, it is possible to obtain accurate and reproducible molar quantities of individual GSL species.

GUIDELINES

This protocol requires the use of some hazardous materials. As such, users must be appropriately trained and hazardous materials stored, used, and disposed of in accordance with your institution's health and safety policies, and approved laboratory policies, risk assessments and codes of practice.

MATERIALS

Protocol status: Working

We use this protocol and it's

Created: Oct 30, 2022

working

Gradient grade for liquid chromatography LiChrosolv® Reag. Ph Eur

https://www.sigmaaldrich.com/catalog/product/mm/100030?

lang=en®ion=GB&gclid=EAIaIQobChMI1JyZu_Xa8AIVcwUGAB10lgVLEAAYASAA

EgJHmfD_BwE

Last Modified: Feb 20, 2023

Anthranilic acid

PROTOCOL integer ID:

72040

A89855 Sigma-Aldrich Reagent grade, ≥98%

Acetonitrile 1.00030

https://www.sigmaaldrich.com/catalog/product/sial/a89855?lang=en®ion=GB

Keywords:

Glycosphingolipids, Glucosyl Ceramide, HPLC, Oligosaccharide Analysis, **ASAPCRN**

Boric acid

ReagentPlus[®], ≥99.5% B0252 Sigma-Aldrich

https://www.sigmaaldrich.com/GB/en/product/sigald/b0252?cm_sp=Insite-_caSrpResults_srpRecs_srpModel_boric%20acid%20reagent%20plus-_-srpRecs3-1

Chloroform

Suitable for HPLC, ≥99.8%, amylene stabilized 34854-M Sigma-Aldrich

https://www.sigmaaldrich.com/catalog/product/sial/34854m?lang=en®ion=GB

Discovery® DPA-6S SPE Tube

https://www.sigmaaldrich.com/GB/en/search/dpa-6s? focus=products&page=1&perPage=30&sort=relevance&term=DPA-6S&type=product_name

Gibco™ DPBS, no calcium, no magnesium

14190094

https://www.thermofisher.com/order/catalog/product/14190094#/14190094

Glycosphingolipid standards

https://www.matreya.com/Departments/Product-Categories/Glycosphingolipids.aspx

Kinesis SPE Columns: TELOS® C18(EC) 100mg/1ml SPE Columns

https://kinesis.co.uk/kinesis-spe-columns-telosr-c18-ec-100mg-1ml-spe-columns-210-100m-001t.html

Ludger - BioQuant 2AA Labelled Chitotriose Standard

Cat. #: BO-CAA-CHI-01 Batch: B37I-02

https://www.ludger.com/docs/products/bg/bg-caa-chi/bg-caa-chi-01-b37i-02cofa.pdf

Ludger - 2-AA Labeled Glucose Homopolymer Ladder

Cat. # CAA-GHP-30

https://www.glycoprofiling.com/docs/products/caa/ludger-caa-ghp-30-guide.pdf

Methanol

34860 Sigma-Aldrich Suitable for HPLC, ≥99.9%

https://www.sigmaaldrich.com/catalog/product/sigald/34860?

lang=en®ion=GB&cm_sp=Insite-_-

caSrpResults_srpRecs_srpModel_methanol%2034860-_-srpRecs3-1

rEGCase I: recombinant Endoglycoceramidase custom-synthesized by GenScript https://www.genscript.com/

Sarstedt 1.5 ml Micro tubes with screw cap and seal

Product reference number 72.692

https://www.sarstedt.com/en/products/laboratory/screw-cap-micro-tubes-reaction-tubes/screw-cap-micro-

tubes/product/dmsarproducts/specificationPdf/Product/72.692.005/

Sodium acetate trihydrate

S1304 Sigma-Aldrich

Meets USP testing specifications

https://www.sigmaaldrich.com/catalog/product/sial/s1304?

<u>lang=en®ion=GB&gclid=CjwKCAjwtJ2FBhAuEiwAlKu19jp2DL5CW1Ke85Y0Ksb5b</u>oyV-5zKyUkzvg3uVVuZXV3abjGZbVruJhoCescQAvD_BwE

Sodium cyanoborohydride

Reagent grade, 95%

156159 Sigma-Aldrich

https://www.sigmaaldrich.com/catalog/product/aldrich/156159?

lang=en®ion=GB

TSKgel® Amide-80 HPLC Column


https://www.sigmaaldrich.com/GB/en/product/supelco/813071

SAFETY WARNINGS

This protocol requires the use of some hazardous solvents, reagents and chemicals. Refer to the Safety Data Sheets (SDS) provided by supplier and applicable Control of Substances Harmful to Health (COSHH). The correct personal protective equipment must be worn, and incidents reported in line with your institution's policy and procedures.

BEFORE START INSTRUCTIONS

Check that you have the required reagents, solvents, chemicals, equipment and PPE.



6 Separate into two phases: by adding 🔼 0.2 mL PBS and then 🚨 0.2 mL chloroform.

Vortex.

8 Centrifuge at 16,000 x g for 00:10:00 at room temperature.

10m



- 9 Remove very carefully the lower phase to a new tube and retain the upper phase.
- Dry down the lower phase under a stream of (oxygen-free) nitrogen in heating block (42 °C).
- 11 When dry, re-suspend the lower phase in \angle 20 μ L chloroform/methanol (1:3).
- 12 Add upper phase to lower phase and vortex.

- Load lower/upper phase mix onto column, let drip through gravity flow.
- Rinse sample tube with A 1 mL water, apply to column to wash.



16 Wash column with 3 x 🚨 1 mL water.



- 17 Elute GSLs into a new tube with:
- 17.2 <u>A 2 mL</u> chloroform/methanol (1:3).
- Vortex and leave Overnight at 4 °C or carry on to enzymatic digestion.



GSL digestion with EGC'ase I

16h

Vortex (\$\textit{L} 5 \text{ mL} \), C18) and dry down samples under a stream of nitrogen in heating block (\$\text{42 °C} \).

- 20 When about $\ \underline{\ \ }\ 150\ \mu L$ sample remaining, transfer to $\ \underline{\ \ }\ 1.5\ mL$ screw-cap tube.
- Rinse sample tube with $200 \, \mu L$ C:M 2:1, vortex and combine with the rest of the sample in the screw-cap tube.
- Rinse sample tube with Δ 200 μL chloroform, vortex and combine with the rest of the sample in the screw-cap tube.
- Dry down, under a slow stream of nitrogen in heating block (42 °C).
- 25 Re-suspend in 🚨 50 µL C:M 2:1, vortex, dry down under a very slow stream of nitrogen.

Note

NB It is essential that ALL the sample is dried in the bottom of the tube.

26 Add $\underline{\Delta}$ 90 μL enzyme/buffer to each sample and vortex:



- 26.1 rEGCase: $\Delta 5 \mu L$ enzyme (stored in freezer) plus $\Delta 85 \mu L$ buffer per sample.
- 26.2 rEGCase buffer (stored at 4 °C): 0.6 % Triton in pH 5.2.
- 27 Incubate at \$ 37 °C for \$ 16:00:00

16h



か

rEGCasel: recombinant Endoglycoceramidase I was custom synthesised by GenScript.

2AA labelling of glycans released from GSLs

1h 30m

28

Note

Note

Labelling mix: 30 mg/mL 2AA in labelling buffer and 45 mg/mL sodium cyanoborohydride **Labelling buffer**: 4% sodium acetate trihydrate and 2% boric acid in methanol

Dissolve 2AA in labelling buffer first with vortexing.

29 Add 2AA in buffer to sodium cyanoborohydride and vortex.



Safety information

Sodium cyanoborohydride is very toxic. Refer to the Safety Data Sheet.

30 Add Δ 310 μL labelling mix to the Δ 90 μL sample (digest) in Δ 1.5 mL screw-cap tube.





31 Incubate in oven at 8 80 °C for 01:00:00 , vortexing at 00:30:00

1h 30m

- Allow to cool to Room temperature
- 33 Add <u>A 1 mL</u> acetonitrile:water (97:3).



- Transfer from screw-cap tube to $\frac{\pi}{2}$ 15 mL tubes.
- Rinse screw-cap tube with 2 x $^{\text{L}}$ 1 mL acetonitrile:water (97:3) and add to $^{\text{L}}$ 15 mL tube.
- Pre-equilibrate 50 mg Discovery® DPA-6S SPE Tube (supplied by Sigma-Aldrich) with:
 - A 1 mL acetonitrile

 2 x A 1 mL water

 2 x A 1 mL acetonitrile
- 37 Apply samples ($\underline{\mathbb{Z}}$ 3.4 mL) to equilibrated DPA-6S columns, let drip through gravity flow.

Wash columns with acetonitrile:water (95:5) \rightarrow add \square 1 mL into the \square 15 mL tubes to wash and add to columns, and then add 3 x \square 1 mL acetonitrile:water (95:5) directly onto the columns.

- 39 Elute with \triangle 600 μ L water into new tubes.
- Water into new tubes.
- Dry down under a slow stream of nitrogen in a heating block (8 42 °C)
- Resuspend in \triangle 60 μ L of 70:30 acetonitrile:water
- 42 For HPLC, inject Δ 50 μL

HPLC protocol

1h

- The HPLC system consists of a Waters Alliance 2695 separations module and an in-line Waters 2475 multi λ -fluorescence detector set at Ex λ 360 nm and Em λ 425 nm.
- Purified 2AA-labeled oligosaccharides are separated and quantified by normal-phase highperformance liquid chromatography (NP-HPLC) as previously described (Neville *et al.*, 2004, see the reference below).

Note

Analysis of fluorescently labeled glycosphingolipid-derived oligosaccharides following ceramide glycanase digestion and anthranilic acid labeling

DCA Neville, V Coquard, DA Priestman, DJM te Vruchte, DJ Sillence, Raymond A Dwek, Frances M Platt, Terry D Butters

Analytical biochemistry 331 (2), 275-282, (2004).

doi: 10.1016/j.ab.2004.03.051.

- The solid phase used is a 4.6 × 250 mm TSK gel-Amide 80 column maintained at (Anachem, Luton, UK).
- The chromatographic flow rate is 0.8 mL/min, and run time was 60 min. The total run time is 0.1:00:00.

Time (min)	0	6	35	37	39	41	42	54	60
Flow (mL/min)	0.8	0.8	0.8	8.0	0.8	0.8	1.2	1.2	0.8
% Solvent A	71.6	71.6	52.8	23.0	23.0	71.6	71.6	71.6	71.6
% Solvent B	8.4	8.4	27.2	57.0	57.0	8.4	8.4	8.4	8.4
% Solvent C	20.0	20.0	20.0	20.0	20.0	20.0	20.0	20.0	20.0

Gradient conditions for Normal Phase HPLC. All chromatography was controlled and data were collected and processed using Waters Empower software.

- The mobile phases are acetonitrile (solvent A), de-ionised water (solvent B) and [M] 100 millimolar (mM) ammonium acetate, (p) 3.85 (solvent C).
- 48 Starting conditions are 71.6% solvent A. 8.4% solvent B and 20% solvent C at 0.8 mL/min constant

6m


flow-rate for 6 mins. 00:06:00

- 49 The gradient is developed from 6 min to 35 min by increasing solvent B from 8.4% to 27.2% with concomitant decrease in solvent A from 71.6% to 52.8 %. (5) 00:29:00
- 29m
- 50 From 35 to 37 min solvent B is increased to 57% and solvent A reduced to 23%, then maintained for 2 min up until 39 min. (5) 00:04:00
- 4m
- 51 Between 39 and 41 min the solvent ratios are returned to the starting conditions and then maintained for 1 min. § 00:03:00
- 3m
- 52 At 42 mins the flow-rate is increased to 1.2 mL/min and maintained for 12 min up to 54 min. **©** 00:12:00
- 12m

- 53 Between 54 and 60 min the flow-rate is gradually returned to 0.8 mL/min. (5) 00:06:00



- 6m
- 54 Individual GSL species are identified by their glucose unit values (GUs), calculated with the HPLC Empower software using a homopolymer dextran ladder (Fig. 2).



- 54.2 Prepare a mixture of commercially available authentic GSLs in a screw-cap tube:
- Matreya Neutral GSL mix Cat No 1505 (Glc-Cer, Lac-Cer, Gb3, Gb4) 1.0

Δ 100 μL Matreya Ganglioside mix Cat No 1510 (Lac-C	Cer, GM3, GD3) 0.5				
Δ 100 μL Matreya Ganglioside mix Cat No 1511 (GA1,	GM1a, GD1a, GD1b,GT1b) 0.5				
mg/ml A 30 µL Matreya GM2 Ganglioside Cat No 1502 (GM2)					
mg/ml A 30 µL Sigma Asialo-GM2 Cat No G9398 (GA2)	1.0 mg/ml				

- 54.3 Dry down the mixture under nitrogen and then digest, label and clean up as done for the CSF samples in steps 25-41.
- 54.4 Aliquot and store at 4 -20 °C
- In order to calculate molar quantities from integrated peaks in the chromatogram, inject a calibration standard containing 2.5 pmol 2AA-labelled chitotriose (Ludger) with each sample set (not shown).

Figures 1- 5: HPLC profile and sugar structures for GSLs in ...

55




Figure 1 HPLC profiles for GSL standards and 2AA-labelled homopolymer dextran ladder. Glucose Units were calculated using the retention times of a 2AA-labelled homopolymer dextran ladder and are shown for each of the GSL standard peaks.




Figure 2 GSL biosynthesis. Biosynthetic enzymes genes are indicated in the blue grid. Ganglioside names are abbreviated according to Svennerholm [1]. [1] L. Svennerholm, Designation and schematic structure of gangliosides and allied glycosphingolipids, Prog Brain Res 101 (1994) XI-XIV.




Figure 3 Typical HPLC profile for 2AA-labelled glycans released from ante-mortem human CSF glycosphingolipids with sugar structures indicated for each peak. In order to calculate values for total GSL concentrations, use the sum of the areas for all peaks from LacCer to GQ1b. GlcCer, which elutes before LacCer, is not included because it is not fully hydrolysed by EGC'ase I.




Figure 4 Typical HPLC profile for 2AA-labelled glycans released from post-mortem human CSF glycosphingolipids with sugar structures indicated for each peak. In order to calculate values for total GSL concentrations, use the sum of the areas for all peaks from LacCer to GQ1b. GlcCer, which elutes before LacCer, is not included because it is not fully hydrolysed by EGC'ase I.




Figure 5 Overlay of typical HPLC profiles for 2AA-labelled glycans released from human: plasma, post-mortem CSF, and ante-mortem CSF. Due to breakdown of the blood-cerebrospinal fluid barrier, post-mortem CSF is characterised by significant increases in GSLs that are prevalent in plasma, notably, GM3 and Gb3.

Note

Figures created with BioRender.com