



MAR 20, 2023

Immunohistochemistry for p53 staining in Breast Cancer Tissue

Freda Halim^{1,2}

¹Doctoral Degree, Faculty of Medicine, Universitas Padjajaran;

²Surgery Department, Faculty of Medicine, Pelita Harapan University



Freda Halim

ABSTRACT

These are protocols used for study of p53 expression differences in Luminal B Her-2 negative patients. The aim of the study is to show p53 expression differences in Luminal B Her-2 negative patients with and without primary endocrine therapy resistance. We used using paraffin sections of 67 samples and stained the slides for p53 expression, using primary antibody Clone DO-7 (DAKO).

OPEN ACCESS

DOI:

dx.doi.org/10.17504/protocols.io.5qpvard7dv4o/v1

Protocol Citation: Freda Halim 2023.

Immunohistochemistry for p53 staining in Breast Cancer Tissue. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.5qpvard7dv4o/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: Mar 13, 2023

Last Modified: Mar 20, 2023

PROTOCOL integer ID:

78640

Deparaffinization and Rehydration

1 Incubate slides in Xylenes for 3 minutes

3m

- 2 Incubate slides in Xylenes for 3 minutes 3m
- 3 Incubate slides in Xylenes for 3 minutes 3m
- 4 Rehydrate slides in 100% Ethanol for 3 minutes 3m
- 5 Rehydrate slides in 96% Ethanol for 3 minutes 3m
- 6 Rehydrate slides in 70% Ethanol for 3 minutes 3m
- 7 Rinse with running tap water and aquadest for 5 minutes 5m

Blockage of Endogenous Peroxidase

- 8 Incubate slides in 3% H₂O₂ for 15 minutes 15m
- 9 Rinse slides with running tap water and aquadest for 5 minutes 5m

Antigen Retrieval

40m

10 Antigen Retrieval

10.1 Antigen Retrieval with Tris EDTA (pH9) with pressure cooker, in 95⁰ Celcius temperature

20m

10.2 Open the lid and cool down in room temperature

15m

10.3 Rinse slides with running water and aquadest for 5 minutes

5m

10.4 Rinse in PBS (Phosphate Buffer Saline) in pH 7.40-7.60

5m

10.5 Excell Block

10m

10.6 Rinse in PBS in pH 7.40-7.60

5m

Primary Antibody

40m

- 11 Wipe excess liquid around the tissue
- 12 apply primary antibody (clone D07, Dako) 120μL
- 13 Incubate for 60 minutes
- 14 Rinse with PBS

1h

5m

Secondary Antibody

40m

- 15 Apply Excell Link as secondary antibody
- 16 Rinse with PBS
- 17 Apply Excell HRP as secondary antibody

15m

5m

20m

Signal Detection/ Histochemistry

27m

- | | | |
|-----------|---|------------|
| 18 | Apply DAB (Diamino-benzidine) 80-100µL for 10 minutes , Rinse with running tap water and aquadest for 5 minutes | 15m |
| | | |
| 19 | Apply Hematoxylline for 1 minutes, Rinse with running tap water and aquadest for 5 minutes | 6m |
| | | |
| 20 | Apply Tatcha's bluing solution and rinse with running tap water and aquadest for 5 minutes | 6m |

Dehydration and Clearing

20m

- | | | |
|-----------|---|------------|
| 21 | Clear excess water from the slides | |
| | | |
| 22 | Dehydrate Slides in 70%, 96%, 100% for 5 minutes each | 15m |
| | | |
| 23 | Incubate slides in Xylenes for 5 minutes | 5m |
| | | |
| 24 | Mount the slides | |