



MAR 28, 2024

Protein Digestion with S-trap Spin Columns using Conditioned Concentrated Media

Forked from [Protein Digestion with S-trap Spin Columns](#)

Joanna Bons¹, J P Rose¹, M A Watson¹, B Schilling¹

¹Buck Institute for Research on Aging

Schilling Lab



M A Watson
Buck Institute

ABSTRACT

Trypsin digestion of isolated proteins using S-trap Spin columns in preparation for downstream proteomic profiling.

For trypsin digestion of proteins from conditioned media Step 1 is modified as we add a volume of media rather than a concentration of protein homogenate. Step 6 is modified because we use S-trap micro spin columns as conditioned concentrated media samples contain $\leq 100 \mu\text{g}$ of protein. In step 14 and 16 we use $2 \mu\text{g}$ sequencing-grade trypsin instead of the usual 1:25 (wt/wt) trypsin/protein ratio.

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.x54v928eml3e/v1

Protocol Citation: Joanna Bons, J P Rose, M A Watson, B Schilling 2024. Protein Digestion with S-trap Spin Columns using Conditioned Concentrated Media. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.x54v928eml3e/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Mar 25, 2024

MATERIALS

Last Modified: Mar 28, 2024

PROTOCOL integer ID: 97346

Keywords: Trypsin Digestion, S-trap spin columns, Proteomics, Mass Spectrometry, Conditioned Media, SASP, Secretome

- S-trap spin columns (Protifi)
- Centrifuge
- Centrifugal vacuum concentrator
- Small hot box/incubator
- HPLC-grade water
- 10% SDS solution
- 1 M triethylammonium bicarbonate (TEAB) solution, pH 8
- 100 mM triethylammonium bicarbonate (TEAB) solution, pH 8 in water
- Dithiothreitol (DTT; 250 mM in 100 mM TEAB, pH 8)
- Iodoacetamide (IAA; 250 mM in 100 mM TEAB, pH 8)
- S-trap buffer (90% methanol in 100 mM TEAB)
- Sequencing-grade trypsin
- 12% phosphoric acid in water

- Digestion buffer/Elution buffer 1: 50 mM triethylammonium bicarbonate (TEAB) solution, pH 8 in water
- Elution buffer 2: 0.5% formic acid (FA) in water
- Elution buffer 3: 50% acetonitrile (ACN), 0.5% formic acid (FA) in water
- 0.2% formic acid (FA) in water

- 1** In a 2-mL microcentrifuge tube add an appropriate proportion of volume from the concentrated condition media, 10% SDS for a final concentration of 4% SDS, 1M TEAB pH 8 solution for a final concentration of 50 mM TEAB, and HPLC-grade water if necessary to bring the volume up to a minimum of 50 μ L.
- 2** Add 250 mM DTT for a final concentration of 20 mM and incubate for 10min at 50°C to reduce the proteins. Then immediately leave for 10 min on the bench at room temperature (RT).
- 3** Add 250 mM IAA for a final concentration of 40 mM and incubate for 30min at RT in the dark to alkylate the proteins.
- 4** Acidify the sample with 12% phosphoric acid for a final concentration of 1.2%.
- 5** Add 7 volumes of S-Trap buffer to the acidified lysate and mix immediately by inversion. Formation of protein colloid may be observed.

- 6 Use S-trap micro spin columns as conditioned concentrated media samples contain ≤ 100 μg of protein. Ensure that the S-Trap spin column is in a 2.0-mL flow-through catch tube.
- 7 Add 100 μL of the acidified lysate/S-Trap buffer mix into the S-Trap spin column. Centrifuge at 4,000 x g for 10 seconds or until all the solution has passed through the column. Discard the flow-through.
- 8 Repeat step 7 until the entire acidified lysate and S-Trap buffer mix has passed through the column.
- 9 Add 200 μL of S-Trap buffer to wash the column. Centrifuge at 4,000 x g for 10 seconds or until all solution has passed through the column. Discard the wash solution.
- 10 Add 200 μL of S-Trap buffer and set aside.
- 11 Prepare the solution of sequencing-grade trypsin.
- 12 Centrifuge the S-trap spin column at 4,000 x g for 10 seconds or until the column is dry.
- 13 Place the S-Trap spin column into a clean 2.0-mL elution tube.

- 14** Add 2 μg of the sequencing-grade trypsin solution to the column. Be careful not to pierce the trap material.
- 15** Loosely cap the S-trap micro spin column or close the S-trap mini spin column and incubate for 1 hour at 47°C with no agitation.
- 16** After 1 hour, add another 2 μg of the sequencing-grade trypsin solution to the column and incubate overnight at 37°C with no agitation.
- 17** After overnight incubation take the samples out of the incubator and elute sequentially in the same 2.0-mL elution tube as follows:
 - 17.1** a. Add 80 μL of elution buffer 1 (50 mM TEAB, pH 8) and centrifuge at 1,000 x g for 1 minute.
 - 17.2** b. Add 80 μL of elution buffer 2 (0.5% FA) and centrifuge at 1,000 x g for 1 minute.
 - 17.3** c. Add 80 μL of elution buffer 3 (50% ACN, 0.5% FA) and centrifuge at 4,000 x g for 1 minute.
- 18** Vacuum dry the eluted peptide solution in a centrifugal vacuum concentrator.

- 19 Reconstitute the dried peptides in 0.2% FA and thoroughly mix the solution before desalting.