



Sep 08, 2021

Collection and preservation of eDNA from marine water samples

Ana Ramón-Laca¹, Abigail Wells¹, Linda Park¹¹NWFSC-NOAA

1 Works for me



Share

dx.doi.org/10.17504/protocols.io.bx37pqrn

Ana Ramón-Laca

DISCLAIMER

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

ABSTRACT

This protocol is designed for water collection from Niskin bottles and filtration at sea using reusable filter cups.

Aim: to collect and filter 2.5 L of water at each depth from each CTD cast and preserved the filter at room temperature.

DOI

dx.doi.org/10.17504/protocols.io.bx37pqrn

PROTOCOL CITATION

Ana Ramón-Laca, Abigail Wells, Linda Park 2021. Collection and preservation of eDNA from marine water samples. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bx37pqrn>

KEYWORDS

eDNA, marine, water, fisheries

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 08, 2021

LAST MODIFIED

Sep 08, 2021

GUIDELINES

Water should be filtered immediately after collection.

Only DNA-free forceps should touch the filter!

MATERIALS TEXT

- DNA away
- Kimwipes
- gloves
- 69 oz Whirlpack bags (2.5 L is at 9.5" of the bag from the bottom)
- 2L pitchers (for whirlpack bags stability)
- utensils for getting membrane from filter (i.e. forceps)
- 500 ml filter cups with a rubber stopper fitted with an adapter for filter cup
- mixed cellulose ester sterile filters (1 µm, 47 mm diameter) (Advantec® Cat. A100H047A)
- pump
- 3-port manifold
- tubing adapters
- wastewater container carboy
- labeled tubes containing 2 ml of Longmire's buffer
- bleach (5%) and bleach bucket
- distilled water
- 2 buckets (for bleach and rinsing)
- mesh bags
- drying racks
- bungees
- absorbent towels
- distilled water from the evaporator onboard and from the lab

Lysis buffer recipe (Longmire et al 1997):

To make 1 liter add in numerical order:

1. 975 ml double-distilled water
2. 100 ml of 1 M Tris-HCL, pH 8.0
3. 200 ml of 0.5 M EDTA, pH 8.0
4. 2 ml of 5 M NaCl
5. 25 ml of 20% SDS (w/v)

Filter the buffer with an autofill PES bottle top filtration device (sterile 500ml, 0.22 µm)

DISCLAIMER:

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Filter cups cleaning and assemblage

- 1 Set up wash station (two buckets, one with 0.5% bleach, one with clean distilled H₂O)

- 2 Replace bleach and water stations at regular intervals (at least every other day, 0.5% bleach = 1.2L bleach in 20L bucket)
- 3 Place all used small items in a small mesh bag and the filter cups in a large mesh bag in bleach for at least 30 minutes
- 4 Rinse the abovementioned items for at least 10 minutes
- 5 Allow to dry on the drying racks
- 6 Wearing a fresh pair of gloves, mount filter cups with new sterile filter

Collection of water from Niskin bottles

- 7 Let the water from the hull ("surface water sample") run for > 3 min while the CTD is being deployed
- 8 Collect 2 samples of 2.5 L each of "surface water" for each CTD cast
- 9 Once the rosette has been brought back to the surface, water needs to be extracted from each Niskin bottle and transported into the lab for filtration. Each Niskin contains 10L of seawater. We want to collect 2.5L of water from each Niskin.

Wipe the spigot on the Niskin with DNAway to remove potential contaminants on the outside of the Niskin
- 10 Flush water (count to 10) from the Niskin to further reduce risk of contamination and collect 2.5 L of water in a whirlpack bag and place in 2L pitcher
- 11 Include a sampling negative control by collecting 2.5L of distilled H₂O (at least daily)

Filtration of water samples

- 12 Wipe working area with DNAway
- 13 Switch on the pump and set vacuum in the middle range between -8 and -12 bars (keep an eye on the gauge & adjust if necessary)

- 14 Pour water from whirlpack bag (supported by pitcher) into the filter cup, repeat until finished (2.5L of water)
- 15 Using DNA-free forceps fold filter and place in pre-labeled tube pre-filled with preservative buffer (2ml, Longmire Buffer in a 5 ml LoBind tube).
- 16 Filter 2L of ship distilled H₂O every time the cleaning distilled H₂O is replaced as a sampling negative control
- 17 Every 2-3 days 2L of lab diH₂O should be filtered as sampling negative controls
- 18 Note tube label in notes. Also note time filtered, membrane type, place filtered (field/lab, which lab), etc.
- 19 Store samples at room temperature away from UV light