



Sep 16, 2020

## Cell culture

In 1 collection

Yingchao Xue<sup>1,2</sup>, Xiping Zhan<sup>3</sup>, Shisheng Sun<sup>4</sup>, Senthilkumar S. Karuppagounder<sup>5,6,7</sup>, Shuli Xia<sup>2,5</sup>, Valina L Dawson<sup>5,6,7,8,9</sup>, Ted M Dawson<sup>5,6,7,8,10</sup>, John Laterra<sup>2,5,8,11</sup>, Jianmin Zhang<sup>1</sup>, Mingyao Ying<sup>2,5</sup>

<sup>1</sup>Department of Immunology, Research Center on Pediatric Development and Diseases, Institute of Basic Medical Sciences, Chinese Academy of Medical Sciences and School of Basic Medicine, Peking Union Medical College, State Key Laboratory of Medical Molecular Biology;

<sup>2</sup>Hugo W. Moser Research Institute at Kennedy Krieger; <sup>3</sup>Department of Physiology and Biophysics, Howard University;

<sup>4</sup>College of Life Sciences, Northwest University; <sup>5</sup>Department of Neurology, Johns Hopkins University School of Medicine;

<sup>6</sup>Neuroregeneration and Stem Cell Programs, Institute for Cell Engineering, Johns Hopkins University School of Medicine;

<sup>7</sup>Adrienne Helis Malvin Medical Research Foundation; <sup>8</sup>Department of Neuroscience, Johns Hopkins University School of Medicine;

<sup>9</sup>Department of Physiology, Johns Hopkins University School of Medicine;

<sup>10</sup>Department of Pharmacology and Molecular Sciences, Johns Hopkins University School of Medicine;

<sup>11</sup>Department of Oncology, Johns Hopkins University School of Medicine

1 Works for me This protocol is published without a DOI.

Neurodegeneration Method Development Community

Tech. support email: [ndcn-help@chanzuckerberg.com](mailto:ndcn-help@chanzuckerberg.com)



Anita Broellochs  
protocols.io

### ABSTRACT

This protocol explains the cell culture and characterization of lines ND1014, N1, and ND27760 from *Synthetic mRNAs Drive Highly Efficient iPS Cell Differentiation to Dopaminergic Neurons*.

### EXTERNAL LINK

<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC6344911/>

### THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Synthetic mRNAs Drive Highly Efficient iPS Cell Differentiation to Dopaminergic Neurons. Xue Y, Zhan X, Sun S, Karuppagounder SS, Xia S, Dawson VL, Dawson TM, Laterra J, Zhang J, Ying M. Stem Cells Transl Med. 2019 Feb;8(2):112-123. doi: 10.1002/sctm.18-0036. Epub 2018 Nov 1. PMID: 30387318

### PROTOCOL CITATION

Yingchao Xue, Xiping Zhan, Shisheng Sun, Senthilkumar S. Karuppagounder, Shuli Xia, Valina L Dawson, Ted M Dawson, John Laterra, Jianmin Zhang, Mingyao Ying 2020. Cell culture . **protocols.io**  
<https://protocols.io/view/cell-culture-9e6h3he>

### MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

Synthetic mRNAs Drive Highly Efficient iPS Cell Differentiation to Dopaminergic Neurons. Xue Y, Zhan X, Sun S, Karuppagounder SS, Xia S, Dawson VL, Dawson TM, Laterra J, Zhang J, Ying M. Stem Cells Transl Med. 2019 Feb;8(2):112-123. doi: 10.1002/sctm.18-0036. Epub 2018 Nov 1. PMID: 30387318

### EXTERNAL LINK

<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC6344911/>

### COLLECTIONS ⓘ



**Protocols for Synthetic mRNAs Drive Highly Efficient iPS Cell Differentiation to Dopaminergic Neurons**

### KEYWORDS

Culture, maintenance, ND1014, N1, ND27760, ipsc, SNCA

#### LICENSE

\_\_\_\_\_ This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

#### CREATED

Nov 17, 2019

#### LAST MODIFIED

Sep 16, 2020

#### OWNERSHIP HISTORY

Nov 17, 2019  Liz Brydon [Protocols.io](#)

Sep 16, 2020  Anita Broellochs [protocols.io](#)

#### PROTOCOL INTEGER ID

29886

#### PARENT PROTOCOLS

Part of collection

[Protocols for Synthetic mRNAs Drive Highly Efficient iPS Cell Differentiation to Dopaminergic Neurons](#)

#### MATERIALS

NAME	CATALOG #	VENDOR
<a href="#">ReLeSR™ 100 mL</a>	5872	<a href="#">Stemcell Technologies</a>
<a href="#">mTeSR™ 1 500 mL Kit</a>	5850	<a href="#">Stemcell Technologies</a>
<a href="#">CytoTune™ iPS 2.0 Sendai Reprogramming Kit</a>	A16517	<a href="#">Thermo Fisher</a>

#### SAFETY WARNINGS

Please refer to the Safety Data Sheets (SDS) for safety and environmental hazards.

#### BEFORE STARTING

Obtain approval to work with human stem cells from an appropriate Institutional Review Board.

#### Cell Culture

- 1 Obtain human iPSC lines from Coriell Cell Repositories or through reprogramming patient cells.  
Step 1 includes a Step case.

#### **Purchase iPSC Reprogram iPSC**

\_\_\_\_\_ step case \_\_\_\_\_

#### **Purchase iPSC**

Obtain human iPSC lines from Coriell Cell Repositories (derived from normal human skin fibroblasts).

- 2 Characterize pluripotency of iPSCs by immunochemistry for pluripotent cell markers (NANOG, OCT4, TRA-1-60, and SSEA-3) and through embryoid body formation assay.
- 3 Maintain iPSCs as feeder-free cultures in mTESR1 medium in 5% CO<sub>2</sub>/95% air condition at **37 °C** and passage using ReLeSR.
- 4 Perform karyotype analysis of G-banded metaphase chromosomes to confirm the chromosomal integrity of iPSCs.

