



JAN 31, 2024

🌐 Optogenetic modulation of dopaminergic neurons

Cristian González-Cabrera¹, Matthias Prigge¹

¹Neuromodulatory Network Group, Leibniz Institute for Neurobiology, Magdeburg.

ASAP Collaborative Research Network

TeamPrigge



priggelab

ABSTRACT

This protocol describes an optogenetic activation of neuromelanine-laden dopaminergic neurons in the SNc-VTA. We induce neuromelanin (NM) in the SNc-VTA through viral expression of a humane tyrosinase virus. Furhtermore we express a red-absorbing opsin in one hemispherend and blue absorbing in the other hemisphere. We detail how we optogenetic stimulateanimals in over four weeks.

This optogenetic deep brain stimulation induces a reduction in neuromelanine levels in dopaminergic neurons, and rescues behavioral phenotypes as evaluate in the grip strength test.

OPEN  ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.dm6gp3kppvzp/v1

Protocol Citation: Cristian González-Cabrera, Matthias Prigge 2024. Optogenetic modulation of dopaminergic neurons. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.dm6gp3kppvzp/v1>

MANUSCRIPT CITATION:

Optical deep brain stimulation can reactivate NM-laden neurons and rescue behavioral phenotypes.

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: Jan 30, 2024

Last Modified: Jan 31, 2024

PROTOCOL integer ID: 94413

Keywords: Optogenetics, optical deep brain stimulation, Neuromelanin

Funders Acknowledgement:

ASAP

Grant ID: 020505

Viral Injection

1 Stereotactic injection of viruses into DAT-Cre (JAX# 006660) animals

Viral Constructs:

1.- pAAV-hSyn1-dlox-ChrimsonR_tdTomato(rev)-dlox-WPRE

2.- pAAV-hSyn1-dlox-hCHR2(H134R)_mcherry(rev)-dlox-WPRE

3.- pAAV_Ef1a_DIO_hTyrHA_minBack

Injections:

- Bilateral injections were performed.
- Right hemisphere: viruses 1 and 3.
- Left hemisphere: viruses 2 and 3.
- Viruses were combined 50:50 and injected a total Volume of 600ul per hemisphere

Stereotactic coordinates:

- Substantia nigra

ML: 1.4mm AP: 3.25mm DV: 4.0 (blunt needle / NF34BL-2)

Injection set-up:

- Hamilton syringe coupled to a WPI injector (UMP3T-1) ([link](#))

- Injection speed 100nl/min

Experimental Timeline

- 2
 - Animals were left to accumulate neuromelanin for 6 or 10 weeks.
 - Animals were stimulated for 30 minutes for 4 weeks. Stimulation during 5 days and two days rest.
 - The stimulation was 635nm, 5ms pulse, 10 Hz at at 10mW. 3 seconds ON 10 seconds OFF.

Optical stimulation

- 3
 - fibers are attached bilaterally to the animals
 - the animal is placed in open field arena
 - stimulation is controlled via PulsePal (OpenEphys [link](#)) on a MRL_III-635L laser