



1 ▼

Nov 30, 2021

dsDNA quantification using Sybr Green I V.1

James JN Kitson¹¹Newcastle University

1



protocol .

Network Ecology Group



James Kitson
Newcastle University

This is a simple protocol that uses Sybr Green 1 and a microplate reader to quantify dsDNA concentrations in unknown samples. This is best suited for situations where a high number of samples need to be quantified (e.g. normalisation of samples before pooling for Illumina sequencing). For smaller sample numbers (e.g. < 50 samples) it will be more efficient to simply use a qubit or nanodrop.

James JN Kitson 2021. dsDNA quantification using Sybr Green I. **protocols.io**
<https://protocols.io/view/dsdna-quantification-using-sybr-green-i-b2f7qbrn>
 James Kitson



The volumes required for n samples in table 2 were unclear, I have updated the table headers to fix this.

Sybr Green 1, dsDNA, quantification

_____ protocol ,

Nov 30, 2021

Nov 30, 2021

55519

10,000X Sybr Green 1

TE buffer

Black flat bottomed microtitre plates

Plate reader

Part 1 - Make a dilution series of a Lambda DNA standard:

- 1 Using [NEB lambda DNA \(item N3011S\)](#) 10 ng/μl dsDNA standard in TE (20ul of 500ng/μl lambda DNA + 980 μl of TE).
 - *TE = 10 mM Tris HCl, 1 mM EDTA, pH 8*

- 2 Serially dilute this as below:

A	B	C	D
	lambda DNA (μL)	TE (μL)	lambda DNA (ng/μl)
1	1000 (from 10 ng/μl stock)	0	10
2	500 (from dilution 1)	500	5
3	500 (from dilution 2)	500	2.5
4	500 (from dilution 3)	500	1.25
5	500 (from dilution 4)	500	0.625
6	500 (from dilution 5)	500	0.313
7	500 (from dilution 6)	500	0.156
8	0 (from dilution 7)	500	0

Table 1: Dilution series for control DNA samples

- 3 Aliquot 25 μl of each of these into PCR strips and store at -20 °c until needed. Thaw one for each run on the plate reader.
 - *If you have previously made your dilution series start at this point.*

Part 2 - Make and distribute Sybr Green 1 working solution.

- 4 Make a 100x Sybr Green 1 solution (5ul [Sybr Green 1](#) + 495 μl of TE).
 - *This can be pipetted out into 50 - 100 μl aliquots and stored at -20 °c. Once thawed, do not reuse.*
 - *96 μl of 100x Sybr Green 1 is required per plate of samples.*
 - *If you have previously made your 100X Sybr Green 1 start at this point.*

- 5 Make a solution of 1x Syber Green 1 as below:

A	B	C
	Volume for one sample	Volume for n + 136 samples
TE	94	
100X Sybr Green 1	1	
Total	95	

Table 2: Illustrative table for calculating the dilution to make enough 1x Sybr Green 1 for n samples

Part 3 - dilute Sybr Green 1 working solution and add DNA.

- 6 Distribute 95 μl of Sybr Green 1 working solution into the first 4 columns of a Nunc black microtiter plate and 100 μl into the last 8 columns.

- *This is best done using a multichannel electronic repeater or 96 well pipette if you are testing a lot of samples.*

7 Using a p10 multichannel pipette, add 5ul of each standard into the first 4 columns of the control plate as shown in figure 1.




Figure 1: Control plate layout for microplate reader when reading DNA samples.

8 Distribute 98 μl of sybr green 1 working solution into n wells of a nunc black microtiter plate.

- *Use multiple plates for more than 96 samples.*

9 Add 2 μl of sample dsDNA into each of the sample wells.

10 Seal with a clear plate seal.

11 Incubate at room temperature in the dark for 10 minutes. ⌚ 00:10:00

10m

Part 4 - load microplate reader and run the assay.

- 12 Run samples on your plate reader according to manufacturers instructions with an excitation wavelength of 485 nM and an emission wavelength of 535 nM.
- 13 Check that your calibration curve has a very high R^2 value (greater than 0.95) and multiply your measured concentrations by 2.5 to get a final concentration.




Figure 2: A typical calibration curve from the above protocol