



Upload image


Jun 10, 2020

Gel Slice Sample Preparation for Proteomics

Jennifer Gin¹, Leanne Chan¹, Christopher Petzold¹¹Lawrence Berkeley National Laboratory

2 Works for me dx.doi.org/10.17504/protocols.io.bgj3juqn

LBNL-omics

 Jennifer Gin

ABSTRACT

This protocol details steps for gel slice sample preparation for proteomic data acquisition. It was adapted from Shevchenko et al. "In-gel digestion for mass spectrometric characterization of proteins and proteomes." *Nature protocols* 1.6 (2006): 2856.

EXTERNAL LINK

<https://doi.org/10.1038/nprot.2006.468>

DOI

[dx.doi.org/10.17504/protocols.io.bgj3juqn](https://doi.org/10.17504/protocols.io.bgj3juqn)

PROTOCOL CITATION

Jennifer Gin, Leanne Chan, Christopher Petzold 2020. Gel Slice Sample Preparation for Proteomics.
protocols.io
[dx.doi.org/10.17504/protocols.io.bgj3juqn](https://doi.org/10.17504/protocols.io.bgj3juqn)



EXTERNAL LINK

<https://doi.org/10.1038/nprot.2006.468>

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

May 18, 2020

LAST MODIFIED

Jun 10, 2020

PROTOCOL INTEGER ID

37211

MATERIALS

NAME	CATALOG #	VENDOR
Eppendorf tubes 1.5 mL uncolored	022363204	Eppendorf Centrifuge
DTT	DTT-RO	Millipore Sigma
Iodoacetamide	I1149	Millipore Sigma
Water LC-MS grade B&J Brand	BJLC365-2.5	VWR Scientific
Ammonium Bicarbonate LC-MS grade	BJ40867-50G	VWR Scientific
Acetonitrile LC-MS grade B&J Brand	BJLC015-2.5	VWR Scientific
Sequencing Grade Modified Trypsin Porcine	V511A	Promega

NAME	CATALOG #	VENDOR
Acetic Acid	818755	Millipore Sigma
Formic Acid	28905	Thermo Fisher Scientific
Eppendorf ThermoMixer C	2231000667	pipette.com
Eppendorf Vacufuge Concentrator	07-748-13	Fisher Scientific

SAFETY WARNINGS

Wear proper PPE (gloves, safety goggle, and lab coat), and prepare solvents in a chemical fume hood.

Store organic solvents in a flammable storage cabinet when not in use.

Discard used solvents and buffers in appropriate waste containers.

Buffer Preparation

30m



Buffers to prepare fresh:

100 ml of 100 mM Ammonium Bicarbonate (AMBIC):	0.79 g Ammonium Bicarbonate (VWR Scientific, Cat.#BJ40867-50G) 100 ml LC-MS grade Water (VWR Scientific, Cat.#BJLC365-2.5)
1 ml of 50mM Ammonium Bicarbonate/Acetonitrile (AMBIC/ACN):	250 µl 100 mM Ammonium Bicarbonate 250 µl LC-MS grade Water 500 µl Acetonitrile (VWR Scientific, Cat.#BJLC015-2.5)
1 ml of 10 mM 1,4-Dithiothreitol (DTT), 25 mM Ammonium Bicarbonate (AMBIC):	10 µl 1M 1,4-Dithiothreitol (Millipore Sigma, Cat.#DTT-RO) 250 µl 100 mM Ammonium Bicarbonate 740 µl LC-MS grade Water
1 ml of 55 mM Iodoacetamide, 25 mM Ammonium Bicarbonate (IAA/AMBIC):	10 mg Iodoacetamide (Millipore Sigma, Cat.#I1149) 250 µl 100 mM Ammonium Bicarbonate 750 µl LC-MS grade Water
0.200 ml of 100 ng/ul of Trypsin:	20 µg Sequencing Grade Modified Trypsin (Promega, Cat.#V511A) 200 µl 50 mM Acetic Acid Millipore Sigma, Cat.#818755)
0.100 ml of 12.5 ng/ml Trypsin/50 mM Ammonium Bicarbonate (Trypsin/AMBIC):	12.5 µl 100 ng/ul Trypsin 50 µl 100 mM Ammonium Bicarbonate 37.5 µl LC-MS grade Water
1 mL of 50% Acetonitrile/5% Formic Acid (ACN/FA):	0.5 ml Acetonitrile 0.05 ml Formic Acid (Thermo Fisher Scientific, Cat.#28905) 0.45 ml LC-MS grade Water
1mL of Buffer A	0.001 ml Formic Acid 1 ml Water



Store DTT in **-20 °C** and Trypsin in **-80 °C**.

IAA is light sensitive. Store in amber tube (Fisher Scientific, Cat.#05-402-31).

Day 1: Washing the gel slices

1h 15m


15m





Wash gel slice with **150 µl Water** in a tube (Eppendorf Centrifuge, Cat.#022363204). Incubate for **00:15:00** then remove liquid.

3  15m

Wash gel slice with  **150 µl 50 mM AMBIC/ACN** and incubate for  **00:15:00** then remove liquid.

 Discard used buffers and solvents in appropriate waste containers.




4  15m

Add  **50 µl Acetonitrile** and incubate for  **00:15:00** then remove liquid.

5  20m




Add  **100 µl 50 mM AMBIC/ACN** and incubate for  **00:20:00** then remove liquid.

6  10m




Add  **50 µl Acetonitrile** to cover gel slice. The gel slice will visibly shrink and turn white. Incubate for  **00:05:00** to  **00:10:00** then remove liquid.

Day 1: Trypsin digestion 3h

7  1h

Add  **50 µl 10 mM DTT, 25 mM AMBIC** and incubate in  **56 °C water bath** for  **01:00:00** then remove liquid.

8  45m

Add  **50 µl 55 mM IAA, 25 mM AMBIC** and incubate in the dark, at  **Room temperature** for  **00:45:00** then remove liquid.

9  15m




Wash gel slice by adding  **50 µl 50 mM AMBIC/ ACN** and incubating for  **00:15:00** then remove liquid.

10 

Repeat wash step 9.






5m

Add  **50 μ l Acetonitrile** to cover gel slice. The gel slice will visibly shrink and turn white. Incubate for  **00:05:00** to  **00:10:00** then remove liquid.



30m

Add  **10 μ l 12.5 ng/ml Trypsin, 50 mM AMBIC** and incubate on  **On ice** for  **00:30:00** , then remove excess liquid, if necessary.



All of the liquid (trypsin) should be absorbed into the gel slice at the end of 30 minutes. The gel slice will look clear/hydrated again. There is usually no liquid leftover because only a small amount of trypsin is added however, if necessary, remove excess trypsin before overnight digestion.



15m

Add  **50 Millimolar (mM) AMBIC** to cover the gel slice.



Wait 15-20 minutes to make sure the gel slice is still completely covered in liquid. Sometimes the gel is not completely rehydrated after the trypsin has been absorbed, and it absorbs a lot of the AMBIC if you only use enough to cover. This is especially a problem with larger gel pieces.



16h

Digest the protein(s) at  **37 °C**  **Overnight** in an incubator or thermocycler.




Day 2: Recovering digested peptides

4h

15 Remove gel slice from  **37 °C** incubator or thermocycler.



15m

Add  **20 μ l 50 mM AMBIC** after tryptic digestion. Incubate at  **Room temperature** for  **00:15:00** .



Collect liquid in a clean tube.

18 

Add  **30 µl 50% ACN/ 5% FA** to cover the gel slice.

19 

20m

Incubate for  **00:20:00** on a benchtop shaker (Pipette.com, Cat.#2231000667) set at  **Room temperature**.

20 

Collect liquid.

21   

Repeat steps 18-20.

22  

5m

Add  **20 µl Acetonitrile**. Incubate for  **00:05:00** then collect liquid.


23 

Use a SpeedVac (Fisher Scientific, Cat.#07-748-13) to dry the liquid completely.

Day 2: Sample Prep For MS

5m

24 

After the tubes are completely dry in Step 24, add  **17 µl Buffer A** then transfer to plastic autosampler vials (Agilent, Cat.#5182-0567, #5182-0564) or 96-well plate (Bio-Rad, Cat.#HSP9655).

25 

Store at  **-20 °C** until ready for LC-MS/MS analysis.