



May 23, 2022

## LRRK2 RCKW Protein Purification

David Snead<sup>1,2</sup>, Yu Xuan Lin<sup>1</sup>, Mariusz Matyszewski<sup>1</sup>

<sup>1</sup>Department of Cellular and Molecular Medicine, University of California, San Diego, La Jolla, CA 92 093;

<sup>2</sup>Department of Biochemistry and Molecular Biology, Johns Hopkins University Bloomberg School o f Public Heath, Baltimore, MD, 21205



dx.doi.org/10.17504/protocols.io.81wgb6693lpk/v1

Mariusz Matyszewski

Protein purification protocol for tag-less LRRK2<sup>RCKW</sup> as done by Leschziner and Reck-Peterson Labs. Same protocol can be used to purify LRRK1<sup>RCKW</sup> as well.

Original Protocol by David Snead. Modified by Yu Xuan Lin and Mariusz Matyszewski for publication.

DOI

dx.doi.org/10.17504/protocols.io.81wgb6693lpk/v1

David Snead, Yu Xuan Lin, Mariusz Matyszewski 2022. LRRK2 RCKW Protein Purification. **protocols.io** 

https://dx.doi.org/10.17504/protocols.io.81wgb6693lpk/v1

	_
ASAP Grant ID: ASAP-000519	
MJFF Grant ID: 18321	

LRRK2, protein purification

\_\_\_\_\_ protocol,

Feb 09, 2022

May 23, 2022



```
57957
```

```
Equipment:
Ultracentrifuge, we used one with a Ti70 rotor
Dounce homogenizer
FPLC with SP, His, and S200I columns.
Also make sure Ni-NTA beads and TEV protease are available.
Buffers (suggested volumes listed):
Day 1:
■300 mL Lysis Buffer:
[M]50 millimolar (mM) HEPES p-7.4
[M] 0.5 Molarity (M) NaCl
[M]20 millimolar (mM) Imidazole pH8.0
[M] 0.5 millimolar (mM) TCEP
[M]5 % volume Glycerol
[M]5 millimolar (mM) MgCl2
[M]20 micromolar (µM) GDP
Plus protease inhibitors (Protein inhibitor tablets and [M] 0.5 millimolar (mM) Pefabloc )
□300 mL Wash Buffer:
[M]50 millimolar (mM) HEPES p+7.4
[M] 0.5 Molarity (M) NaCl
[M] 0.5 millimolar (mM) TCEP
[M]5 % volume Glycerol
[M]5 millimolar (mM) MgCl2
[M]20 micromolar (µM) GDP
■100 mL Elution Buffer:
[M]50 millimolar (mM) HEPES p-7.4
[M] 0.5 Molarity (M) NaCl
[M]300 millimolar (mM) Imidazole p+8.0
[M] 0.5 millimolar (mM) TCEP
[M]5 % volume Glycerol
[M]5 millimolar (mM) MgCl2
[M]20 micromolar (µM) GDP
```

mprotocols.io

```
■300 mL Dilution Buffer :
[M]50 millimolar (mM) HEPES p-7.4
[M] 0.5 millimolar (mM) TCEP
[M]5 % volume Glycerol
[M]5 millimolar (mM) MgCl2
[M]20 micromolar (µM) GDP
■500 mL Buffer A (250 mM NaCl) :
[M]50 millimolar (mM) HEPES p+7.4
[M]250 millimolar (mM) NaCl
[M] 0.5 millimolar (mM) TCEP
[M] 5 % volume Glycerol
[M]5 millimolar (mM) MgCl2
[M]20 micromolar (µM) GDP
■300 mL Buffer B (2.5 M NaCl) :
[M]50 millimolar (mM) HEPES |p+7.4
[M]2.5 Molarity (M) NaCl
[M] 0.5 millimolar (mM) TCEP
[M]5 % volume Glycerol
[M]5 millimolar (mM) MgCl2
[M]20 micromolar (µM) GDP
Day 2:
□500 mL Buffer 0:
[M]50 millimolar (mM) HEPES p+7.4
[M] 0.5 Molarity (M) NaCl
[M] 0.5 millimolar (mM) TCEP
[M]5 % volume Glycerol
[M]5 millimolar (mM) MgCl2
[M]20 micromolar (µM) GDP
■300 mL Buffer 300 :
[M]50 millimolar (mM) HEPES p-7.4
```



```
[M]0.5 Molarity (M) NaCl

[M]300 millimolar (mM) Imidazole p⊦8.0

[M]0.5 millimolar (mM) TCEP

[M]5 % volume Glycerol

[M]5 millimolar (mM) MgCl2

[M]20 micromolar (µM) GDP

□300 mL Storage Buffer :

[M]20 millimolar (mM) HEPES p⊦7.4

[M]0.7 Molarity (M) NaCl

[M]0.5 millimolar (mM) TCEP

[M]5 % volume Glycerol

[M]2.5 millimolar (mM) MgCl2

[M]20 micromolar (µM) GDP
```

Day 1: Protein Pellet Lysis

1h 28m

- Ready the Ultracentrifuge and cool it down to 8 4 °C
- 2 Make Lysis Buffer (check Materials section)
- 3 Thaw protein pellets & On ice
- 4 Resuspend each pellet in **■40 mL Lysis Buffer**

Might need to use a spatula and a pipette to resuspend fully.

5 Homogenize each pellet with Dounce homogenizer & On ice

protocols.io

15 plunges loose, followed by 17 plunges tight. 5.1 Gentle when homogenizing, don't completely pull douncer out of liquid, avoid bubbles Pour lysate into rotor tubes (pre-chilled, clean, check for cracks; we used Ti70 tubes) 6.1 Make sure they are balanced. Use Lysis buffer to balance if needed. 1h 28m 7 Spin lysate in ultracentrifuge \$\mathbb{G}\$50000 rpm, 4°C, 01:28:00 , Ti70 rotor Prep the rest of Day 1 buffers while waiting 8 Day 1: Ni-NTA gravity column 2m Start equilibrating 12 mL of Ni-NTA beads with Lysis buffer when there is about 30 mins left on the centrifuge 9.1 We equilibrated by resuspending **a** 6 **m** L in a 15 mL falcon tube, twice. They were spun down at **31000 rpm, 4°C, 00:02:00, 3 times**, each time getting rid of liquid and doing a 50:50 resuspension. 10 After ultracentrifuge is done, pour the supernate into 3x 50 mL falcon tubes.

Add **4 mL equilibrated Ni-NTA beads** to 40 mL of supernate. Bring up to **50 mL with Lysis buffer** 

protocols.io

- 12 Incubate in cold room, while rotating "hot-dog style".

  © rpm, 4°C, 01:00:00 medium-slow rotation
- 13 After incubation, add onto a gravity column. Wash with remaining **Lysis Buffer**, followed by **200 mL Wash Buffer**
- 14 Elute with 50 mL Elution Buffer
  Make sure to resuspend the beads with stopper closed on the column. Might need multiple repeats of resuspension for best results. Leave column dry afterwards.
- Dilute the solution to [M]250 millimolar (mM) NaCl by adding □50 mL Dilution Buffer
- 16 Syringe filter the solution.

## Day 1: FPLC Cation Exchange and Overnight TEV cleavage

- 17 Use a 5 mL HiTrap SP FF column. Pump wash FPLC with **Buffer A** and **Buffer B**. Equilibrate the column with **Buffer A**.
- 18 Apply sample onto the column with a sample pump
- Run a gradient program, [M]250 millimolar (mM) NaCl to [M]2.5 Molarity (M) NaCl , 0% B to 100% B.
- 20 Combine protein fractions into a 50 mL Falcon Tube. Check protein concentration; should be around [M]1.1 micromolar (μM)
  - 20.1 Afterwards dilute to [M]600 millimolar (mM) NaCl

## 20.2 Bring up to **50 mL with Wash Buffer**

Add TEV, and incubate overnight at § 4 °C while rotating "hot-dog style".

Our final TEV concentration was around [M]0.2 micromolar (µM)

## Day 2: FPLC HisTrap Column

- 22 Make Day 2 buffers. Make sure to degas the Storage Buffer.
- Put the 5 mL HisTrap column on the FPLC. Pump wash the FPLC with **Buffer 0** and **Buffer 300**. Equilibrate the column with **Buffer 0**.
- 24 Run a Step Gradient Program

0% for 4 CV
6.7% for 5 CV
13.3% for 5 CV
26.7% for 5 CV
50% for 4 CV
100% for 6 CV (Column Wash)
0% for 5 CV (Column Re-equilibration)

25 Combine protein elution into a 15 mL tube before concentrating for Size Exclusion Column

Day 2: SEC Column 2m



26	Equilibrate a S200I 10/300 column with the <b>Storage Buffer</b> .
27	Concentrate the protein to <b>□1 mL</b> using 100 kD cutoff concentrators. Should be done nearing the end of the column equilibration.
	We concentrated using 2 concentrators to save time.
	27.1 After concentrating, filter the protein using a 0.1 uM spin filter. Spin at max speed for © 00:02:00 at § 4 °C
28	Inject filtered sample and run elution program.
29	Collect and concentrate the protein to about $\  \  \  \  \  \  \  \  \  \  \  \  \ $
	About 3-6 mL of protein initially.
30	Aliquot and freeze protein. Make sure to check pureness of protein as well.
	2, 5, and/or 10 uL aliquots recommended.