



FEB 01, 2024

STC-1 Passaging and PD Associated Microbe Stimulation Protocol

Malu G Tansey¹

¹College of Medicine |University of Florida



Senthilkumar Karuppagounder

ABSTRACT

STC-1 Passaging and PD Associated Microbe Stimulation Protocol

MATERIALS

Complete Medium
DMEM, high glucose, no glutamine
4mM L-glutamine
1mM sodium pyruvate
5% FBS

Other Reagents
0.05% Trypsin-EDTA
1X D-PBS without calcium and magnesium
75cm² Flask
12 well plates

OPEN ACCESS



DOI:
dx.doi.org/10.17504/protocols.io.261gedmn7v47/v1

Protocol Citation: Malu G Tansey 2024. STC-1 Passaging and PD Associated Microbe Stimulation Protocol.
protocols.io
<https://dx.doi.org/10.17504/protocols.io.261gedmn7v47/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Feb 01, 2024

Last Modified: Feb 01, 2024

PROTOCOL integer ID: 94578

Keywords: ASAPCRN

Funders Acknowledgement:

Aligning Science Across
Parkinson's (ASAP) Collaborative
Research Network

Passaging & Plating Protocol

- 1 Cells are passaged when they reach approximately 70% confluence

Warm trypsin, 1XPBS, and complete medium in water bath for 15 min.
- 2 Aspirate off culture medium
- 3 Rinse cells with 10mL D-PBS
- 4 Add 3mL Trypsin to flask. Place flask back into 37°C/5% CO₂ incubator for 7 minutes.
- 5 Observe cells under microscope to ensure they have detached from the bottom of the flask.

a. These cells are particularly sticky, so giving them a nice tap may be necessary to detach them.
- 6 Add 9mL of complete medium to flask and transfer cell suspension to a 15mL tube.

- 7 Spin down at 130g for 5 min
- 8 Aspirate off supernatant and resuspend cell pellet in 1mL culture medium.
- 9 Count cells
- 10 Prepare cell suspension at a concentration of 3×10^5 cells/mL
- 11 Add 1mL cell suspension per well in 12 well plate(s)
- 12 For PD associated microbe stimulation, allow cells to rest for approximately 24 hours,
- 13 Treat STC-1s with 10 colony forming units (CFU) of PD associated microbe per STC-1 cell

- 14** Incubate for desired time and then collect conditioned media and/or cells for downstream analyses

LPS Pre-Treatment

- 15** Following plating, allow cells to rest for 24 hours
- 16** Treat cells with 10ng/mL LPS for 24 hours
- 17** 24 hr after LPS treatment, stimulate cells with 10 CFU of PD associated microbes per STC-1 cell
- 18** Incubate for desired time and then collect conditioned media and/or cells for downstream analyses