

Aug 30, 2024

Preparing 10X Genomics CytAssist FFPE Sections

DOI

dx.doi.org/10.17504/protocols.io.n92ldme3ol5b/v1

YuHong Fu^{1,2}, Glenda Halliday^{1,2}

¹Brain and Mind Centre & Faculty of Medicine and Health School of Medical Sciences, The University of Sydney, Sydney, NSW 2050, Australia;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815



YuHong Fu

USTD

OPEN  ACCESS



DOI: **dx.doi.org/10.17504/protocols.io.n92ldme3ol5b/v1**

Protocol Citation: YuHong Fu, Glenda Halliday 2024. Preparing 10X Genomics CytAssist FFPE Sections. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n92ldme3ol5b/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 11, 2024

Last Modified: August 30, 2024

Protocol Integer ID: 93359

Keywords: FFPE, 10X, Sectioning, ASAPCRN



Funders Acknowledgement:

**The Michael J. Fox
Foundation for Parkinson's
Research (MJFF) and the
Aligning Science Across
Parkinson's (ASAP) initiative
Grant ID: ASAP-020529
The Michael J. Fox
Foundation for Parkinson's
Research (MJFF) and the
Aligning Science Across
Parkinson's (ASAP) initiative
Grant ID: ASAP-020505**

Abstract

Designing and properly mounting FFPE sections on large-window CytAssist slides is important for downstream analysis and data quality. This protocol describes the key steps for this preparation.



Set-up before cutting FFPE sections

- 1 Clean the surface of the bench area, microtome (HistoCore MULTICUT R, Leica), Water bath (HI1210, Leica), and tools (including forceps and blades) with 80% ethanol. 10m
- 2 Spray RNase ZapTM (Thermo Fisher Scientific) on the surface of the bench area, instruments, and tools and leave them to completely dry before use. 30m
- 3 Filled the Water bath with RNase/DNase free water and set the temperature at 40°C. 30m

Cut FFPE sections

1d 4h 15m

- 4 Cut FFPE sections at 6µm thickness.
Note: The time cost indicated for this step is for one FFPE block. 15m
- 5 Mount sections in the center region (< 11mm x 11mm) of the poly-prep slide according to the 10X genomics analysis requirement.
Note: The time cost for this step is experiment task-dependent. 4h
- 6 Incubated mounted sections at 42°C for 3 h and dried overnight at room temperature in a desiccator. 1d
- 7 Ready for downstream CytAssist procedure.