

FEB 06, 2023

OPEN ACCESS

Protocol Citation: Andreas Sagen 2023. General initiation protocol for HEK-Blue cells.

protocols.io

<https://protocols.io/view/general-initiation-protocol-for-hek-blue-cells-cna4vagw>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development
We are still developing and optimizing this protocol

Created: Jan 25, 2023

Last Modified: Feb 06, 2023


PROTOCOL integer ID:
75836

Keywords: HEK, Endotoxin, TLR, TOLL-Like receptor

General initiation protocol for HEK-Blue cells

Andreas Sagen¹

¹The National Institute of Occupational Health in Norway



Andreas Sagen

University of Oslo, The National Institute of Occupational H...

ABSTRACT

HEK-Blue is a product from Invivogen, which provide reporter cells for endotoxin-testing among others. Here is a generalized protocol for initial seeding of cells with reporter characteristics.

GUIDELINES

The protocol has been proven to work with different HEK-Blue cells, but not all. Take time to optimize.

MATERIALS

Centrifuge
Laminar flow cabinet
CO2 incubator

SAFETY WARNINGS



Use appropriate PPE when working with selective antibiotics.

BEFORE START INSTRUCTIONS

Prepare complete initiation medium

Complete initiation medium (CIM)

- Initially, we want to expand the cells quickly and make them healthy. This require double the amount of FBS and no selection antibiotics in the medium for the first 2-3 passages.

1.1 In  400 mL DMEM, mix  100 mL FBS

5m

Note

Always use ultra-low endotoxin (>0.1 EU/mL) mediums, FBS and supplements if possible

Materials:

 DMEM high glucose GlutaMAX **Gibco - Thermo Fischer Catalog #31966021**

 ULTRA-LOW ENDOTOXIN FETAL BOVINE SERUM (FBS) **BioWest Catalog #S1860**



1.2 Aliquot CIM in  50 mL units and store frozen at  -20 °C

Note


Create descriptive labels containing the following information: HEK-Blue initiation medium, DMEM w/GlutaMAX, 20% FBS, Date

Thawing from frozen stock

2 The cells are fragile in the beginning. Take extra good care of the cells when thawing.


2.1 Add  11 mL warm CIM in a T-75 flask and incubate for  00:30:00 in a 5% CO₂ incubator

30m


2.2 Get a vial of HEK-Blue cells, and thaw in a warm water bath for approx.  00:02:00

2m

2.3 Transfer vial content to a centrifuge tube with  5 mL warm CIM

2.4 Centrifuge with  200 rcf, 00:05:00

5m

2.5 Discard supernatant and resuspend pellet in  1 mL warm CIM

2.6 Add all content from centrifugation tube to T-75 and incubate in a 5% CO₂ incubator



Note

Incubate for 2-3 days, then change medium. Passage when cells reach 80-90% confluence

Subculture of early culture

40m

3 The cells are fragile the first 2-3 passages. To maximize growth, it is important to use a special medium formulation and take extra care when working.



3.1 Add  11 mL warm CIM in a T-75 flask and incubate for  00:30:00 in a 5% CO₂ incubator

30m


3.2 Remove supernatant from T-75 and wash with PBS

Materials:

 PBS pH 7.2 **Gibco - Thermo Fischer Catalog #20012019**


3.3 Add  3 mL PBS and incubate for  00:05:00

5m

3.4 Add  3 mL CIM and transfer content to centrifugation tube

3.5 Centrifuge with  200 rcf, 00:05:00

5m

3.6 Discard supernatant and resuspend pellet in  1 mL warm CIM

3.7 Split between 1:2 to 1:6 of resuspended suspension from centrifugation tube to T-75 and incubate in a 5% CO₂ incubator

Note

Incubate for 2-3 days, then change medium. Passage when cells reach 80-90% confluence

Repeat this procedure once, then change to the general subculture protocol

Protocol



NAME

General subculture protocol for HEK-Blue cells

CREATED BY

Andreas Sagen

PREVIEW