

NOV 14, 2022



WORKS FOR ME

Luciferase Assay pH-dependent Fox-activity

DOI

dx.doi.org/10.17504/protocols.io.n2bvj833wgk5/v1

kyle.kisor1

<sup>1</sup>UCSF



COMMENTS 0

**ABSTRACT** 

Working protocol for how to perform luciferase assay. Worked for Fox reporter

DOI

dx.doi.org/10.17504/protocols.io.n2bvj833wgk5/v1

PROTOCOL CITATION

kyle.kisor 2022. Luciferase Assay pH-dependent Fox-activity. **protocols.io** <a href="https://dx.doi.org/10.17504/protocols.io.n2bvj833wgk5/v1">https://dx.doi.org/10.17504/protocols.io.n2bvj833wgk5/v1</a>

**KEYWORDS** 

**luciferase** 

LICENSE

This is an open access protocol distributed under the terms of the <u>Creative Commons</u>

<u>Attribution License</u>, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Nov 14, 2022



1

Citation: kyle.kisor Luciferase Assay pH-dependent Fox-activity <a href="https://dx.doi.org/10.17504/protocols.io.n2bvj833wgk5/v1">https://dx.doi.org/10.17504/protocols.io.n2bvj833wgk5/v1</a>

LAST MODIFIED

Nov 14, 2022

PROTOCOL INTEGER ID

72739

## MATERIALS TEXT

- Generic tissue culture materials
- Cell line of choice
- Luciferase plasmid and Renilla control plasmid
- Dual-Glo luciferase kit (promega)

SAFETY WARNINGS

## No safety warnings

**BEFORE STARTING** 

Read entire protocol before starting

## **Plate Cells**

- 1 1. Plate 500K cells/well for each condition
- 2 Incubate o/n before transfection

## **Transfection**

- For each transfection condition create the following tube mixtures. (Volumes are for 1 transfection per condition and can be scaled up)
  - a. Tube 1: 125uL Optimem + 4uL lipofectamine 3000
  - b. Tube 2: 125uL Optimem + 2uL p3000 reagent + 1 ug total DNA
- 4 flick and incubate 5 minutes
- 5 Mix together tube 1 and tube 2 and incubate 15 min



6	Add dropwise to cells and incubate o/n	
7	Incubate 8 hours	
8	Change media to fresh media or media + 10mM EIPA	
	48 hours incubation	
9	Incubate cells for 48 hours after transfection. Change media at 24 hours post transfection with fresh media and fresh media + EIPA	
	Luciferase Assay	
10	Slowly thaw luciferase reagent in room temp water	
11	Wash gently each well with PBS	
12	Add 500uL of luciferase reagent directly to plate and nutate for 10 min RT	
- <del>-</del>		
13	Scrape each well with cell scraper and put into fresh microfuge tube	
protoc	protocols in	

14	Spin down 13,000 rpm 5 min at RT
15	Add 100uL of supernatant to 4wells of white 96well plate. Wait 10 minutes
16	Read luciferase on plate reader
17	For Renilla reading, calculate reagent needed for 100uL each well using a 1:100 dilution of reagent: buffer
18	Add 100uL of renilla reagent to each well and wait 10 min before reading