

Sep 11, 2024

In Vivo Metabolomics

DOI

dx.doi.org/10.17504/protocols.io.kqdg326b7v25/v1

Yoshitaka J Sei¹, Szu-Chi Liao¹, Johanna ten Hoeve-Scott², Ken Nakamura¹

¹Gladstone Institutes; ²UCLA Metabolomics Center

ASAP Collaborative Rese...



Haru Yamamoto

UCSF

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.kqdg326b7v25/v1

Protocol Citation: Yoshitaka J Sei, Szu-Chi Liao, Johanna ten Hoeve-Scott, Ken Nakamura 2024. In Vivo Metabolomics. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kqdg326b7v25/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 25, 2024

Last Modified: September 11, 2024

Protocol Integer ID: 98820

Keywords: ASAPCRN

Funders Acknowledgement:

ASAP

Grant ID: 020529

Abstract

This protocol describes steps for targeting metabolites with [U-13C]glucose in mice for metabolomics analysis.

Tail vein infusion and tissue harvest

1 **Materials**

- Surgical tools – scissors, tweezers, forceps, spatula, razor
- Heating pad (K&H Thermo-Peep Heated Pad)
- Tail illuminator (Braintree Scientific: MTI STD)
- Glue (Loctite 1699233)
- 30g needles (BD: 305106)
- PE tubing (BD: PE#10)
- Syringe pump (WPI: SP220I)
- 1mL syringes (BD: 309628)
- Vapor splitter (VetEquip: CX-IM 5 place manifold)
- Dry ice
- Sterile 1X PBS or saline
- [U-13C]glucose solutions (Cambridge Isotope Laboratories: CLM-1396-1)
- Glucometer (Zoetis: AlphaTrak 2)
- Luer adapter

2 **Day Prior to Procedure**

- 2.1 Weigh mice.
- 2.2 Fast mice overnight (14-16 hours).

3 **Day of Procedure**

- 3.1 Make [U-13C]glucose solutions with PBS/saline.
 - Bolus: 200 μ L of 0.4 mg/g.
 - Infusion: 75 μ L of 0.72 mg/g.
- 3.2 Anesthetize mice in an induction chamber with 2% isoflurane with flowrate set to 1.2 L/min.
- 3.3 Move a mouse to the tail illuminator, cannulate the tail vein and glue needle in place.



- 3.4 Test that the catheter works by infusing 0.1 mL of saline – the saline should flow up the tail vein if it was cannulated correctly.
- 3.5 Gently move mouse and catheter to the isoflurane-splitter.
- 3.6 Repeat steps 3-4 to fill up the remaining spaces on the isoflurane-splitter (up to 5 mice).
- 3.7 Prick the paw of each mouse with a needle (18G) to measure the pre-infusion blood glucose level.
- 3.8 Inject the bolus glucose solution into each mouse and transfer the luer adapter to the infusion syringes.
- 3.9 Start infusion at a flow rate of 150 $\mu\text{L}/\text{h}$ for 75 μL (30 min).
- 3.10 After infusion is done, euthanize mice by cervical dislocation.
- 3.11 Quickly dissect out the brain and flash freeze in isopentane cooled to -80°C or liquid nitrogen.
- 3.12 Store samples in -80°C until further processing.

Metabolite Extraction

- 4 Protocol adapted from [UCLA Metabolomics Center](#).

Materials

- Aluminum foil
- Metal Pad
- Dry ice
- Hammer
- Microcentrifuge tubes
- 80% methanol
- Vortex mixer
- Ice



- Centrifuge
 - Screw-cap glass vials
- 5 Wrap frozen brains in pre-cooled aluminum foil and place on a metal pad on ground dry ice.
- 6 Pulverize tissue with a hammer.
- 7 Transfer the ground tissue into a pre-cooled microcentrifuge tube on dry ice.
- 8 Add 1mL 80% MeOH at -80°C and vortex vigorously for 30 s.
- 9 Incubate samples at -80°C for an hour.
- 10 Allow samples to warm up a bit on ice.
- 11 Vortex samples again for 30 s.
- 12 Centrifuge at top speed (16,000 g) for 15 min at 4°C.
- 13 Transfer the supernatant to a screw-cap glass vial and store at -80°C until shipping to UCLA Metabolomics Core.

Metabolomics Analysis

- 14 Protocol from [UCLA Metabolomics Center](#).

Materials

- Q Exactive mass spectrometer (Thermo Scientific)
- Vanquish Flex UHPLC (Thermo Scientific)
- Ion Chromatography System ICS-5000+
- Luna NH2 3um 100A (150 × 2.0 mm) column (Phenomenex)
- LC-MS grade acetonitrile (ACN)
- LC-MS grade water

- 5 M ammonium acetate
 - ammonium hydroxide
- 15 Dilute the MeOH extract with 50% ACN to 50 mg tissue/ml concentration (1:5-1:10).
 - 16 Spin 10 min at full speed at 4C to clarify the samples.
 - 17 Transfer the supernatant to glass vials with glass insert.
 - 18 Combine 10 ul of each sample to make a pool sample used to condition the column prior loading the individual samples.
 - 19 Place samples in autosampler 4C.
 - 20 Load 5 ul onto a Luna NH2 column using a Vanquish Flex.
 - 21 Separate the metabolites at a flow rate of 200 µl/min with mobile phases A (5 mM NH₄AcO pH 9.9) and B (100% ACN) using a linear gradient from 15% A to 95% A over 18 min followed by a 7 min isocratic wash at 95% A.
 - 22 Re-equilibrate the column to 15% A for 15 min prior to loading the next sample.
 - 23 Acquire data with a Q Exactive mass spectrometer using polarity switching in full scan mode using a range of 70-975 m/z and 70,000 resolution.
 - 24 Convert the RAW files to mzXML format using the msConverter (part of open source ProteoWizard).
 - 25 Extract metabolite intensities with Maven (v 8.1.27.11) using a targeted list of polar metabolites of central carbon metabolism with expected retention time and accurate mass measurements (< 5 ppm) for identification.
 - 26 Perform data analysis using R scripts functions (or other tools). Correct isotopologue data for the natural abundance of C¹³ using AccuCor.