



Apr 30, 2024

Recipe for 50x TAE buffer



Forked from [Recipe for 50x TAE buffer](#)

This protocol is a draft, published without a DOI.

Anna Behle¹, Alice Pawlowski¹

¹Institute for Synthetic Microbiology

Axmann Lab

iGEM Duesseldorf 2018

1 more workspace



Skylar Carlson

University of the Pacific

OPEN  ACCESS



Protocol Citation: Anna Behle, Alice Pawlowski 2024. Recipe for 50x TAE buffer. [protocols.io](https://protocols.io/view/recipe-for-50x-tae-buffer-dbh82j9w) <https://protocols.io/view/recipe-for-50x-tae-buffer-dbh82j9w>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: March 29, 2024

Last Modified: April 30, 2024

Protocol Integer ID: 97568

Abstract

Stock solution for 50x TAE. TAE buffer is a solution made up of Tris base, acetic acid and EDTA (Tris-acetate-EDTA). It is a common buffer for DNA separation using standard agarose gel electrophoresis.



Materials

MATERIALS



EDTA



Acetic acid, glacial **Sigma Aldrich Catalog #537020**



Tris (Tris Base) **Gold Biotechnology Catalog #T-400**



Components required

1	A	B	C	D	E	F	G
		Component	mass	VOLUME	Molarity in 50x	Final in 1x	Molecular weight
	1	EDTA disodium salt	93.05 g		50 mM	1 mM	372.24 g/mol
	na	Dissolve	400 mL	100 mL	pH = 8.0		
		AUTOCLAVE					
	2	Tris-base	242 g		2 M	40 mM	121.14 g/mol
	na	Dissolve	700 mL	700 mL			
	3	glacial / acetic acid		57.1 mL	1 M	20 mM	60.05

Preparation of a 0.5 M EDTA stock solution

- 2 For 500 ml:
- weigh out 93.05 grams of EDTA disodium salt (MW=372.24 g/mol)
 - Dissolve in 400 milliliter deionized water and adjust the pH with solid sodium hydroxide (NaOH) plates, EDTA will not go completely into solution until the pH is adjusted to about 8.0!
 - Top up the solution to a final volume of 500 milliliters
 - autoclave

Preparation of 50 x TAE buffer

- 3
- weigh out 242 grams of Tris-base (MW = 121.14 g/mol) and dissolve in approximately 700 milliliters of deionized water
 - Carefully add 57.1 milliliters of 100 % glacial acid (or acetic acid) and 100 milliliters of 0.5 M EDTA (pH 8.0)
 - adjust the solution to a final volume of 1 liter
 - the pH of this buffer is not adjusted and should be about 8.5
 - store stock solution at room temperature

Preparation of 1 x TAE working solution

- 4
- 20 mL 50 x TAE
add 920 mL DI H₂O

final concentration in gel / running buffer: 40 mM Tris, 20 mM acetic acid, 1 mM EDTA