




Upload image

Jul 25, 2020

Porcine Circovirus 3 (PCV3) Complete Genome Sequence PCR

PLOS One

Chew Yee Tan¹¹Universiti Putra Malaysia**1** Works for me dx.doi.org/10.17504/protocols.io.bdd9i296 Chew Yee Tan
Universiti Putra Malaysia

ABSTRACT

This protocol describes four separate PCR runs to amplify the complete genome sequence of Porcine circovirus 3 (PCV3) as four overlapping segments.

EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0235832>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Tan CY, Opaskornkul K, Thanawongnuwech R, Arshad SS, Hassan L, Ooi PT (2020) First molecular detection and complete sequence analysis of porcine circovirus type 3 (PCV3) in Peninsular Malaysia. PLoS ONE 15(7): e0235832. doi: [10.1371/journal.pone.0235832](https://doi.org/10.1371/journal.pone.0235832)

DOI

dx.doi.org/10.17504/protocols.io.bdd9i296

PROTOCOL CITATION

Chew Yee Tan 2020. Porcine Circovirus 3 (PCV3) Complete Genome Sequence PCR. **protocols.io**
dx.doi.org/10.17504/protocols.io.bdd9i296

MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

Tan CY, Opaskornkul K, Thanawongnuwech R, Arshad SS, Hassan L, Ooi PT (2020) First molecular detection and complete sequence analysis of porcine circovirus type 3 (PCV3) in Peninsular Malaysia. PLoS ONE 15(7): e0235832. doi: [10.1371/journal.pone.0235832](https://doi.org/10.1371/journal.pone.0235832)

EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0235832>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Mar 09, 2020

LAST MODIFIED

Jul 25, 2020

PROTOCOL INTEGER ID

33953

MATERIALS

NAME	CATALOG #	VENDOR
QIAgen DNeasy Blood and Tissue Kit, 50 rxn	69504	Qiagen
MyTaq™ Red Mix 2X master mix		Bioline
RedSafe™ nucleic acid staining solution		iNtRON Biotechnology, South Korea
GelPilot 100 bp Plus Ladder		Qiagen

- 1 DNA extraction was performed using DNeasy Blood & Tissue Kit extraction kit (Qiagen, Germany) in accordance to manufacturer's instructions. Samples involved include lung, inguinal lymph node, spleen, tonsil, kidney, mesenteric lymph node, heart, liver and brain tissues.
- 2 Four separate conventional PCR runs were performed to amplify the complete nucleotide sequence of PCV3 as four overlapping segments.
- 3 Primers, their respective cycling conditions and expected PCR product band size are tabulated below:

Primer Pair	Nucleotide Sequence (5'–3')	Product Length (bp)	PCR Cycling Condition (Temperature / Time)						Reference
			Initial Denaturation	Number of Cycle	Denaturation	Annealing	Extension	Final Extension	
P1F	CACCGTGTGAGTGGATA TAC	854	94°C / 4 min	35	94°C / 20 s	55°C / 30 s	72°C / 30 s	72°C / 5 min	Palinski R, Piñeyro P, Shang P, Yuan F, Guo R, Fang Y, et al. A novel porcine circovirus distantly related to known circoviruses is associated with porcine dermatitis and nephropathy syndrome and reproductive failure. J Virol. 2017 Jan 1;91(1):e01879-16.
P1R	CAAAACCCACCTTAACA G								
P2F	GTCGTCCTGGAGCCAAG TG	807	95°C / 5 min	35	94°C / 30 s	62°C / 30 s	72°C / 1 min	72°C / 7 min	
P2R	CGACCAAAATCCGGGTAA GC								
KF	TTACTTAGAGAACGGAC TTGTAACG	649	94°C / 5 min	35	94°C / 30 s	55°C / 30 s	72°C / 1 min	72°C / 10 min	Ku X, Chen F, Li P, Wang Y, Yu X, Fan S, et al. Identification and genetic characterization of porcine circovirus type 3 in China. Transbound Emerg Dis. 2017 Jun;64(3):703-8.
KR	AAATGAGACACAGAGCT ATATTACG								
V4F	GAAAACGCGGAAGCTT GTG	806	95°C / 5 min	35	94°C / 30 s	56°C / 30 s	72°C / 1 min	72°C / 10 min	Designed in this study based on PCV3 strain KY075986
V4R	CCACTTCTGGCGGGAAC TAC								

- 4 Type, concentration and volume of reagents used in the four separate PCR runs are the same: 12.5 µL of MyTaq™ Red Mix 2X master mix (Bioline, United Kingdom) and 1.0 µM of each primer in a pair, in a 25 µL total PCR reaction volume.
- 5 PCR products were stained using RedSafe™ nucleic acid staining solution (iNtRON Biotechnology, South Korea) and analysed by agarose electrophoresis using 2.5% agarose gel with GelPilot 100 bp Plus Ladder as marker (Qiagen, Germany).