

May 13, 2024

# Protocol for performing PINK1 siRNA knockdown in mouse embryonic fibroblasts (MEFs)

DOI

[dx.doi.org/10.17504/protocols.io.kxygx343zg8j/v1](https://dx.doi.org/10.17504/protocols.io.kxygx343zg8j/v1)

Enrico Bagnoli<sup>1</sup>, Miratul Muqit<sup>1</sup>

<sup>1</sup>Medical Research Council Protein Phosphorylation and Ubiquitylation Unit, School of Life Sciences, University of Dundee, Dow Street, Dundee DD1 5EH, UK

ASAP Collaborative Rese...

Enrico Protocols



Francesca Tonelli

MRC-PPU at The University of Dundee

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.kxygx343zg8j/v1](https://dx.doi.org/10.17504/protocols.io.kxygx343zg8j/v1)

**Protocol Citation:** Enrico Bagnoli, Miratul Muqit 2024. Protocol for performing PINK1 siRNA knockdown in mouse embryonic fibroblasts (MEFs). [protocols.io](https://dx.doi.org/10.17504/protocols.io.kxygx343zg8j/v1) <https://dx.doi.org/10.17504/protocols.io.kxygx343zg8j/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** December 13, 2023

**Last Modified:** May 13, 2024

**Protocol Integer ID:** 92308

**Keywords:** ASAPCRN

**Funders Acknowledgement:**

**ASAP**

**Grant ID:** ASAP-000463

## Abstract

This protocol details the siRNA knockdown in mouse embryonic fibroblasts (MEFs) for PINK1, but applicable for any other target.

## Attachments



[926-2384.docx](#)

17KB

## Guidelines

### 72-Hr siRNA knockdown in MEFs in a 6-well plate

- If using Immortalised MEFs – seed 100,000 cells per well.
- If using Primary MEFs – seed 200,000 cells per well.


## Materials

### Key Reagents

- siRNA as purchased from Horizon Discovery as ON-TARGETplus siRNA Reagents (Dharmacon<sup>TM</sup> siRNA solutions).


#### Note


Protocol for 5 nmol quantities of purchased siRNA

-  Lipofectamine<sup>®</sup>; RNAiMAX Transfection Reagent **Thermo Fisher Catalog #13778150**
- Oligomycin
- Antimycin
- MLI2




## Day 1 - Cell Seeding and siRNA preparation

1 For Primary MEFs seed 200,000 cells per well in a 6-well plate at a total volume of  2 mL .

2 For Immortalised MEFs seed 100,000 cells per well in a 6-well plate at a total volume of  2 mL .

3

## Day 1 - Cell Seeding and siRNA preparation

4 Resuspend the dried siRNA – protocol for 5 nmol of purchased siRNA to reconstitute at [  10 micromolar ( $\mu\text{M}$ ) ].

4.1 Add  400  $\mu\text{L}$  of RNA-free water.



4.2 Add  100  $\mu\text{L}$  of 5x siRNA Buffer.




4.3 Incubate in hood for  00:05:00 , vortex vigorously and store at  -20 °C .

5m



### Note

The protocol can be amended if the nmol quantity of siRNA purchased is higher/lower by adjusting the resuspension to ensure a final [  10 micromolar ( $\mu\text{M}$ ) ] concentration.














## Day 2 - siRNA Knock-Down

5m

5 Prepare mastermix according to the number of wells. Four tubes have to be prepared, 2 with OptiMem and lipofectamine and one each with the PINK1 and scramble siRNA respectively.

**Note**

For PINK1 knockdown a final concentration of 25 nM of siRNA in each well is used.

- 5.1 Tube1 PINK1 siRNA:  
Per well -  5  $\mu$ L of PINK1 siRNA at [  10 micromolar ( $\mu$ M) ] diluted in  100  $\mu$ L OPTI-MEMTube 3
- 5.2 Tube 2:  
Per well -  10  $\mu$ L of Lipofectamine diluted in  100  $\mu$ L of siRNA-OPTI-MEM
- 5.3 Tube 3 scramble siRNA  
Per well -  5  $\mu$ L of scramble siRNA at [  10 micromolar ( $\mu$ M) ] diluted in  100  $\mu$ L OPTI-MEM
- 5.4 Tube 4 (same as tube 2)  
Per well -  10  $\mu$ L of Lipofectamine diluted in  100  $\mu$ L of siRNA-OPTI-MEM
- 6 Vortex slowly and combine the content of Tube 1 with tube 2 (PINK1 siRNA) and similarly with tube 3 and 4.
- 7 Vortex again slowly and incubate for  00:05:00 at  Room temperature .
- 8 Add drop by drop,  200  $\mu$ L of PINK1 siRNA mix to each well for PINK1 KD and similarly of the scramble siRNA for controls.

5m



**Note**

For an experiment with 4x 6 well plates (WT and mutant cells treated with scramble or PINK1)  
the following volumes can be used:

- Tube1 and tube 3: 60  $\mu\text{L}$  siRNA + 1200  $\mu\text{L}$  OPTiMEM
- Tube 2 and 4: 120  $\mu\text{L}$  of lipofectamine + 1200  $\mu\text{L}$  of otpimem
- This gives 2580  $\mu\text{L}$  of each mix ( 180  $\mu\text{L}$  spare)

**Day 4 - Mitochondrial depolarization**

5m

9 Make a 500x of Oligomycin/antimycin solution. The final concentration in the well is:

- oligo: 1 micromolar ( $\mu\text{M}$ )
- antimycin: 10 micromolar ( $\mu\text{M}$ )

9.1 Add 4  $\mu\text{L}$  of this 500x solution will be to each well.

**Note**

For a 100  $\mu\text{L}$  of 500x solution, we use:

A	B
Compound	Quantity
Antimycin A 50 mM solution	10 $\mu\text{L}$
Oligomycin 6.4 mM solution	7.8 $\mu\text{L}$
DMSO	82.2 $\mu\text{L}$



- 10 Treat cell with Oligomycin/antimycin A for ⌚ 24:00:00 . OA should be added ⌚ 48:00:00 after the addition of the siRNA. Include DMSO control.

3d



## Day 5 - Cell lysis

1h 50m

- 11 Prepare working stock of MLI2 at a concentration of [M] 100 micromolar ( $\mu\text{M}$ ) in DMSO.

- 12 Treat cells with MLI2( [M] 10 nanomolar (nM) ) for ⌚ 01:30:00 (2 ul/well).

1h 30m



- 13 Lyse cells using 🧪 50  $\mu\text{L}$  of Lysis buffer/well



- 14 Lysate can be precleared by centrifugation 🌀 17000 x g for ⌚ 00:15:00 at 🌡 4 °C .

15m



- 15 Perform protein estimation and subject lysate to immunoblotting.