



Easy Covid LAMP V.2

ahmet ¹, ayten ¹

¹Hayat

Version 2 ▾

Dec 05, 2020

1

Works for me

dx.doi.org/10.17504/protocols.io.bqfzmt6

Hayat

ahmet

DOI

dx.doi.org/10.17504/protocols.io.bqfzmt6

PROTOCOL CITATION

ahmet , ayten 2020. Easy Covid LAMP. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bqfzmt6>

Version created by ahmet

WHAT'S NEW

We added a step for fluorimetric assessment of the results to make it a more robust test.

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Dec 05, 2020

LAST MODIFIED

Dec 05, 2020

PROTOCOL INTEGER ID

45273

- 1 Collect sample into 1.5ml eppendorf tube.
- 2 Insert sample transfer plastic into the sample tube, mix for 5 seconds, then insert into pre-filled 0.2ml master mix PCR tube, mix for 5 seconds, discard.
- 3 Place PCR tube on the pre-heated heat block device and run . 🔥 **65 °C (Heat block is pre-set)**
🕒 **00:29:00**
- 4 Check the color of the PCR tube with naked eye: positive if yellow, negative if pink or orange.

- 5 Optional: Check the fluorescence signal using handheld fluorescent viewer: positive if bright green, negative if matte green.