





Sep 20, 2022

RNA extraction, cDNA synthesis and Taqman qPCR

[gurvir.virdi](#)¹¹UCL Institute of Neurology

1 Works for me

 Sharedx.doi.org/10.17504/protocols.io.36wgqj5rkvk5/v1 gurvir.virdi

ABSTRACT

RNA extraction, cDNA synthesis and TaqMan qPCR on samples.

DOI

dx.doi.org/10.17504/protocols.io.36wgqj5rkvk5/v1

PROTOCOL CITATION

gurvir.virdi 2022. RNA extraction, cDNA synthesis and Taqman qPCR.
protocols.io
<https://protocols.io/view/rna-extraction-cdna-synthesis-and-taqman-qpcr-cgvstw6e>



LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 20, 2022

LAST MODIFIED

Sep 20, 2022

PROTOCOL INTEGER ID

70290

RNA extraction

- 1 Cell pellets are snap-frozen using dry ice.

- 2 RNA is extracted using the Maxwell® RSC simply RNACells kit (Promega), and the Maxwell® RSC 48 instrument, following manufacturer instructions.
- 3 After RNA extraction, RNA concentration and quality are measured and assessed using a nanodrop.

cDNA synthesis

- 4 Up to 1 µg of RNA per sample is reverse-transcribed into cDNA using the High-Capacity cDNA Reverse Transcription kit (ThermoFisherScientific).

quantitative PCR

- 5 qPCR is performed using TaqMan™ Gene Expression Assay (Thermo Fisher Scientific) according to the manufacturer's instructions.
- 6 For each gene, TaqMan™ probes were used along with the TaqMan™ master mix, and sample cDNA following the manufacturer's protocol.
- 7 Samples, along with a minus reverse transcriptase control (-RT) were ran for each gene on the QuantStudio 6 Flex Real-Time PCR System (Applied Biosystems).
- 8 The -RT served as a negative control, and the gene expression levels were normalised to the housekeeping gene GAPDH following the delta-delta Ct method. Gene expression values were expressed as the normalisation to their respective control sample.