

VERSION 1

MAR 30, 2023

OPEN ACCESS

DOI:
dx.doi.org/10.17504/protocols.io.q26g7yzz3gwz/v1

Protocol Citation: Antoine Champie, Amélie De Grandmaison 2023. HTTM : DNA Extraction. [protocols.io](https://dx.doi.org/10.17504/protocols.io.q26g7yzz3gwz/v1) <https://dx.doi.org/10.17504/protocols.io.q26g7yzz3gwz/v1> version created by Antoine Champie

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
 We use this protocol and it's working

Created: Oct 27, 2022


Last Modified: Mar 30, 2023

PROTOCOL integer ID:
 71938

HTTM : DNA Extraction V.1

Antoine Champie¹, Amélie De Grandmaison¹

¹Université de Sherbrooke



Antoine Champie

ABSTRACT

Part two of the HTTM protocol. A low-cost and high-throughput Tn-seq protocol. This part cover the DNA extraction from cell pellets of transposon insertion mutants and subsequent silica columns regeneration.

MATERIALS

■ Homemade DNA lysis Buffer :

A	B
Component	Amount for 1000 ml of solution
CTAB 2%	20 g
1.5 M Guanidine HCl	143.2 g
10 mM Tris HCl	1.57 g

Mix well and adjust volume to 1 l with water and adjust pH to 8.0.

■ Homemade wash solution :

A	B
Component	Amount for 1000 ml of solution
Ethanol 100%	800 ml
Tris HCl 1 M pH 8.0	10 ml
NaCl 4 M	25 ml
EDTA 0.5 M	2 ml

Mix well and adjust volume to 1 l with water and adjust pH to 8.0.

■ Elution Buffer (Low TE Buffer): 10 mM Tris-HCl (pH 8.0) + 0.1 mM EDTA

Solutions for plate regeneration, from this protocol :

(1)<https://doi.org/10.1016/j.ab.2008.10.021>.

■ NaOH 1N + Triton X-100 0.15% (v/v)

A	B
Component	Amount for 1000 ml of solution
Water	960 ml
NaOH	40 g
Triton X-100	1.5 ml

Mix well and store in a base resistant container.

■ HCl 1.5N + Triton X-100 0.15% (v/v)

A	B
Component	Amount for 1000 ml of solution
Water	873.5 ml
HCl Stock (37%)	125 ml
Triton X-100	1.5 ml




Mix well and store in an acid resistant container.



Silica columns array come from the following commercially available kit :

96-Well Plate Bacteria Genomic DNA Miniprep Kit from Biobasic. CAT#: SK1295

DNA extraction

2h 5m

- 1 Prepare the lysis solution by adding  165 µL of proteinase K to  66 mL of homemade lysis buffer and mix well.
- 2 Add  600 µL of lysis solution to each well of the deep-well plate and resuspend the pellet.

3 Cover with an adhesive aluminum cover and incubate at  55 °C for  01:00:00 . 1h


4 While still warm, add  260 µL of ethanol 100%, without overmixing.

Note

Overmixing will result in DNA agglomeration and difficulty with the extraction.

5 Transfer immediately to a deep-well plate fitted with an array of silica columns.





6 Centrifuge twice at  4000 x g, 00:10:00 . 10m

7 Discard flowthrough and add  500 µL of wash solution.

8 Centrifuge at  3000 x g, 00:10:00 . 10m


8.1 Repeat steps 7 and 8.

9 Discard flowthrough.

- 10 Centrifuge at  3000 x g, 00:10:00 to eliminate traces of wash solution. 10m
- 11 Discard flowthrough.
- 12 Add a collector plate between the silica column array and the deep-well plate.
- 13 Add  50 µL of low TE to the silica matrix in each well.
- 14 Cover with an adhesive aluminum foil and incubate at  55 °C for  00:15:00 . 15m
- 15 Centrifuge at  3000 x g, 00:10:00 . 10m

Silica array regeneration (Optional)

1 h 5m

- 16 Put the contaminated silica array on an empty deep-well plate.
Add  150 µL of 1N NaOH + 0.15%(v/v) Triton X-100 to each well.

- 17 Incubate at  Room temperature for  00:05:00 . 5m
- 18 Centrifuge  3000 x g, 00:05:00 . 5m
- 19 Add  200 μ L of 1.5N HCl+ 0.15% (v/v) Triton X-100 to each well.
- 20 Incubate at  Room temperature for  00:30:00 . 30m
- 21 Centrifuge  3000 x g, 00:05:00 . 5m
- 22 Add  150 μ L of 1N NaOH + 0.15%(v/v) Triton X-100 to each well.
- 23 Incubate at  Room temperature for  00:05:00 . 5m
- 24 Centrifuge  3000 x g, 00:05:00 . 5m
- 24.1 Collect the flowthrough in a beaker. Neutralize pH if needed and dispose of the flow through.

25 Add  200 µL of ddH₂O to each well.

26 Centrifuge  3000 x g, 00:05:00 .

5m

26.1 Repeat steps 25 and 26.

27 Silica columns array are ready to be reused.