



Sep 12, 2022

# Final QC, Pooling and Sequencing

Katarina A Cohen<sup>1,2</sup>, Khai-Minh H Nguyen<sup>1,2</sup>, Oksana Polesskaya<sup>1,2</sup>,  
Abraham Palmer<sup>1,2</sup>

<sup>1</sup>UC San Diego; <sup>2</sup>Palmer Lab

1 Works for me

Share

[dx.doi.org/10.17504/protocols.io.yxmvmnw29g3p/v1](https://dx.doi.org/10.17504/protocols.io.yxmvmnw29g3p/v1)

CGORD

Palmer Wet-Lab Protocols



Khai-Minh H Nguyen

## ABSTRACT

This protocol is conducted after a set of libraries are completed and ready to quantify and pool. This protocol outlines the final steps before submitting for sequencing.

## DOI

[dx.doi.org/10.17504/protocols.io.yxmvmnw29g3p/v1](https://dx.doi.org/10.17504/protocols.io.yxmvmnw29g3p/v1)

## PROTOCOL CITATION

Katarina A Cohen, Khai-Minh H Nguyen, Oksana Polesskaya, Abraham Palmer  
2022. Final QC, Pooling and Sequencing . **protocols.io**  
<https://protocols.io/view/final-qc-pooling-and-sequencing-b8i8ruhw>



## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

May 01, 2022

## LAST MODIFIED

Sep 12, 2022

## PROTOCOL INTEGER ID

61760

## MATERIALS TEXT

### Equipment

- PC Running Excel
- Pipette 1-10uL
- Pipette 20-200uL
- Bioanalyzer
- Qubit Fluorometer

### Reagents

-  [Ultrapure Distilled, Nuclease Free Water](#) **Contributed by users** In 2 steps

### Consumables

- Pipette Tips
- 1.5mL Eppendorf Tubes
- Qubit Assay
- TapeStation D1000

## BEFORE STARTING

Ensure all libraries that will be pooled are uniquely indexed.

### Library QC

- 1
  - Quantify purity and concentration of library with Nanodrop and a Qubit Assay
  - Obtain average fragment size of library with TapeStation (D1000 Assay)

- 1.1
  - Libraries should have an average fragment size between 420bp - 650bp.
  - 260/280 should be around 1.80 - 2
  - 260/230 should be around 2-2.2.
  - We have been able to get good data from libraries with relatively poor nanodrop purities.
  - Qubit concentrations can widely range. We get a range from 10ng/ul - 60ng/ul

### Pooling

- 2 Download  [Pooling Template.xlsx](#)

- 3
  - Increase "Target Vol (uL) per Sample" if any Sample Vol is lower than 1ul

- Use [☒ Ultrapure Distilled, Nuclease Free Water Contributed by users](#) when adding water to final pool

4 Enter library name in Column A "Sample Name"

5 Enter qubit concentrations in column B

6 Enter average fragment size of library in column C for each library.

7 Increase "Target Vol (uL) per Sample" if any Sample Vol is lower than 1ul

8 Label new tube as Riptide Pool ## , and Date. Add calculated volume of water (shown in cell I6). Add calculated volumes of sample shown in column E.

- Use [☒ Ultrapure Distilled, Nuclease Free Water Contributed by users](#) when adding water to final pool

#### Checking Pooling

### 9 RECOMMENDED OPTIONAL STEP

Check pooling with an illumina MiSeq run

- % Reads Identified for each library shouldn't vary more than 30% from each other.

#### Sample-Barcode List

10

[📄 YYYY-MM-DD-Flowcell Sample-Barcode list.xlsx](#)

- Download the file above
- Check Twist Bioscience site for updates to the sample barcodes used for the Twist 96-Plex Kit.

11 Open the library file created in the "EPMotion - Normalization and Randomization".

- Go to the "Sample\_Randomization" Tab

- 12 When the "Sample\_Randomization" tab is opened, copy and paste the randomized Transponder ID's into Column A of the Flowcell Sample-Barcode list.
  - NOTE: If less than 96 samples are processed in a library, delete the unfilled rows within that 96 library set. (Essentially, you want to delete the sample barcodes that are not associated with a sample).
- 13 Transfer any comments from the library file into Column G of the Sample barcode list
  - Ensure that the comment matches with the correct sample ID.
- 14 Enter the PCR index barcode used for that particular library (ex. 1,2,3,4,5...) in Column D
- 15 Enter Library name in Column E

#### Submitting for Sequencing

- 16 Submit 30ul of Pool
  - Platform: NovaSeq S4
  - Run Type: PE150

You can find the updated inline i7 and i5 index sequences on Twist Bioscience Site