



Jul 27, 2021

# genomic DNA-extraction

**Michael Burgis**<sup>1</sup><sup>1</sup>Philipps University Marburg, Zentrum für synthetische Mikrobiologie

1 Works for me

Share

[dx.doi.org/10.17504/protocols.io.bwx3pfqn](https://dx.doi.org/10.17504/protocols.io.bwx3pfqn)**Michael Burgis**

## ABSTRACT

DNA isolation protocol for quick and rapid extraction of genomic DNA from *Chlamydomonas reinhardtii* cells. (also works for cyanobacteria)

This protocol is based on a paper published by Newman et al 1990 (PMID: [1981764](#)).

The protocol is designed to be used for single colonies or 1ml of cell culture but can be used for a larger amount of biomass if scaled correctly. Here I describe the protocol for 10ml of overnight culture.

## DOI

[dx.doi.org/10.17504/protocols.io.bwx3pfqn](https://dx.doi.org/10.17504/protocols.io.bwx3pfqn)

## PROTOCOL CITATION

Michael Burgis 2021. genomic DNA-extraction . **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bwx3pfqn>



## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Jul 27, 2021

## LAST MODIFIED

Jul 27, 2021

## PROTOCOL INTEGER ID

51931

## MATERIALS TEXT

TEN buffer (10mM Tris-HCl, 10mM EDTA [Ethylenediaminetetraacetic acid], 150mM NaCl); SDS-EB (2% SDS, 400mM NaCl, 40mM EDTA, 100mM Tris-HCl, pH8.0), Nuclease free water, phenol/Chloroform:isoamyl alcohol(1:1), Chloroform:isoamyl alcohol, 70% Ethanol, 100% Ethanol, TE Buffer (10mM Tris-HCl pH 8.0, 1mM EDTA pH 8.0)

### Buffer preparation

- 1
  - TEN Buffer (10 mM Tris-HCl, 10 mM EDTA, 150 mM NaCl)
  - SDS-EB Buffer (2% SDS, 400 mM NaCl, 40 mM EDTA, 100 mM Tris-HCl, pH 8.0)
  - TE Buffer (10mM Tris-HCl pH 8.0, 1mM EDTA pH 8.0) [autoclaved]

### Resuspension and lysis

- 2 Pellet the 10ml overnight culture of your algae or cyanobacteria in a 15ml falcon at 5.000g.

- 3 Resuspend in 5ml of TEN Buffer by vortexing, spin for 10 sec. and aspirate off supernatant.
- 4 Resuspend cells in 150 ul H<sub>2</sub>O on ice and add 300 ul of SDS-EB buffer, vortex to mix.

#### Phenol Chloroform washing

- 5 Extract once with 3.5 ml phenol/Chloroform:isoamyl alcohol(1:1) for a few min by vortexing, separate phases by centrifugation for 5 min., transfer aq. phase to a new falcon tube.
- 6 Extract once with 3 ml Chloroform:isoamyl alcohol (24:1), transfer aq. to a new tube.
- 7 Repeat step in order to maximize purity.

#### Precipitation and Resuspension

- 8 Add 2 volumes abs. ethanol, incubate on ice for 30 min., centrifuge for 10 min., wash pellet once with 200 ul 70% ethanol.
- 9 Dry at 60°C until all ethanol evaporated and resuspend in desired volume of TE Buffer.
- 10 Incubate at 50 °C for 1 to 2 h to ensure complete re-suspension. (Can be extended to overnight incubation)
- 11 Measure the quality and quantity of your DNA using a NanoDrop.