

Sep 19, 2024

Hepatitis E and Enterovirus Multiplex Nested RT-PCR

DOI

dx.doi.org/10.17504/protocols.io.rm7vzjd74lx1/v1

Shannon Fitz¹, Alex Shaw¹, Dilip Abraham², michael Owusu³

¹Imperial College London; ²CMC Vellore; ³Kwame Nkrumah University of Science and Technology

GEMS - Genomic Environ...



Shannon Fitz

Imperial College London

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.rm7vzjd74lx1/v1

Protocol Citation: Shannon Fitz, Alex Shaw, Dilip Abraham, michael Owusu 2024. Hepatitis E and Enterovirus Multiplex Nested RT-PCR. [protocols.io](https://dx.doi.org/10.17504/protocols.io.rm7vzjd74lx1/v1) <https://dx.doi.org/10.17504/protocols.io.rm7vzjd74lx1/v1>

Manuscript citation:

Hepatitis E primers extracted from: Wang B, Harms D, Papp CP, Niendorf S, Jacobsen S, Lütgehetmann M, Pischke S, Wedermeyer H, Hofmann J, Bock CT. Comprehensive Molecular Approach for Characterization of Hepatitis E Virus Genotype 3 Variants. J Clin Microbiol. 2018 Apr 25;56(5):e01686-17. doi: 10.1128/JCM.01686-17. PMID: 29514938; PMCID: PMC5925713.

Enterovirus primers: designed by Dr Alex Shaw, unpublished.

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 19, 2024

Last Modified: September 19, 2024

Protocol Integer ID: 107998

Keywords: PCR, Hepatitis E, Enterovirus, Multiplex PCR, Nested PCR

Funders Acknowledgement:

The Gates Foundation

Grant ID: INV-049092 PA4273

Abstract

For multiplex, nested amplification of Hepatitis E (HEV) and Enterovirus (EV) samples.

Primer sets and target amplicon size

Outer HEV: HEV_38_f and HEV_39_r (538 bps)

HEV_38_f: GARGCYATGGTBGAGAARG

HEV_39_r: GCCATRTTCCARACRGTRTTCC

Inner HEV: HEV_37_f and HEV_27_r (306 bps)

HEV_37_f: GGTTYCGYGCYATTGARAARG

HEV_27_r: TCRCCRGARTGYTTCTTCC

Outer EV: For1_5NCR and Rev7_5NCR (535 bps)

For1_5NCR: CGGTACCYTTGTRCGCCTG


Rev7_5NCR: ATTGTCACCATAAGCAGCC


Inner EV: For2_5NCR and Rev6_5NCR (388 bps)

For2_5NCR: CAAGCACTTCTGTTWCCC

Rev6_5NCR: CCAAAGTAGTCGGTCCGC

Materials

 DreamTaq PCR Master Mix (2X) **Thermo Fisher Catalog #K1072**

 SuperScript III One-Step RT-PCR System with Platinum Taq **Invitrogen - Thermo Fisher Catalog #12457-026**

Approximate total cost per sample: £8.01

8 primers required (cost excluded from estimate as primers do not need to be ordered each time)

Superscript cost per sample: ~£7.56


Dreamtaq cost per sample: ~£0.45 per sample


Extra equipment required:

Vortex, mini centrifuge, thermocycler



Protocol materials

 DreamTaq PCR Master Mix (2X) **Thermo Fisher Catalog #K1072** Materials, Step 1.3


 SuperScript III One-Step RT-PCR System with Platinum Taq **Invitrogen - Thermo Fisher Catalog #12457-026**

Materials, Step 1.1

Hepatitis E and Enterovirus Multiplex Nested PCR

1 First round PCR for combined enterovirus and hepatitis E samples

1.1

 SuperScript III One-Step RT-PCR System with Platinum Taq **Invitrogen - Thermo Fisher Catalog #12457-026**

Prepare master mix on ice per number of samples/controls + 10%. Flick master mix gently to mix, and spin down. Aliquot 20 ul into each well, then add 5 ul of template.

A	B
Reagent	Volume for 1 rxn (ul)
2x reaction mix	12.5
HEV_38_f primer	0.5
HEV_39_r primer	0.5
For1_5NCR primer	0.5
Rev7_5NCR primer	0.5
Nuclease free water	4.5
SSIII taq	1
Template	5 ul
Total volume	25 ul

1.2 Cycling conditions for first round PCR

Flick samples gently to mix, spin down tubes. Amplify samples using the following conditions:



A	B	C
Number of cycles	Temperature (°C)	Time
1	50	30 minutes
1	94	5 minutes
40	94	30 seconds
	52	30 seconds
	68	45 seconds
1	68	5 minutes
	10	Hold

1.3 Second round PCR

 DreamTaq PCR Master Mix (2X) **Thermo Fisher Catalog #K1072**

Prepare on ice per number of samples/controls + 10%. Vortex briefly, and spin down. Aliquot 24 ul of master mix into each well, then add 1 ul of template.

A	B
Reagent	Volume per 1 reaction (ul)
Dreamtaq	12.5
HEV_37_f	0.5
HEV_27_r	0.5
For2_5NCR	0.5
Rev6_5NCR	0.5
Nuclease free water	9.5
Template (amplicon from first round PCR)	1
Total	25 ul

1.4 Cycling conditions for second round PCR

Vortex samples briefly, spin down, and amplify using the following conditions:



A	B	C
Number of cycles	Temperature (oC)	Time
1	95	10 minutes
40	95	30 seconds
	52	30 seconds
	72	45 seconds
1	72	5 minutes
	10	Hold

1.5 Check all samples and controls on an agarose gel or a tapestation.

Protocol references

Wang B, Harms D, Papp CP, Niendorf S, Jacobsen S, Lütgehetmann M, Pischke S, Wedermeyer H, Hofmann J, Bock CT. Comprehensive Molecular Approach for Characterization of Hepatitis E Virus Genotype 3 Variants. J Clin Microbiol. 2018 Apr 25;56(5):e01686-17. doi: 10.1128/JCM.01686-17. PMID: 29514938; PMCID: PMC5925713.