



May 25, 2022

# GUV preparation

 GUV preparation and assayLiv Jensen<sup>1</sup>, Chunmei Chang<sup>2</sup><sup>1</sup>Hurley Lab, University of California, Berkeley; <sup>2</sup>Team Hurley, University of California, Berkeley

1

[dx.doi.org/10.17504/protocols.io.3byl4b398vo5/v1](https://dx.doi.org/10.17504/protocols.io.3byl4b398vo5/v1)

Liv Jensen

GUV preparation for membrane tube assay application

DOI

[dx.doi.org/10.17504/protocols.io.3byl4b398vo5/v1](https://dx.doi.org/10.17504/protocols.io.3byl4b398vo5/v1)

Liv Jensen, Chunmei Chang 2022. GUV preparation. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.3byl4b398vo5/v1>

[GUV preparation and assay, Chunmei Chang](#)

ASAPCRN

 protocol ,

May 02, 2022

May 25, 2022

61817

## 1 GUV Preparation

1.1 Clean the coverglass.

1.2 Coat cleaned coverslips with 100  $\mu$ L 5% (w/w) polyvinyl alcohol (PVA) with a

molecular weight of 145,000 (Millipore).

- 1.3 Place the coated coverslip in a heating incubator at 60 °C to dry the PVA film for 30 min.
- 1.4 Spread a lipid mixture with a molar composition of 70% DOPC, 20% DOPE, 5% DO-PI(3)P, 5% DOPS, 0.3% Atto647N DOPE, 0.01% PEG2000-biotin-DSPE at 1 mg/ml uniformly onto the PVA film.
- 1.5 Put the lipid-coated coverslip under vacuum overnight to evaporate the solvent.
- 1.6 Use 100 µL 320 mOsm sucrose solution for swelling for 1 h at room temperature
- 1.7 Harvest the GUVs and use them with 12 h.