


MAR 20, 2023

## 🌐 Synaptic immunohistochemistry - wholemount via acetone permeabilization

Anya Suppermpool<sup>1</sup>, FishFloorUCL<sup>1</sup>

<sup>1</sup>University College London



Anya Suppermpool

### ABSTRACT

Whole-mount Immunohistochemistry – acetone permeabilization

Works with anti-MAGUK ab.

Modified from M Westerfield protocol.

Lavinia Sheets: T. Nicolson Lab September, 2010

Revised August 2011

Anya Suppermpool Synaptic (MAGUK) Version 2018

Anya's Benchling Protocol Archive at:

<https://benchling.com/s/prt-nLJhuyARGTPcTi887wL1?m=slm-FGIOfYoYTo8bWHpVnUZ0v>

### OPEN ACCESS

#### DOI:

[dx.doi.org/10.17504/protocols.io.261ge38q7l47/v1](https://dx.doi.org/10.17504/protocols.io.261ge38q7l47/v1)

**Protocol Citation:** Anya Suppermpool, FishFloorUCL 2023. Synaptic immunohistochemistry - wholemount via acetone permeabilization. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.261ge38q7l47/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**Created:** Mar 08, 2023

**Last Modified:** Mar 20, 2023

**PROTOCOL integer ID:**  
78330

## MATERIALS

### PO<sub>4</sub>buffer

8 parts 0.1M NaH<sub>2</sub>PO<sub>4</sub>

2 parts 0.1M Na<sub>2</sub>HPO<sub>4</sub>

### BT buffer (make fresh)

1.0g sucrose

18.75µL 0.2M CaCl<sub>2</sub>

PO<sub>4</sub>buffer to 15mL

### BT fix (make fresh)

1.2mL PO<sub>4</sub>buffer

4.8mL BT buffer

2mL 16% PFA (used bought Thermofisher fix)

### PBS/BSA/DMSO

50mL PBS

0.5g BSA

500µL DMSO

## Reagents

- 1** PO<sub>4</sub>buffer  
8 parts 0.1M NaH<sub>2</sub>PO<sub>4</sub>  
2 parts 0.1M Na<sub>2</sub>HPO<sub>4</sub>
- 2** BT buffer (make fresh)  
1.0g sucrose  
18.75µL 0.2M CaCl<sub>2</sub>  
PO<sub>4</sub>buffer to 15mL
- 3** BT fix (make fresh)  
1.2mL PO<sub>4</sub>buffer  
4.8mL BT buffer  
2mL 16% PFA (used bought Thermofisher fix)
- 4** PBS/BSA/DMSO  
50mL PBS  
0.5g BSA  
500µL DMSO

## Day 1

- 5 Tricaine, Dechorionated beforehand (2dpf)
- 6 Fix 1.5-2h in BT fix (\*use bought fix) at 4C
- 6.1 (less fix time higher SNR but mushy) 5h fixation works for 5dpf – 9dpf larvae, time may need to be adjusted for different ages
- 7 Replace fixative with PO4 buffer.  
(Optional) Store at 4C overnight (shaking not necessary)
- 8 Wash for 5' RT with shaking if continuing with permeabilization and blocking in same day
- 9 Chill a small volume (~50mL glass bottle) acetone at -20C at least 20' prior to perm steps.  
(optional) if rushed/forgotten, put acetone in -80C for a few minutes, but this could lead to over-permeabilization
- 10 Transfer larvae to glass vials using glass Pasteur pipette.  
Larvae will be sticky so be careful while transferring between tubes!
- 11 Wash with dH2O 5' RT with rocking.  
\*Timing is critical for permeabilization steps\*—process tubes in same order for each step and adjust volume of solutions to allow for quick handling

- 12 Wash with cold acetone 5' at -20C (put in freezer)
- 13 Wash with dH2O 5' RT with rocking.  
This step can go slightly longer (~10') if necessary
- 14 Wash with PO4, 5' at RT (Lavinia's original protocol extra step)
- 15 Block >2h RT in PBS/BSA/DMSO + 2% goat serum.  
(optional) can go overnight at 4C
- 16 Incubate in primary antibodies in PBS/BSA/DMSO 4C overnight.  
Use ~1mL antibody mix per tube

## Day 2

- 17 Remove primary antibody mix.  
Add 1:1000 20%NaN<sub>3</sub> to mix if saving for future use
- 18 Wash 5X+ with PBS/BSA/DMSO >20' RT with rocking.  
(optional) any washes after antibody incubation can go overnight at 4C if necessary
- 19 Incubate with secondary antibodies in PBS/BSA/DMSO 2-3h at RT or overnight at 4C in the dark

## Day 3

**20** Wash 3X with PBS/BSA/DMSO >20' RT with rocking  
(Wash with PBS/BSA/DMSO + DAPI (1:2000 dilution) >20' RT with rocking)

**21** Wash with PBS >20' RT with rocking like 5x

**22** Move to glycerol – 20%, 40%, 60%, 80%

**23** Mount and image using glycerol lens (did 20x before)

**24**

Reagents	B	C
Anti-MAGUK	1:500	MABN72 Anti-pan-MAGUK, clone K28/86
Fix	#28906	ThermoFisher (from Ana Faro)
tRFP	1:500	Anti-tRFP Rabbit Polyclonal (AB233 - evrogen)