



VERSION 1

JUN 01, 2023

OPEN  ACCESS

Protocol Citation: Lakme Caceres 2023. Nuclei Isolation and Immunoprecipitation for 10X Sequencing . **protocols.io** <https://protocols.io/view/nuclei-isolation-and-immunoprecipitation-for-10x-s-cue6wthe>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development
We are still developing and optimizing this protocol.

Created: May 18, 2023

Last Modified: Jun 01, 2023

PROTOCOL integer ID:
82110

Nuclei Isolation and Immunoprecipitation for 10X Sequencing V.1

Lakme
Caceres¹

¹Princeton Neuroscience Institute



Lakme Caceres

ABSTRACT

This protocol is for assessing the MIT AAVs from pooled injections.

GUIDELINES

Keep tissue/nuclei on ice as much as possible.

MATERIALS

Dynabeads - 4C
anti-GFP - 4C (roche)
Triton X 100 - 4C
DTT - 4C
RNAsin -20C

Prepare Stock Solutions

- 1 Make 20 mL **10% BSA** by combining 2 mL of BSA with 18 mL of MilliQ water in a 50 mL falcon tube. (4°C - 2 weeks)

- 2 Make 20 mL **10% Triton X-100** by combining 18 mL MilliQ water with 2 mL Triton X-100 in a 50 mL tube. Vortex and then incubate at room temperature for 20 minutes. Filter it through a 0.22 μ m filter with a syringe into a clean 50 mL tube. (4°C - 1 month)
- 3 Make 250 mL **Nuclear Isolation Media** by filling a 250 mL bottle with 200 mL of MilliQ water and then adding 2.5 mL 1M Tris, 6.26 mL 1M KCl, 1.25 mL 1M MgCl₂, and 21.45 g Sucrose. Shake until sucrose is dissolved then fill to 250 mL with MilliQ water. (4°C - 2 weeks)
- 4 Make 50 mL **Citrate-Phosphate Buffer** (pH 5) by adding 0.48 g Citric Acid and 0.91 g Dibasic Sodium Phosphate Dihydrate to a falcon tube and then filling it to 50 mL with MilliQ water. Titrate the pH with NaOH using the pH meter. (4°C - 1 month)

Prepare Fresh Solutions

- 5 Make 3 mL **Homogenization Buffer** per sample by adding 2.895 mL Nuclear Isolation Media (filtered via syringe) to a 5 mL eppendorf. Then add 3 μ L 100 mM DTT and 30 μ L 10% Triton X-100. Add 15 uL RNAsin and invert to mix. Store on ice.
- 6 Make 200 uL **Blocking Buffer** per sample by dividing your total desired volume of blocking buffer by 10 to get the amount of 10% BSA in uL. Add this amount to a tube and then fill the remainder with 1X PBS.

Homogenization

- 7 Clean dounce, scalpel, and forceps using MilliQ water, ethanol, RNase Zap, then MilliQ again. The red-tape forceps are for unfixed tissue.
- 8 Get tissue sample from -80°C freezer and place on dry ice. Weigh it on a sterile, tared weigh boat.
- 9 Add tissue to dounce and push it to the bottom using 1 mL of Homogenization Buffer and the pestle. Homogenize the tissue without creating bubbles. Then add the remaining 2 mL of the

Homogenization Buffer and continue to dounce until homogenized.

- 10 Pass all of the nuclei suspension through three FlowMi filters, 1000 μ L at a time into a new 5 mL eppendorf.
- 11 Centrifuge at 900 g/rcf for 10 minutes at 4°C.
- 12 Discard the supernatant and resuspend the pellet in 210 μ L Blocking Buffer. Incubate for 10 minutes on ice.

Cell Count

- 13 Add 9 μ L of sample to a PCR tube and then add 1 μ L of acridine orange.
If sample is very concentrated, instead add 2 μ L sample to 2 μ L of acridine orange and 16 μ L 1X PBS.
- 14 Pipette mix and then add the total volume to a three-chamber cell counting chip and make note of the channels used (A, B, or C).
- 15 On the cell counter, select Fluorescence Cell Counting -> Cell Lines -> Advanced -> Protocol -> and then choose "Nuclei" from the list of protocols.

Then go to Settings and choose the 3 Channel option. Select all of the channels you will be using.

Then hit "Count" and then "Start Count."
- 16 When the cell count is complete, you will get a reading in cells/mL. Convert this to cells/ μ L by dividing this number by 1,000.
- 17 Save approximately 40,000 nuclei from the original sample to use as our unpurified population

for sequencing. (If the concentration was 1,000 cells/uL, then save 40 uL).

- 18** If there are channels left on the cell counter chip, mark the used channels on the back and place it back in the drawer for future use.

Immunoprecipitation

- 19** Vortex stock of Dynabeads Protein G and then add 5 uL to a 1.5 mL eppendorf. Place it on the mag rack for a minute and discard the supernatant remaining at the bottom of the tube without disturbing the beads.
- 20** Remove the tube from the magnet and resuspend the beads with 500 uL of Citrate Phosphate Buffer. Place the tube back on the mag rack for a minute and discard the supernatant. Repeat this step once more.
- 21** Add 100 uL of 1X PBS and 2 uL anti-GFP to the washed Dynabeads. Add this to the rotator in the 4°C fridge and let it incubate for 1 hour.
- 22** After incubating, place the tube on the magnet for 2 minutes and discard the supernatant without disturbing the beads.
- 23** Wash with 100 uL 1X PBS Buffer, enough to submerge the beads while they are on the rack. Discard the supernatant.
- 24** Add the total volume of nuclei sample to the beads. Incubate on the rotator in the 4°C fridge for 1 hour.
- 25** After incubation, place the samples on the mag rack for 2 minutes. Collect the supernatant in appropriately labeled eppendorf tubes. This will be our supernatant population for sequencing.



- 26 Wash the bead tubes with enough 1X PBS to submerge the beads (don't resuspend, gently) and then discard the supernatant.
- 27 Remove the tubes from the rack and then resuspend the beads in 50 uL 1X PBS.

Cell Count

- 28 Take another cell count of the supernatant and bead-bound samples. Reverse pipette to record their volumes.
- 29 Dilute all of the samples (unpurified, supernatant, and bead-bound) to 1,000 cells/uL. The bead-bound population likely won't need dilution.
- 30 Follow the 10X Chromium manual for the rest of the protocol to proceed with library prep and sequencing.