

Aug 22, 2024

## Human Dorsal Root Ganglion spatial ATAC-seq protocol

DOI

**[dx.doi.org/10.17504/protocols.io.36wgqn5n5gk5/v1](https://dx.doi.org/10.17504/protocols.io.36wgqn5n5gk5/v1)**

Úrzula Franco-Enzástiga<sup>1</sup>, Theodore Price<sup>1</sup>

<sup>1</sup>University of Texas at Dallas



Úrzula Franco-Enzástiga

University of Texas at Dallas

OPEN  ACCESS



DOI: **[dx.doi.org/10.17504/protocols.io.36wgqn5n5gk5/v1](https://dx.doi.org/10.17504/protocols.io.36wgqn5n5gk5/v1)**

**Protocol Citation:** Úrzula Franco-Enzástiga, Theodore Price 2024. Human Dorsal Root Ganglion spatial ATAC-seq protocol. protocols.io **<https://dx.doi.org/10.17504/protocols.io.36wgqn5n5gk5/v1>**

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** April 11, 2024

**Last Modified:** August 22, 2024

**Protocol Integer ID:** 98090

**Keywords:** ATAC-seq, spatial ATAC-seq, human DRG, chromatin, epigenomics

**Funders Acknowledgement:**

**National Institutes of Health**

**Grant ID:** U19NS130608

**National Institutes of Health**

**Grant ID:** R01NS111929



## Abstract

In this protocol, we describe how to perform spatial ATAC-seq on fresh-frozen human dorsal root ganglion.

## Materials

### *Key materials*

- OCT - Fisher HealthCare (Cat. No. 23-730-571)
- SuperFrost Plus charged slides - Fisher Scientific (Cat. No. 12-550-15)

### *Service*

- Spatial ATAC-seq technology product (AXO-0310); Microfluidic barcoding system.

### *Additional materials and equipment required*

- Dry ice
- Stainless steel base mold - Ted Pella (Cat. No. 27276-6)
- Cryostat Leica CM1950

## Before start

Required PPE: All work must be done wearing appropriate PPE including lab coat and gloves. Proper safety precautions should be taken into account while working with and disposing human tissue.

## Tissue harvesting and storage

- 1 Surgically excised lumbar L4 or L5 human dorsal root ganglia (hDRGs) are obtained from organ donors at Southwest Transplant Alliance (STA) at 2 h post cross-clamp.
- 2 Right after dissection, tissue is fresh frozen on powdered dry ice and stored at -80 °C until processing.

## Tissue sectioning

- 3 Frozen hDRGs are embedded in OCT (Optimal Cutting Temperature embedding medium) in a cryomold by adding small volumes of OCT over dry ice to avoid thawing.
- 4 Tissues are cryostat sectioned at 10 µm onto SuperFrost Plus charged slides.
- 5 Slides are briefly warmed to allow the sections to adhere to the slides but are immediately returned to the -20 °C cryostat chamber until sectioning is complete.
- 6 The slides are removed from the cryostat and stored in a -80 °C freezer until shipment to AtlasXomics for further processing. Samples are shipped on dry ice.

## Transposition for Spatial ATAC-seq (AtlasXomics: AXO-0310)

- 7 Spatial ATAC-seq assay is performed at AtlasXomics (Product code: AXO-0310; <https://www.atlasxomics.com/howitworks>).
- 7.1 In brief, tissue is fixed with 0.2% formaldehyde in PBS for 5 min and quenched with 1.25 M glycine for 5 min.
- 7.2 After fixation, tissue is washed with PBS and nuclease-free water and left to dry.
- 7.3 Tissue sections are permeabilized for 15 min with a buffer containing 10 mM Tris-HCl pH 7.4, 10 mM NaCl, 3 mM MgCl<sub>2</sub>, 0.01% Tween-20, 0.1% NP40, 0.001% digitonin and 1%BSA, and washed with a solution containing 10 mM Tris-HCl pH 7.4, 10 mM NaCl, 3 mM MgCl<sub>2</sub>, 0.1% Tween-20 and 1%BSA.



- 7.4 Tn5 transposition on the fixed hDRG sections is carried out at 37°C for 30 min. Transposition reaction buffer contains tagmentation buffer, 0.1% Tween-20, 0.01% digitonin, PBS and Tn5). Transposition reaction is stopped with 40 µM EDTA for 5 min.
- 8 Spatial barcoding is performed with microfluidics system (Product code: AXO-0310).
- 8.1 To correlate spatially barcoded accessible chromatin with tissue morphology, images of hDRG sections are acquired (Keyence BZ-X810 at 10X magnification).
- 8.2 The tissue is lysed using a reverse cross-linking solution containing 200 nM NaCl, 50 mM Tris-HCl pH 8.0, 1 mM EDTA, 1% SDS and 0.4 mg/mL proteinase K at 58°C for 2 h to release DNA fragments, which are then purified and amplified through PCR for library preparation using standard Next Generation Sequencing methods.
- 8.3 The PCR product is mixed with SPRI magnetic beads and incubated for 10 min. A magnetic stand is used to collect the beads, thereby purifying the final PCR product, which is eluted in nuclease-free water.
- 8.4 150X150 paired-end sequencing is performed on a NextSeq 2000 with 15% PhiX. A total of 300 million reads per sample are targeted.

## Spatial ATAC-seq analysis (AtlasXomics: AXO-0310)

- 9 Processing of spatial ATAC-seq is performed by AtlasXomics.
- 9.1 Use *Chromap* to trim sequencing adapters, align sequenced reads to the human reference genome GRCh38/hg38 and remove duplicates from fastq files.
- 9.2 Fragments files containing coordinates of sequenced DNA mapping to barcodes are generated for downstream analysis.
- 9.3 Use *AtlasXbrowser* to designate pixels on tissue based on microscopy images and to create Seurat-compatible metadata files.
- 9.4 Use *ArchR* to process fragments files and remove pixels not on tissue.
- 9.5 At this point of the processing, the .tsv fragments file contains the filtered aligned reads which can be used to proceed with peak calling. We have provided a tissue\_positions\_list.csv file containing a table with rows corresponding to the 2500 spots on the array, high- and low-

resolution images of the tissue, a `scale_factors.json` file, and experiment metadata. Together, these files can help for further processing of spatial ATAC-seq analysis using *ArchR*.

A brief description of some of the provided files is shown below.

### ***tissue\_positions\_list.csv***

Columns on the file `tissue_positions_list.csv` correspond to:

- Column 1: barcode
- Column 2: `in_tissue` binary status indicating if the spot falls on a tissue spot
- Column 3: array row
- Column 4: array column
- Column 5: the row pixel coordinate of the center of the spot in the full resolution image
- Column 6: the column pixel coordinate of the center of the spot in the full resolution image

### ***scale\_factors.json***

Parameters in the `scale_factors.json` file correspond to the relative scales of the supplied image and the array barcode features. Four aspects are provided:

- `fiducial_diameter_fullres`: Number of pixels spanning the diameter of a fiducial spot in the original, high-resolution image
- `spot_diameter_fullres`: Number of pixels spanning the diameter of a spot in the original, high-resolution image
- `tissue_highres_scalef`: Scaling factor to convert pixel positions in the original, full-resolution image to pixel positions in the provided high resolution image (.png file)
- `tissue_lowres_scalef`: Scaling factor to convert pixel positions in the original, full-resolution image to pixel positions in the provided low-resolution image (.png file)

## Protocol references

Deng Y, *et al.* (2022) Spatial profiling of chromatin accessibility in mouse and human tissues. *Nature* 609(7926):375-383.

Granja JM, *et al.* (2021) ArchR is a scalable software package for integrative single-cell chromatin accessibility analysis. *Nat Genet* 53(3):403-411.

Zhang H, *et al.* (2021) Fast alignment and preprocessing of chromatin profiles with Chromap. *Nat Commun* 12(1):6566.