



May 29, 2020

NEBExpress Ni Spin Column Reaction Protocol (NEB #S1427)

New England Biolabs¹¹New England Biolabs

1 Works for me

This protocol may be deleted by the owner

Alicia Bielik

ABSTRACT

Ni resin can be used for the purification of His-tagged fusion proteins under native or denaturing conditions

- The binding capacity of NEBExpress[®] Ni Spin Columns is ≥ 1 mg per column. The binding capacity can vary depending on the size of the target protein, binding conditions and the accessibility of the His-tag. An exact protocol may need to be optimized by the user.
- It is recommended to estimate the expression level of the His-tagged protein of interest by first running a sample on an SDS-PAGE gel.

EXTERNAL LINK

<https://www.neb.com/protocols/2019/08/28/nebexpress-ni-spin-column-reaction-protocol-neb-s1427>

GUIDELINES

Chemical Compatibility

Reagent	Tolerance (up to)
EDTA	10 mM*, 100 mM**
DTT	5 mM
b-mercaptoethanol	20 mM
TCEP	5 mM
Triton X-100	2 %
Tween 20	2 %
NP-40	2 %
Cholate	2 %
CHAPS	1 %
Tris-HCl, HEPES, MOPS	100 mM
Urea	8 M
Guanidine-HCl	6 M

* If reagents contain 10 mM EDTA, do not mix the sample and the resin for more than 24 hours before washing and eluting.

** If reagents contain 100 mM EDTA, do not mix the sample and the resin for more than 2 hours before washing and eluting.

MATERIALS

NAME	CATALOG #	VENDOR
NEBExpressNi Spin Columns – 25 columns	S1427L	New England Biolabs
NEBExpressNi Spin Columns – 10 columns	S1427S	New England Biolabs

MATERIALS TEXT

Sodium Phosphate
NaCl
Imidazole, pH 7.4
H₂O
Centrifuge
2ml microcentrifuge tubes

SAFETY WARNINGS

Please refer to the Safety Data Sheets (SDS) for health and environmental hazards.

BEFORE STARTING

Ni resin can be used for the purification of His-tagged fusion proteins under native or denaturing conditions

- The binding capacity of NEBExpress® Ni Spin Columns is ≥ 1 mg per column. The binding capacity can vary depending on the size of the target protein, binding conditions and the accessibility of the His-tag. An exact protocol may need to be optimized by the user.
- It is recommended to estimate the expression level of the His-tagged protein of interest by first running a sample on an SDS-PAGE gel.

Buffer Preparation for Ni Spin Columns

- 1 Prepare buffers for Ni Spin Columns as follows:

	Lysis/Binding Buffer:	Wash Buffer:	Elution Buffer:
	20 mM sodium phosphate, 300 mM NaCl	20 mM sodium phosphate, 300 mM NaCl, 5 mM Imidazole	20 mM sodium phosphate, 300 mM NaCl, 500 mM Imidazole
2X IMAC Buffer	7.5 ml	5.0 ml	2.5 ml
2M Imidazole	-	0.025 ml	1.25 ml
H2O	7.5 ml	5 ml	1.25 ml
Total	15.0 ml	10.0 ml	5.0 ml



Notes:

1. When the recommended protocol is followed, each isolation requires the following volumes of buffer: **0.75 ml Lysis/Binding Buffer** , **0.75 ml Wash Buffer** and **0.4 ml Elution Buffer** . An excess volume of each concentrated buffer is provided for preparation of cell lysates or to optimize the concentration of imidazole in the washes.
2. Crude lysate should be prepared with a lysis buffer without imidazole. To further minimize contaminants in the eluate, the concentration of imidazole in the wash buffer can be increased to 10 mM (≥ 10 mM can reduce the isolated yield but may result in increased purity). imidazole concentration should be determined empirically.
3. Refer to the Chemical Compatibility table prior to including other additives.

- 2 To prepare buffers for Ni Spin Columns under **Native** Conditions: bring all three buffers to a final pH of 7.4
To prepare buffers for Ni Spin Columns under **Denaturing** Conditions: bring all three buffers to a final concentration of **[M]8 Molarity (M) Urea** or **[M]6 Molarity (M) Guanidine** .

Sample Preparation

- 3 Harvest cells by centrifugation at **4000 x g 00:15:00** , store the pellet at **-20 °C** or process immediately.



Note: it is recommended to pre-weigh the vessel prior to addition of cell suspension, in order to determine the mass of cell pellet used.

- 4 Resuspend cell pellet in Lysis Buffer and lyse using method of choice (use approximately **5 ml lysis buffer** per **1 g cell paste**).



Note: Cells can be lysed by standard methods including sonication, repeated freeze-thaw cycles, French press, etc. Other commercially available lysis reagents can also be used, following manufacturer's instructions. It is recommended that imidazole be omitted from any lysis buffer.

- 5 Centrifuge sample at **12000 x g 00:15:00** to pellet cellular debris. Remove the clarified protein lysate supernatant and transfer to a new microcentrifuge tube on ice, **retain a 2 µl aliquot of the clarified lysate** for SDS-PAGE analysis.



Note: A standard isolation typically employs 0.5 ml of clarified lysate.

Column Preparation

- 6 Remove the bottom tab of the column by twisting, loosen the top cap and place the column in the collection tube provided.
- 7 Centrifuge column at **800 x g 00:01:00** to remove the storage buffer, discard the buffer.
- 8 Add **250 µl Lysis/Binding buffer** to the column.
- 9 Centrifuge column at **800 x g 00:01:00** , discard the Lysis/Binding buffer.
- 10 Place the column in a new 2 ml microcentrifuge tube.

Lysate Binding

- 11 Add **up to 500 µl protein sample lysate** to the column.



Note: If sample volume is greater than 500 µl multiple applications can be performed; collect the flow through in separate microcentrifuge tubes.

- 12 Tap the column to mix the lysate with the resin and allow binding for **00:02:00** .



Note: Binding of some His-tagged proteins can be increased with longer mixing times. Cap the column and seal the bottom using the plug provided. Mix end-over-end at **4 °C** for desired time (typically **00:05:00** - **00:15:00**). Prolonged mixing may result in more non-specific binding.

13 Centrifuge column at **800 x g 00:01:00**, reserve flow through.

14 Place the column in a new 2 ml microcentrifuge tube.

Column Wash

15 Add **250 µl wash buffer** to the column and centrifuge at **800 x g 00:01:00**.

16 Repeat wash step twice, collect each wash in a separate 2 ml microcentrifuge tube: Add **250 µl wash buffer** to the column and centrifuge at **800 x g 00:01:00**. (1/2)

17 Repeat wash step twice, collect each wash in a separate 2 ml microcentrifuge tube: Add **250 µl wash buffer** to the column and centrifuge at **800 x g 00:01:00**. (2/2)

Protein Elution

18 Place column in a new 2 ml microcentrifuge tube.

19 Add **200 µl elution buffer** to the column. Mix the resin with the elution buffer thoroughly.



Note: Elution volume can be reduced to **100 µl** if a more concentrated protein sample is desired, total target protein yield may be lower.

20 Centrifuge at **800 x g 00:01:00**, save eluted sample.

21 Place column in a new 2 ml microcentrifuge tube and repeat elution step. **Typically, >90% of the bound protein is eluted following the second elution.**

22 Analyze the clarified cell lysate (load), flow through, washes and eluates by SDS-PAGE.