



MAR 14, 2024

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.3b4l4qmoovo5/v1

Protocol Citation: Xiaobo Mao, Ramhari Kumbhar, Hanseok Ko, Valina L. Dawson, Ted Dawson 2024. Grip Strength. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.3b4l4qmoovo5/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Mar 14, 2024

Grip Strength

Xiaobo Mao^{1,2,3,4}, Ramhari Kumbhar^{1,2,3,4}, Hanseok Ko^{1,2,4}, Valina L. Dawson^{1,2,4,3,5,6}, Ted Dawson^{1,2,4,3,5,6,7}

¹Neuroregeneration and Stem Cell Programs, Institute for Cell Engineering, Johns Hopkins University School of Medicine, Baltimore, MD 21205, USA;

²Department of Neurology, Johns Hopkins University School of Medicine, Baltimore, MD 21205, USA;

³Adrienne Helis Malvin Medical Research Foundation, New Orleans, LA 70130-2685, USA;

⁴Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815, USA;

⁵Department of Physiology, Johns Hopkins University School of Medicine, Baltimore, MD 21205 USA;

⁶Solomon H. Snyder Department of Neuroscience, Johns Hopkins University School of Medicine, Baltimore, MD 21205, USA;

⁷Institute for NanoBioTechnology, Johns Hopkins University, Baltimore, MD, USA

ASAP Collaborative Research Network

Kapliitt Protocols



Eileen Ruth Torres
Weill Cornell Medicine

ABSTRACT

Use to measure rodent grip strength

MATERIALS

Bioseb Rodent Grip Strength BIO-GS4

Last Modified: Mar 14, 2024

PROTOCOL integer ID: 96717

Keywords: ASAPCRN, mouse model, behavior

Funders Acknowledgement:
Aligning Science for Parkinson's
Grant ID: 020608

- 1 Set up grip strength meter on table.
- 2 Hold the mouse gently by the tail, and place either its forepaws or all paws on the grip strength grid.
- 3 Gently pull the mouse back until it releases its grasp.
- 4 Record the maximum holding force displayed on the grip strength meter.