



JAN 30, 2024

🌐 Immunohistochemistry free-floating rat brain cryosections

📁 In 1 collection

mariangela.massarocenere^{1,2,3}

¹Department of Experimental Neuroscience, Santa Lucia Foundation IRCCS, Rome, Italy;

²Department of Systems Medicine, University of Rome Tor Vergata, Rome, Italy;

³Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, United States



mariangela.massarocenere

ABSTRACT

Protocol for immunohistochemistry on rat brain cryosections

MATERIALS

Reagents:

- Phosphate buffer solution (PB)1X.
- TritonX-100
- Endogenous peroxidase blocking solution: PB 1x + 3% H₂O₂
- Blocking Buffer: PB 1X + 10% NGS
- 1st Ab : Diluted in1X PBS +5% NGS+ 0.3% TX
- 2nd Ab : Diluted in1X PBS +5% NGS+ 0.3% TX
- 3,30-diaminobenzidine (DAB)
- Xilen
- Ethanol : 100%, 95%, 70%, 50%

OPEN  ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.3b4l4q4zrvo5/v1

Protocol Citation: mariangela.massarocenere 2024.

Immunohistochemistry free-floating rat brain cryosections .

protocols.io

<https://dx.doi.org/10.17504/protocols.io.3b4l4q4zrvo5/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: Jan 25, 2024

1. Sections selection

- 1 Collect the cryosections needed for a caudo-rostral representation of each brain region (every fifth or sixth sections depending on the section's thickness, brain region, and animal species) into a 24-well-plate (3-4 sections per well)

1.1 Wash in PB1X 3 times x 5 min (500 ul) at RT

2. Inactivation of endogenous peroxidase

10m

- 2 Incubate sections with endogenous peroxidase blocking solution: PB1X + 3% H₂O₂ (500 µl/well) for at RT 10m

2.1 Wash in PB1X 2 times x 5 min (500 ul) at RT

3. Blocking

- 3 Incubate with blocking solution (500 µl/well) : PB1X + 10%NGS (serum from the same species as the host of the secondary antibody) + 0.3% TritonX-100, at RT 1h

4. Primary antibody incubation

- 4 Remove the blocking solution and incubate sections with PB1X+ 5%NGS+0.3% TritonX-100+ primary Ab (500 µl/well) for 24/72 h, depending on the Ab, at +4°C

4.1 Wash in PB1X 3 times x 5 min (500 ul) at RT

5. Secondary biotinylated antibody incubation

2h

5 Incubate sections with adequate secondary biotinylated Ab diluted in PB1X+ 0.3% TritonX-100, at RT

2h

5.1 Wash n PB1X 3 times x 5 min (500 ul) at RT

6. Extravidin-peroxidase reaction

1h 30m

6 Incubate with Extravidin-peroxidase buffer solution diluted in PB1X+0.3% TritonX-100, at RT

1h 30m

6.1 Wash n PB1X 3 times x 5 min (500 ul) at RT

7. Substrate preparation

7 Allow each DAB tablet to reach room temperature before use and Dissolve a DAB tablet in 15 mL of Tris-buffered saline, pH 7.6

7.1 Add 12 µL of fresh 30% hydrogen peroxide prior to use

7.2 Filter the solution through a 0.2 µm filter immediately before use

8. Development

3m

8 Add 500 µl/ well of substrate solution in the dark

3m

8.1 Remove and clean all material with bleach

8.2 Wash 2 times x 5 min in PB1x

9. Section Mounting

9 Mount sections on SuperFrost plus slides and let it dry overnight

10. Dehydration

1m

10 Dehydrate in consecutive steps of ethanol solutions (50%, 70%, 95%, 100%)

1m

10.1 Incubate slides in 2 times x 5 min in Xylene

11. Mount Coverslips

11 Coverslip slides with mounting medium, remove bubbles if any, and let dry