



2 ▼

Feb 23, 2022

Cross-linking of IgG to Protein A or G Beads (S1425/S1430) V.2

[New England Biolabs¹](#)¹New England Biolabs

1

dx.doi.org/10.17504/protocols.io.bhkwj4xe**New England Biolabs (NEB)**Tech. support phone: **+1(800)632-7799** email: **info@neb.com****New England Biolabs**
New England Biolabs

This protocol consists of an IgG purification step followed by covalent cross-linking of the IgG to the Protein A/G solid support. For IgG that has been previously purified, proceed directly to the cross-linking protocol.

DOI

dx.doi.org/10.17504/protocols.io.bhkwj4xe<https://www.neb.com/protocols/0001/01/01/cross-linking-of-igg-to-protein-a-or-g-beads>

New England Biolabs 2022. Cross-linking of IgG to Protein A or G Beads (S1425/S1430). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bhkwj4xe>

Julia Rossmanith



Protein A, Protein G, cross-linkage, cross linked, covalent cross-linking of the IgG to protein A/G, igG purification to beads

_____ protocol ,

Jun 17, 2020

Feb 23, 2022

Feb 23, 2022



New England Biolabs

New England Biolabs

38262

Overview

This protocol consists of an IgG purification step followed by covalent cross-linking of the IgG to the Protein A/G solid support. For IgG that has been previously purified, proceed directly to the cross-linking protocol.

MATERIALS




 [Protein A Magnetic Beads - 1 ml New England](#)

Biolabs Catalog #S1425S



 [Protein G Magnetic Beads - 1 ml New England](#)

Biolabs Catalog #S1430S

Materials Needed:

- Protein A ([NEB #S1425S](#)) or Protein G ([NEB #S1430S](#)) Magnetic Beads
- Elution Buffer: 0.1 M glycine-HCl (pH 2.5)
- Binding Buffer: 0.1 M NaPhosphate Buffer (pH 8.0)
- Dimethyl pimelidate dihydrochloride (Sigma, D-8388) dissolved at 25 mM in Cross-linking Buffer.
- Cross-linking Buffer: 0.2 M triethanolamine (pH 8.2)
- Blocking Buffer: 0.1 M ethanolamine (pH 8.2)
-  [10x Phosphate Buffered Saline Gibco,](#)
- **ThermoFisher Catalog #70011044**
-  [Tween 20 Sigma](#)
- **Aldrich Catalog #P9416-50ML**
-  [Sodium azide Sigma](#)
- **Aldrich Catalog #71289**
- Immunoglobulin in Binding Buffer
- 100 µl PBS, 0.1% Tween 20, 0.02% sodium azide

For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet).

The IgG Purification protocol is for the binding of  **20 µg purified IgG** or isolation of  **20 µg IgG from serum**.

IgG Immobilization

1 

The IgG immobilization protocol is for the binding of **20 µg purified IgG** or isolation of **20 µg IgG from serum**.

Vortex and thoroughly resuspend Protein A Magnetic Beads.

2 

Aliquot **100 µL bead suspension** to a sterile microcentrifuge tube.

3 


Add **500 µL 0.1 M NaPhosphate Buffer (pH 8.0)** and vortex to resuspend. Apply magnet for **00:00:30**, to pull beads to the side of the tube and remove supernatant. (Wash 1/2)

4 

Repeat wash: Add **500 µL 0.1 M NaPhosphate Buffer (pH 8.0)** and vortex to resuspend. Apply magnet for **00:00:30**, to pull beads to the side of the tube and remove supernatant. (Wash 2/2)

5 

Add to the beads **80 µL 0.1 M NaPhosphate Buffer (pH 8.0)**.

6 

Add **15 µL - 25 µL serum OR 20 µg purified IgG** in a maximum volume of **30 µL**.

7 

Mix thoroughly and incubate at **4 °C** with agitation for **00:30:00**.

8 Apply magnet and remove supernatant.

9 

Add **500 µL 0.1 M NaPhosphate Buffer (pH 8.0)** and vortex to resuspend. Apply magnet for **00:00:30**, to pull beads to the side of the tube and remove supernatant. (Wash 1/3)

10 

Add **500 µL 0.1 M NaPhosphate Buffer (pH 8.0)** and vortex to resuspend. Apply magnet for **00:00:30**, to pull beads to the side of the tube and remove supernatant. (Wash 2/3)

11 

Add **500 µL 0.1 M NaPhosphate Buffer (pH 8.0)** and vortex to resuspend. Apply magnet for **00:00:30**, to pull beads to the side of the tube and remove supernatant. (Wash 3/3)

IgG Cross-linking to Protein A/G Magnetic Beads

30m


12 

Add **1 mL Cross-linking Buffer (0.2 M triethanolamine, [pH 8.2])** to the Protein A/G immobilized antibody. (Wash 1/2)

At this point the purified IgG can be eluted from the beads or used directly for immunoprecipitation of target proteins. The purified IgG can also be cross-linked to the Protein A beads (see cross-linking protocol) to create a reusable immunoprecipitation bead which prevents the co-elution of antibody with target protein.

13 

Vortex to resuspend. (Wash 1/2)


14 Apply magnet for  **00:00:30** , to pull beads to the side of the tube and remove supernatant.
(Wash 1/2)

15 

Add  **1 mL Cross-linking Buffer (0.2 M triethanolamine, [pH 8.2])** to the Protein A/G immobilized antibody. (Wash 2/2)

16 

Vortex to resuspend. (Wash 2/2)

17 Apply magnet for  **00:00:30** , to pull beads to the side of the tube and remove supernatant.
(Wash 2/2)


18 

Resuspend in  **1 mL Cross-linking Buffer** containing
[M]25 Milimolar (mM) DMP (6.5 mg DMP/ml of buffer) .

19  

Mix thoroughly and incubate at  **Room temperature** for  **00:45:00** with agitation.
















20 Apply magnet for  **00:00:30** , to pull beads to the side of the tube and remove supernatant.


21 

Add  **1 mL Blocking Buffer (0.1 M ethanolamine, [pH 8.2])** .

22 

Vortex to resuspend.

- 23 Apply magnet for  **00:00:30** , to pull beads to the side of the tube and remove supernatant.
- 24  Add  **1 mL Blocking Buffer** .
- 25  Vortex to resuspend.
- 26  30m Incubate for  **00:30:00** at  **Room temperature** with agitation.
- 27 Apply magnet for  **00:00:30** , to pull beads to the side of the tube and remove supernatant.
- 28  Add  **1 mL PBS** . (Wash 1/3)
- 29  Vortex to resuspend. (Wash 1/3)
- 30 Apply magnet for  **00:00:30** , to pull beads to the side of the tube and remove supernatant. (Wash 1/3)
- 31  Add  **1 mL PBS** . (Wash 2/3)
- 32  Vortex to resuspend. (Wash 2/3)


33 Apply magnet for  **00:00:30** , to pull beads to the side of the tube and remove supernatant.
(Wash 2/3)

34 

Add  **1 mL PBS** . (Wash 3/3)

35 

Vortex to resuspend. (Wash 3/3)

36 Apply magnet for  **00:00:30** , to pull beads to the side of the tube and remove supernatant.
(Wash 3/3)

37 

Add  **1 mL Elution Buffer (0.1 M glycine-HCl [pH 2.5])** .

38 

Vortex to resuspend.

39 Apply magnet for  **00:00:30** , to pull beads to the side of the tube and remove supernatant.

This elutes bound antibody that is not cross-linked with DMP.

40 

Add  **1 mL PBS** . (1/2)

41 

Vortex to resuspend. (1/2)

42 Apply magnet for  **00:00:30** , to pull beads to the side of the tube and remove supernatant. ^{30s}


(1/2)

43 

Add  **1 mL PBS** . (2/2)

44 

Vortex to resuspend. (2/2)

45 Apply magnet for  **00:00:30** , to pull beads to the side of the tube and remove supernatant. ^{30s}
(2/2)

46 

Resuspend and store beads in  **100 µL PBS** , 0.1% Tween 20, 0.02% sodium azide.