



NOV 09, 2022



SHARE

WORKS FOR ME 1



Organoid Immunohistochemistry on Ventana Benchmark XT

Forked from [Organoid Immunohistochemistry on Ventana Benchmark XT](#)**This protocol is published without a DOI.**[gustavo.parfitt](#)¹, [kristen.whitney](#)¹¹Icahn School of Medicine at Mount Sinai

Ahfeldt Lab



gustavo.parfitt

COMMENTS 0

ABSTRACT

Organoid Immunohistochemistry on Ventana Benchmark XT

PROTOCOL CITATION

[gustavo.parfitt](#), [kristen.whitney](#) 2022. Organoid Immunohistochemistry on Ventana Benchmark XT.
protocols.io
<https://protocols.io/view/organoid-immunohistochemistry-on-ventana-benchmark-ci33ugqn>



FUNDERS ACKNOWLEDGEMENT

Kristen Whitney

FORK NOTE

FORK FROM

[Forked from Organoid Immunohistochemistry on Ventana Benchmark XT, kristen.whitney](#)

LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Nov 09, 2022

LAST MODIFIED

Nov 09, 2022

PROTOCOL INTEGER ID

72539

Immunohistochemistry of organoid sections

- 1 Fix organoids in 10% NBF for 15-30 minutes.
- 2 Prepare paraffin-embedded blocks on tissue processor.
- 3 Cut 5µm sections from FFPE (formalin-fixed paraffin-embedded) blocks using microtome.
- 4 Mount sections on positively charged slides and bake overnight at 70° C.
- 5 Prepare primary antibody dilution kit (depending on antibody concentration) for Ventana Benchmark XT automatic staining platform (Roche Diagnostics, Indianapolis, IN) according to the manufacturer's protocol.
- 6 Load slides onto Ventana and run through alcohol, antigen retrieval, primary antibody incubation and visualization steps.
- 7 Coverslip slides and let dry overnight.

Digital Histopathologic Analysis in QuPath

- 8 Scan slides at 40x on Aperio VERSA 8 (Leica Biosystems, Wetzlar Germany) digital slide scanner.
- 9 Manually segment organoids to generate annotations in QuPath (version 0.2.3, <https://QuPath.github.io/>).
- 10 Perform color deconvolution and determine appropriate settings for positive-pixel-based analysis.