



Aug 05, 2020

# De-salting of tryptic digested peptides

avinash.kale<sup>1</sup><sup>1</sup>domnic colvin

1

Works for me

[dx.doi.org/10.17504/protocols.io.bjdqki5w](https://dx.doi.org/10.17504/protocols.io.bjdqki5w)

UM-DAE Centre for Excellence in Basic Sciences



avinash.kale

DOI

[dx.doi.org/10.17504/protocols.io.bjdqki5w](https://dx.doi.org/10.17504/protocols.io.bjdqki5w)

DOCUMENT CITATION

avinash.kale 2020. De-salting of tryptic digested peptides. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bjdqki5w>

LICENSE

This is an open access document distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 05, 2020

LAST MODIFIED

Aug 05, 2020

DOCUMENT INTEGER ID

40080

## De-salting of tryptic digested peptides

The desalting of tryptic digested samples were optimized and performed using octadecyl carbon chain (C18) pipette tips (ZipTip® Pipette Tips, Millipore). Briefly, for the C18 tips, the desalting pipette tip was pre-equilibrated using a solution of 50% Acetonitrile (SRL) in LC-MS grade water and then equilibrated twice in 0.1% formic acid in LC-MS water. After carefully and gently pipetting the entire sample up and down for at least 10 times, the membrane was washed with 0.1% formic acid in LC-MS water. The peptides were eluted in 70% acetonitrile in LC-MS grade water and vacuum dried and stored for further analysis(26).