



Sep 06, 2022

In vitro assembling of RNP for nucleofection of hPSCs

In 1 collection

Hanqin Li¹, Oriol Busquets², Steven Poser², Dirk Hockemeyer¹, Frank Soldner²

¹University of California, Berkeley; ²Albert Einstein College of Medicine

1 Works for me Share

dx.doi.org/10.17504/protocols.io.b4pfqvjn



ABSTRACT

This protocol describes the procedure for the in vitro assembly of ribonucleoprotein (RNP) which can be delivered into human pluripotent stem cells (hPSCs) using nucleofection.

General notes

1. Throughout this protocol, the term hPSC is used to collectively refer to both hiPSCs and hESCs. All described procedures have been tested and work equally well for hiPSCs and hESCs.

DOI

dx.doi.org/10.17504/protocols.io.b4pfqvjn

PROTOCOL CITATION

Hanqin Li, Oriol Busquets, Steven Poser, Dirk Hockemeyer, Frank Soldner 2022. In vitro assembling of RNP for nucleofection of hPSCs. **protocols.io** https://dx.doi.org/10.17504/protocols.io.b4pfqvjn

FUNDERS ACKNOWLEDGEMENT

Aligning Science Across Parkinson's

Grant ID: ASAP-000486



COLLECTIONS (i)

Nucleofection (Amaxa) and electroporation (Biorad) of hPSCs

KEYWORDS

ASAPCRN

LICENSE

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Feb 03, 2022

LAST MODIFIED

Sep 06, 2022

PROTOCOL INTEGER ID

57799

PARENT PROTOCOLS

Part of collection

Nucleofection (Amaxa) and electroporation (Biorad) of hPSCs

MATERIALS TEXT

	Item	Vendor	Catalog #
	Cas9, purified protein, 40 μM	Macrolab, QB3 UC	
		Berkeley	
	Synthetic pegRNAs	IDT or Synthego	
	Synthetic sgRNAs	Synthego	

- 1 Thaw Cas9 or PE2 protein § On ice . We don't usually reuse leftover protein, since freezing and thawing compromises protein activity.
- 2 In each nucleofection reaction, prepare 10 μl RNP

3 For regular CRISPR/Cas9 editing, use:

5 μl, autoclaved H2O

2 μl, 40 μM purified Cas9 protein

3 μl, 100 μM synthetic sgRNA

4 For prime editing PE2 strategy, use:

7-x µl, autoclaved H2O

x μl (90 pmol), purified PE2 protein

3 μl, 100 μM synthetic pegRNA

5 For prime editing PE3 strategy, use:

7-x µl, autoclaved H2O

x μl (90 pmol), purified PE2 protein

2 μl, 100 μM synthetic pegRNA

1 μl, 100 μM synthetic ngRNA

- 6 Pipet the proper amount of each component in the order showing above into a microcentrifuge tube. Mix using a P10/20 tip.
- 7 Incubate at & Room temperature for © 00:10:00 to assemble RNPs. Once assembled, RNPs are stable at room temperature for 30 min, or at least 2h on ice.
- 8 If HDR template is desired, mix 100 pmol of template ssODN or 1 µg targeting vector with assembled RNPs right before mixing with cells.