

Sep 12, 2022

DNA Extraction

klaus.stark¹, Kira Julia Stanzick¹¹Department of Genetic Epidemiology, University of Regensburg

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.5qpvojbw7l4o/v1

Kira Julia Stanzick

ABSTRACT

DNA extraction from frozen blood clots is challenging. Here we applied QIAGEN Clotspin Baskets and the Gentra Puregene Blood Kit for DNA extraction from 5.5 ml whole blood without anticoagulating additives.

DOI

dx.doi.org/10.17504/protocols.io.5qpvojbw7l4o/v1

PROTOCOL CITATION

klaus.stark, Kira Julia Stanzick 2022. DNA Extraction. **protocols.io**
<https://protocols.io/view/dna-extraction-cfiutkew>



LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 19, 2022

LAST MODIFIED


















Sep 12, 2022

PROTOCOL INTEGER ID



68916







Blood clot preparation and red blood cell lysis 40m 19s



- 1 Frozen blood clots in 15 ml tubes were transferred from -20°C to a warming cabinet ^{10m} 55 °C for 00:10:00 and thereafter immediately placed On ice

- 2 The tube was inverted to loosen the clot. The blood clot was completely poured with  **5 mL**  **RBC Lysis Solution Qiagen** into the Clotspin Basket placed on a 50 ml tube
- 3 To disperse the clot, centrifugate  **2000 x g, 00:05:00** 5m
- 4 The remaining clot material from the Clotspin basket was transferred through the basket to ^{10m} the filtrate with  **00:10:00**  **RBC Lysis Solution Qiagen** and the basket was discarded.
- 5 To completely disperse the clotted material, the filtrate was **vortexed** vigorously for 5m 3s
 **00:00:03** and placed for  **00:05:00** at  **Room temperature** on a circulating **shaker** (250 1/min).
- 6 The tubes were again **vortexed** vigorously for  **00:00:03** and centrifugated 5m 3s
 **2000 x g, 00:05:00** at 2000 x g for 5 min
- 7 The supernatant was carefully discarded, taking care that the pellet remains in the tube. If no pellet was visible, about  **0.5 mL** of the supernatant was kept in the tube
- 8 The tube was **vortexed** rigorously for  **00:00:10** and additional  **5 mL** 5m 13s
 **RBC Lysis Solution Qiagen** was added to the pellet, followed by **vortexing** for
 **00:00:03** and incubation on a circulating shaker (250 1/min) for  **00:05:00** at
 **Room temperature**






White blood cell lysis 1d 0h 5m 30s




- 9 Centrifugation  **2000 x g, 00:05:00** to pellet the DNA-containing white blood cells 5m
- 10 Discard supernatant carefully, leaving about 0.2 ml of residual liquid. Vortex rigorously for 10s
 **00:00:10**

11 Adding of  **5 mL**  **Cell Lysis Solution Qiagen** and  **25 µL** of  **Proteinase K Qiagen**  **20 mg/mL** as followed by rigorously **vortexing** for  **00:00:10** 10s



12 For complete lysis of DNA containing cells, incubate  **Overnight** at  **55 °C** 10s


Protein precipitation 17m 20s

13 samples were cooled down  **On ice** for  **00:05:00** and  **1.7 mL**  **Protein Precipitation Solution Qiagen** was added, followed by rigorous **vortexing** for  **00:00:20** 5m 20s

14 After centrifugation  **2000 x g, 00:10:00** , incubate the samples for  **00:02:00**  **On ice** . 12m

DNA precipitation 3m

15 For precipitation of DNA, the supernatant was **carefully** transferred in a 50 ml tube, containing  **5 mL**  **Isopropanol Contributed by users**

16 The samples were mixed by gently **inverting the tube 50 times** and centrifuged  **2000 x g, 00:03:00** 3m
























17 The supernatant was **carefully** discarded, and the tube was drained on a clean piece of absorbent paper without losing the pellet. The DNA pellet was inspected visually and with two following options
Step 17 includes a Step case.
big red to brown pellet
white pellet

DNA washing 4h 23m 20s








step case

big red to brown pellet

Additional washing step to remove remaining erythrocytes

- 18  **2 mL** of  **Cell Lysis Solution Qiagen** and  **10 µL** of  **Proteinase K (20 mg/ml) Contributed by users** were added, followed by incubation at  **55 °C** for  **02:00:00** in a warming cabinet **or** at  **Room temperature**  **Overnight** on a circulating shaker (**250 1/min**) 4h
- 19 samples were cooled down  **On ice** for  **00:05:00** and  **1 mL**  **Protein Precipitation Solution Qiagen** was added, followed by rigorous **vortexing** for  **00:00:20** 5m 20s
- 20 After centrifugation  **2000 x g, 00:10:00**, the samples were incubated for  **00:02:00**  **On ice** 12m
- 21 DNA precipitation was done with  **3 mL**  **Isopropanol Contributed by users**. The samples were mixed by gently **inverting the tube 50 times** and centrifugated  **2000 x g, 00:03:00** 3m
- 22 The supernatant was carefully discarded. The tube was drained on a clean piece of absorbent paper, taking care that the pellet remained in the tube
- 23  **5 mL** of  **70 % (v/v)**  **ethanol Contributed by users** were added immediately and the tubes were inverted until the pellet was detached
- 24 After centrifugation  **2000 x g, 00:03:00** at 2000 x g for 3 min, the supernatant was carefully discarded, and the DNA pellet was air dried at room temperature for 10 min or until the pellet got glassy 3m

DNA hydration 1h 2m

- 25 Addition of  **0.5 mL**  **DNA Hydration Solution Qiagen** for a **large** pellet **or**  **0.3 mL** for a **smaller** pellet was followed by incubation at  **65 °C** in a warming cabinet for  **01:00:00**. To fully dissolve the DNA, the sample were put on a **shaker**  **Overnight** at  **Room temperature** 1h 2m

26 Samples were centrifugated briefly and the solved DNA was transferred to a 2 ml cup and stored at -20°C .