



Dec 05, 2020

Rapid column Ag test

Anne Berking¹¹Berking Biotechnology

1

Works for me

dx.doi.org/10.17504/protocols.io.bqdjms4n

Coronavirus Method Development Community

XPRIIZE Rapid Covid Testing



Berking Biotechnology

ABSTRACT

The "Rapid Ag test based on Alpaca Nanobodies" is a test developed with single domain antibodies. This type of antibody is capable of effectively recognize the Spike protein on the surface of SARS-CoV2 leading a signal that is 10 times stronger than conventional antibodies.

It takes only 16 minutes from sample to result. The procedures for using this test were designed to ensure an easy use and interpretation, making it a quick and simple test.

DOI

dx.doi.org/10.17504/protocols.io.bqdjms4n

PROTOCOL CITATION

Anne Berking 2020. Rapid column Ag test. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bqdjms4n>



LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Dec 03, 2020

LAST MODIFIED

Dec 05, 2020

PROTOCOL INTEGER ID

45195

A	B	C
Equipment	Supplier	Catalog #
NucleoVac 24 Vacuum Manifold	Macherey-Nagel	740299
Manifold for 24 samples	KNF	046209/049540

A
Kit Content
Sample collection tube
Collection Buffer
Pasteur pipette
Binding buffer
Detection Column
Washing Buffer 1
Washing Buffer 2
Developing solution

SAFETY WARNINGS

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

Please, wear the correct personal protection equipment (i.e. gloves) and wash your hands immediately after removing the gloves.

When working with human saliva and other human bodily fluids, there may pathogens present.

DISCLAIMER:

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

ABSTRACT

The "Rapid Ag test based on Alpaca Nanobodies" is a test developed with single domain antibodies. This type of antibody is capable of effectively recognize the Spike protein on the surface of SARS-CoV2 leading a signal that is 10 times stronger than conventional antibodies.

It takes only 16 minutes from sample to result. The procedures for using this test were designed to ensure an easy use and interpretation, making it a quick and simple test.

BEFORE STARTING

Please follow this protocol in order to use "The Alpaca Nanobodies Rapid Ag test" correctly and get the accurate results

Sample Collection		5m 30s
1	Add 3 mL of collection buffer to sample collection tube and collect saliva in sample collection tube (3-5mL) wait 5 minutes	5m 30s
Sample Analyse		6m 50s
2	Mix binding buffer tube, load 1 drop of binding buffer into Column connected to a vacuum device	30s
3	Load sample into column, turn on vacuum device, and let sample go through column wait 5 minutes	5m 30s
4	Add 12mL of washing buffer 1, turn on vacuum device, and let washing buffer go through column	15s
5	Add 6 mL of washing buffer 2, turn on vacuum device, and let washing buffer go through column	15s
6	Add 12 mL of washing buffer 1, turn on vacuum device, and let washing buffer go through column	15s
7	Add 1 drop of developing solution	5s
Analysis of the Result		5m
8	Wait 5 min to visualize reaction blue color = positive, white color = negative	5m