



Jul 25, 2022

Tube Radius Calculation

Liv Jensen¹¹Hurley Lab

1 Works for me

 Sharedx.doi.org/10.17504/protocols.io.x54v9y3rqg3e/v1 Liv Jensen

DISCLAIMER

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

ABSTRACT

This protocol details Tube Radius Calculation.

ATTACHMENTS

[416-898.pdf](#)

DOI

dx.doi.org/10.17504/protocols.io.x54v9y3rqg3e/v1

PROTOCOL CITATION

Liv Jensen 2022. Tube Radius Calculation. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.x54v9y3rqg3e/v1>

KEYWORDS

Tube Radius Calculation, Fluorescence ratio calibration, ASAPCRN

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

May 03, 2022

LAST MODIFIED

Jul 25, 2022

OWNERSHIP HISTORY

May 03, 2022  maria.s

May 25, 2022  Liv Jensen

PROTOCOL INTEGER ID

61846

MATERIALS TEXT

Materials:

- Image stack of multichannel confocal fluorescence images
- Pressure data over time for GUV aspiration pipette
- Force data over time for bead held in optical trap
- FIJI
- Image segmentation and quantification script (Python)

DISCLAIMER:

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Fluorescence ratio calibration:

1 Measure the radius of the GUV and of the aspiration pipette from imaging data.

2 Calculate radius from:

$$R = \frac{F}{4\pi \left(\frac{\Delta P * r_{pip}}{2(1 - \frac{r_{pip}}{r_{GUV}})} \right)}$$

3 In FIJI, create new image stacks from two rectangular selections: one containing a section of the membrane tube, and one containing the approximately horizontal section of the guv.

4 Save these new stacks as .tif files and paste paths into segmentation and quantification script.

5 Use the output of the segmentation and quantification script to calculate the ratio of fluorescence intensity of the membrane tube and the GUV surface (I_{tub} / I_{GUV}).

6 Plot tube radius as calculated in step 2 vs. the fluorescence ratio calculated in step 5. Perform a linear fit of the data, recording the slope, K_{tub} , for use in subsequent experiments.

Calculating tube radius with fluorescence images and K_{tub} :

7 In FIJI, create new image stacks from two rectangular selections: one containing a section of the membrane tube, and one containing the approximately horizontal section of the guv.

8 Save these new stacks as .tif files and paste paths into segmentation and quantification script.

9 Input K_{tub} as a parameter in the quantification function to normalize fluorescence ratio to tube radius.

