

VERSION 2  
SEP 27, 2023

OPEN ACCESS



**DOI:**  
[dx.doi.org/10.17504/protocols.io.14egn2pqzg5d/v2](https://dx.doi.org/10.17504/protocols.io.14egn2pqzg5d/v2)

**Protocol Citation:** Devansh Agarwal, Kevin W. Mazo, Karl Wahlin 2023. Differentiation of RGC Induced Neurons (RGC-iNs). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.14egn2pqzg5d/v2> Version created by [Devansh Agarwal](#)

## 🌐 Differentiation of RGC Induced Neurons (RGC-iNs) V.2

Kevin W. Karl  
Devansh Agarwal<sup>1,2</sup>, Mazo<sup>3</sup>, Wahlin<sup>2</sup>

<sup>1</sup>Department of Bioengineering, University of California San Diego, La Jolla, United States;

<sup>2</sup>Department of Ophthalmology, University of California San Diego, La Jolla, United States;

<sup>3</sup>Department of Biological Sciences, University of California San Diego, La Jolla, United States

Devansh Agarwal: Correspondence: [d4agarwa@ucsd.edu](mailto:d4agarwa@ucsd.edu);

Kevin W. Mazo: Correspondence: [kewima99@gmail.com](mailto:kewima99@gmail.com);

Karl Wahlin: Correspondence: [kwahlin@ucsd.edu](mailto:kwahlin@ucsd.edu)



Devansh Agarwal

### ABSTRACT

This protocol is designed to convert human induced pluripotent stem cells (PSCs) into retinal ganglion cell induced neurons (RGC-iNs) using a doxycycline-inducible polycistronic transcription factor gene cassette containing human NEUROG2, ATOH7, ISL1, and POU4F2. The TetO-driven transgene cassette is integrated into the CLYBL safe harbor site using a CRISPR-Cas12a ribonucleoprotein. The process of generating neurons is greatly enhanced by the inclusion of the BMP blocker LDN-193189.

### GUIDELINES

Apart from observation under the microscope, counting, and centrifugation, all steps should be carried out in a sterile biological safety cabinet.

### MATERIALS

#### Table: Key resources or reagents required.

A	B	C
REAGENT or RESOURCE	SOURCE	IDENTIFIER
40µm Cell Strainer, EASYstrainer (for generating single cell suspension during plating)	Greiner	Cat# 542040
Accutase (single cell passaging of hPSCs)	Sigma-Aldrich	Cat# A6964
B27 vitamin (neural supplement)	Thermo Fisher Scientific	Cat# 17504044

**MANUSCRIPT CITATION:**

Devansh Agarwal, Nicholas Dash, Kevin W. Mazo, Manan Chopra, Maria PA. Garcia, Amit Patel, Ryan M. Wong, Cairang Jia, Hope Do, Jie Cheng, Colette Chiang, Shawna L. Jurlina, Mike Perry, Jong Rho, Risa Broyer, Cassidy Lee, Robert N. Weinreb, Cezar Gavrilovici, Nick W. Oesch, Derek S. Welsbie, Karl J. Wahlin. Human Retinal Ganglion Cell Neurons Generated by Synchronous BMP inhibition and Transcription Factor Mediated Reprogramming. *npj Regen Med* (2023). DOI: 10.1038/s41536-023-00327-x

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Sep 26, 2023

**Last Modified:** Sep 27, 2023

**PROTOCOL integer ID:**  
88422

A	B	C
BDNF (growth factor for RGC growth and survival)	Qkine	Cat# Qk050
Blebbistatin (ROCK inhibitor for improving cell survival)	Sigma-Aldrich	Cat# B0560
BrainPhys Neuronal Medium (basal media for supporting long-term growth of neurons)	StemCell Technologies	Cat# 05790
CultureOne supplement (for enhancing neural conversion)	Thermo Fisher Scientific	Cat# A3320201
DMEM (basal media )	Thermo Fisher Scientific	Cat# 11965
DMEM/F12 50:50 (basal media)	Thermo Fisher Scientific	Cat# 11330
Doxycycline hyclate (antibiotic for transgene induction)	Sigma-Aldrich	Cat# D5207
F12 (basal media)	Thermo Fisher Scientific	Cat# 11765
GDNF (growth factor for enhancing RGC growth and neuronal survival)	Qkine	Cat# Qk051
Insulin-Human Recombinant (N2 supplement component)	Roche	Cat# 11376497001
L-ascorbic acid (N2 supplement component)	Sigma-Aldrich	Cat# A8960
LDN-193189 (pre-patterning BMP pathway inhibitor)	Sigma-Aldrich	Cat# SML0559
Matrigel® Growth Factor Reduced (GFR) Basement Membrane Matrix (cell attachment/differentiation of hPSCs)	Corning	Cat# 354230
mTeSR1 (maintenance and propagation of hPSCs)	Stem Cell Technologies	Cat# 85850
N-2 Supplement (neural supplement)	Thermo Fisher Scientific	Cat# 17502048
NEAA (non-essential amino acids for supporting neuronal growth)	Thermo Fisher Scientific	Cat# 11140
Nicotinamide (NIC) (vitamin B3 supplement to enhance differentiation)	Sigma-Aldrich	Cat# 72340
Poly-L-ornithine (PLO) hydrobromide (for neural attachment)	Sigma Aldrich	Cat# P3655
Progesterone (N2 supplement component)	Sigma-Aldrich	Cat# P8783
Putrescine dihydrochloride (N2 supplement component)	Sigma-Aldrich	Cat# P5780
Sodium selenite (N2 supplement component)	Sigma-Aldrich	Cat# S5261
Thiazovivin (alternate ROCK inhibitor for cell survival)	LC Labs	Cat# T-9753
holo Transferrin Human (N2 supplement component)	Sigma-Aldrich	Cat# T0665

A	B	C
GlutaMAX Supplement (auxiliary energy source for cells)	Thermo Fisher Scientific	Cat# 35050061

## BEFORE START INSTRUCTIONS

For media/reagent recipes see the last section of the protocol.

## PSC expansion step

- Grow PSCs in mTeSR1 under hypoxia ( $[M] 5\% (v/v)$   $O_2$  /  $[M] 10\% (v/v)$   $CO_2$ ) or normoxia ( $[M] 20\% (v/v)$   $O_2$  /  $[M] 5\% (v/v)$   $CO_2$ ) at  $37^\circ C$ .  
**12-well plate ( $3.5\text{ cm}^2$ ):** Plate 5,000 PSCs into each well of 12-well plates in mTeSR1 in the presence of  $[M] 5\text{ micromolar } (\mu M)$  blebbistatin (blebb; 2,000x stock). Feed daily in mTeSR1 (without blebb) for ~4 more days. Typically, we get ~200,000 – 500,000 cells per well when cells are ready for passaging.  
**6-well plate ( $9.6\text{ cm}^2$ ):** Plate 15,000 PSCs into each well of 6-well plates in mTeSR1 in the presence of blebb. Feed daily in mTeSR1 (without blebb) and grow for ~4 more days. Typically, we can get ~750,000 – 1,000,000 cells per well when cells are ready for passaging.

### Note

For PSC expansion make sure that colonies are not overgrown (70-80 % confluent) and not touching.

Alternatively,  $[M] 2\text{ micromolar } (\mu M)$  thiazovivin (10 mM or 5000x stock) can be used instead of blebb as a ROCK inhibitor.

## Day -1: Priming of stem cells for neural induction

- One day prior to neural induction, pre-coat TC plates overnight with  $[M] 0.1\text{ mg/mL}$  poly-L-ornithine hydrobromide (PLO). To do this, add  $\text{1 mL}$   $[M] 2\text{ mg/mL}$  PLO (20x) to  $\text{19 mL}$  cell culture grade H<sub>2</sub>O, use  $\text{1 mL}$  per well for 6-well plates,  $\text{0.5 mL}$  per well for 12 well plates and  $\text{0.25 mL}$  per well for 24 well plates. Incubate overnight at  $37^\circ C$ . 1d

Pre-treat stem cells by replacing mTeSR1 media with fresh mTeSR1 supplemented with

**1**  $1\ \mu\text{g/mL}$  doxycycline (dox; 1,000x stock) and **1** 100 nanomolar (nM) LDN-193189 (10,000x stock).

#### Note

For long-term experiments, you need better adhesion of cells so you can use a higher concentration of PLO. To do this, dilute **1** 1 mL **1** 2 mg/mL PLO (20x) into **1** 9 mL cell culture grade H<sub>2</sub>O for a final concentration of **1** 0.2 mg/mL.

## Day 0: Neural Induction



















- 3** Wash overnight PLO-coated plates >3 times with culture grade H<sub>2</sub>O. Let plates dry completely in the back of the TC hood for **01:00:00**, then coat with **1** 1 % (v/v) Matrigel. For Matrigel coating, use **1** 1 mL per well for 6-well plates, **1** 0.5 mL per well for 12 well plates and **1** 0.25 mL per well for 24 well plates. Incubate for > **03:00:00** at **37 °C**. 4h
- 4** Prepare Neural Induction Medium (NIM) initiation cocktail with **1** 2  $\mu\text{g/mL}$  dox (500x stock), **1** 100 nanomolar (nM) LDN (10,000x stock), 1x CultureOne (100x stock) and **1** 5 micromolar ( $\mu\text{M}$ ) blebb.

#### Note

Dox is light-sensitive, so keep aliquots cool after thawing and dark when not in use (up to ~1 month at **4 °C**).

In addition, it is critical to prevent the drying of matrigel-coated dishes after the matrigel has been added.





- 5** Incubate the cells with prewarmed Accutase (~ **00:05:00**) in the incubator for **00:12:00**. The volume of Accutase to use is 1/2 the volume that you maintain the cells in (e.g., **1** 1 mL per well of a 6 well plate, **1** 0.5 mL per well of a 12 well plate and **1** 0.25 mL per well of a 24 well plate). 17m
- 6** Gently rinse the wells using a P1000 micropipette and pipet up and down 3 times to further break up the cell clumps into single cells.

- 7 Put the cells into a conical  1.5 mL or  5 mL tube with 2 times the volume of prewarmed mTeSR1 with  5 micromolar ( $\mu\text{M}$ ) blebb (e.g.  1 mL Accutase +  2 mL mTeSR1) to quench the Accutase, then pellet the cells for  00:05:00 at  80 x g. 5m
- 8 Aspirate the supernatant and then resuspend the cell pellet in prewarmed NIM +  5 micromolar ( $\mu\text{M}$ ) blebb. 2m
- 9 Filter cells with a  40  $\mu\text{m}$  cell strainer. Count the cells with a hemocytometer. 5m
- 10 Plate cells onto PLO/matrigel-coated plates. 5m
- To make a 6-well plate (**9.6  $\text{cm}^2$ /well; 57.6  $\text{cm}^2$  total; 7,000 cells/ $\text{cm}^2$** ): Add 403,200 cells into  12 mL (67,200 cells per well) of NIM initiation cocktail (NIM + dox, LDN, CultureOne, blebb) in a  15 mL conical tube, mix well and distribute across the wells ( 2 mL per well).
- To make a 12-well plate (**3.5  $\text{cm}^2$ /well; 42  $\text{cm}^2$  total; 7,000 cells/ $\text{cm}^2$** ): Add 294,000 cells into  12 mL (24,500 cells per well) of NIM initiation cocktail in a  15 mL conical tube, mix well and distribute across the wells ( 1 mL per well).
- To make a 24-well plate (**1.9  $\text{cm}^2$ /well; 45.6  $\text{cm}^2$  total; 7,000 cells/ $\text{cm}^2$** ): Add 319,200 cells into  12 mL (13,300 cells per well) of NIM initiation cocktail in a  15 mL conical tube, mix well and distribute across the wells ( 0.5 mL per well).

## Day 1: Maintenance

- 11 Do nothing.

## Day 2: Feed - add 1/3 of media

- 12 For 6 well plate: Add  1 mL **NIM** +  1 µg/mL **dox** + **1xCultureOne** to each well. 5m
- For 12 well plate: Add  0.5 mL **NIM** + **dox** + **CultureOne** to each well.
- For 24 well plate: Add  0.25 mL **NIM** + **dox** + **CultureOne** to each well.









#### Note

Do this very carefully by adding media to the sides of the dish. If you are not very careful cells will detach.

## Day 3: Maintenance

- 13 Do nothing.

## Day 4: Feed - exchange 1/3 of media

- 14 For 6 well plate: Remove  1 mL media and replace with fresh  1 mL **NIM** +  1 µg/mL **dox** + **1xCultureOne** + **NIC** ( 10 millimolar (mM) ). 5m
- For 12 well plate: Remove  0.5 mL media and replace with fresh  0.5 mL **NIM** + **dox** + **CultureOne** + **NIC**.
- For 24 well plate: Remove  0.25 mL media and replace with fresh  0.25 mL **NIM** + **dox** + **CultureOne** + **NIC**.

**Beyond day 4, plates need to be fed every other day.**







#### Note

This is the last day to add CultureOne.

#### **For all subsequent exchanges:**

Neurons tend to easily dissociate from the dish, so be very careful when aspirating. Take care to aspirate and dissociate by tilting the dish so that the medium accumulates on one side. Then, aspirate/dispense with the pipette directed toward the wall of the dish (i.e., away from the cells at the bottom).







## Day 6: Feed - exchange 1/3 of media

- 15 For 6 well plate: Remove  1 mL media and replace with fresh  1 mL **BrainPhys + B27 (1x)** + **[1 μg/mL dox + BDNF ( [50 ng/mL] ) + GDNF ( [10 ng/mL] ) + NIC ( [10 millimolar (mM)] )**. 5m
- For 12 well plate: Remove  0.5 mL media and replace with fresh  0.5 mL **BrainPhys + B27 + dox + BDNF + GDNF + NIC**.
- For 24 well plate: Remove  0.25 mL media and replace with fresh  0.25 mL **BrainPhys + B27 + dox + BDNF + GDNF + NIC**.

### Note

This is the last day to add dox.

## Day 8: Feed - exchange 1/3 of media

- 16 For 6 well plate: Remove  1 mL media and replace with fresh 1  1 mL **BrainPhys + B27 (1x) + BDNF ( [50 ng/mL] ) + GDNF ( [10 ng/mL] ) + NIC ( [10 millimolar (mM)] )**. 5m
- For 12 well plate: Remove  0.5 mL media and replace with fresh  0.5 mL **BrainPhys + B27 + BDNF + GDNF + NIC**.
- For 24 well plate: Remove  0.25 mL media and replace with fresh  0.25 mL **BrainPhys + B27 + BDNF + GDNF + NIC**.
- 17 For long-term experiments >1 week, continue feeding by 1/3 media exchange every other day with **BrainPhys + B27 + BDNF + GDNF + NIC**. 5m

## Media/Reagent Recipes

- 18 **Neural Induction medium (NIM medium):**


A	B	C
Component	Catalog #	Volume
DMEM/F12, HEPES	Thermo Fisher Scientific # 11330	485 ml
N2 supplement (100x)	Thermo Fisher Scientific # 17502048 or recipe below	5 ml


A	B	C
NEAA (non-essential amino acids, 100x)	Thermo Fisher Scientific # 11140	5 ml
GlutaMAX supplement (100x)	Thermo Fisher Scientific # 35050061	5 ml
		Total 500 ml



## 19 N2 supplement (100X) recipe - If making in-house add the following:

A	B	C	D	E	F	G
Component	Catalog #	100 ml	150 ml	200 ml	250 ml	500 ml
Transferrin	Sigma-Aldrich # T0665	1 g	1.5 g	2 g	2.5 g	5.0 g
Insulin	Roche # 11376497001	50 mg	75 mg	100 mg	125 mg	250 mg
Progesterone	Sigma-Aldrich # P8783	63 µg	94.5 µg	126 µg	157 µg	315 µg
Putrescine	Sigma-Aldrich # P5780	161 mg	241.5 mg	322 mg	402.5 mg	805 mg
Sodium Selenite	Sigma-Aldrich # S5261	50.2 µg	75.3 µg	100.4 µg	125.5 µg	251 µg
DMEM/F12	Thermo Fisher Scientific # 11330	to 100 ml	to 150 ml	to 200 ml	to 250 ml	to 500 ml

## 20 Reagent Stock Dilutions:

**Recombinant Human BDNF Protein** (Qkine, Cat# Qk050): [1M] 50 µg/mL stock in [1M] 10 millimolar (mM) HCl with [1M] 0.1 % (v/v) BSA; 1,000X; use at [1M] 50 ng/mL. Store aliquots in  -80 °C.

**Recombinant Human GDNF Protein** (Qkine, Cat# Qk051): [1M] 10 µg/mL stock in cell culture grade H2O with [1M] 0.1 % (v/v) BSA; 1,000X; use at [1M] 10 ng/mL. Store aliquots in  -80 °C.

**Doxycycline (dox)** (Sigma-Aldrich, Cat# D5207): Make [1M] 1 mg/mL stock (working concentration is [1M] 0.5-2 µg/mL) in cell culture grade ddH2O and filter sterilize; 1,000X; use at [1M] 1 µg/mL. Store  1 mL aliquots at  -20 °C. Protect from light.

**LDN-193189 hydrochloride (LDN)** (Sigma Cat# SML0559-5MG): [1M] 1 millimolar (mM) (10,000x):





$g = \text{Molecular Weight (g/mol)} * \text{Molarity (M)} * \text{Volume (L)}$



$0.005 \text{ g} = (406.48 \text{ g/mol}) * (0.001 \text{ M}) * (\text{Volume})$

$0.005 \text{ g} = 0.41 \text{ g/L} * \text{Volume}$




Volume =  0.012 L or  5 mg LDN in  12 mL DMSO; Store aliquots in  -80 °C .

**Matrigel (GF reduced)** ( 1 % (v/v) ) (Corning, Cat# 354230):


Thaw stock 10ml bottle overnight  On ice (wet) before aliquoting. Always keep matrigel and tubes on ice and never let it come to room temperature or it will gel.

Make  200 µL aliquots. Store aliquots in  -80 °C .




Add  200 µL matrigel to  24 mL ice-cold DMEM/F12 (~  1 % (v/v) final).

Add  1 mL per well of 6-well plate,  0.5 mL per well of 12-well plate or  0.25 mL per well of 24-well plate.

**Nicotinamide (NIC)** (Sigma Cat# 72340):  1 Molarity (M) (100x) stock solution (



 10 millimolar (mM) working solution). Soluble in water to ~ 1g/10ml.



$g = \text{Molecular Weight (g/mol)} * \text{Molarity (M)} * \text{Volume (L)}$ ;  $g = (122.12 \text{ g/mol}) * (1 \text{ M}) * (0.05 \text{ L})$

=  6.11 g NIC in  50 mL of DMEM/F12 (or water); filter sterilize and store at  4 °C .

**poly-L-ornithine hydrobromide (PLO)** - (Sigma Cat# P3655-500MG; mol wt 30,000-70,000 Da):

 2 mg/mL (20x) stock:

Add  500 mg PLO to  250 mL cell culture ddH2O.

Make  1 mL aliquots. Store aliquots in  -80 °C .

**Blebbistatin (blebb)** ( 10 millimolar (mM) or 2,000x stock, Sigma Cat# B0560-1MG):

$g = \text{Molecular Weight (g/mol)} * \text{Molarity (M)} * \text{Volume (L)}$

$0.001 \text{ g} = (292.33 \text{ g/mol}) * (0.01 \text{ M}) * (\text{Volume})$

$0.001 \text{ g} = 2.92 \text{ g/L} * \text{Volume}$

Volume =  $0.000342 \text{ L} = 0.342 \text{ mL} = 342 \text{ µL}$

Dissolve  1 mg blebbistatin into  342 µL DMSO, divide into  20 µL aliquots (store at  -80 °C ). Working concentration is  5 micromolar (µM) .