



OCT 20, 2023

OPEN ACCESS



**DOI:**

[dx.doi.org/10.17504/protocols.io.ewov1q31ogr2/v1](https://dx.doi.org/10.17504/protocols.io.ewov1q31ogr2/v1)

**Protocol Citation:** Lauren D Walter, Benjamin D Cosgrove 2023. LDW\_BDC Muscle scRNAseq. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.ewov1q31ogr2/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Oct 20, 2023

**Last Modified:** Oct 20, 2023

## 🌐 LDW\_BDC Muscle scRNAseq

Lauren D Walter<sup>1</sup>, Benjamin D Cosgrove<sup>1</sup>

<sup>1</sup>Cornell University

SenNet JMN-MSMP



Lauren D Walter  
Cornell University

### ABSTRACT

single-cell RNA sequencing of muscle

## Skeletal muscle injury and collection


- 1 In a given mouse, both tibialis anterior (TA) muscles were injected with 10 µl of notexin (10 µg/ml; Latoxan, France).
- 2 At various days post injury, the mouse was sacrificed and the TAs were collected. Each TA was processed independently to generate single cell suspensions.

## Generate single-cell suspensions

- 3 Muscles were enzymatically digested with 8 mg/ml Collagenase D (Roche, Basel, Switzerland) and 10 U/ml Dispase II (Roche, Basel, Switzerland) and then manually dissociated to generate cell suspensions.
- 4 Myofiber debris was removed by filtering the cell suspensions through a 100 µm and then a 40 µm filter (Corning Cellgro # 431752 and # 431750).
- 5 After filtration, erythrocytes were removed by incubating the cell suspension in erythrocyte lysis buffer (IBI Scientific # 89135-030).
- 6 After digestion, the single-cell suspensions were washed and resuspended in 0.04% BSA in PBS at a concentration of  $10^6$  cells/ml. A hemocytometer was used to manually count the cells to determine the concentration of the suspension.

## Single-cell RNA library preparation and sequencing

- 7 Single-cell RNA-sequencing libraries were prepared using the Chromium Single Cell 3' reagent kit v3 (10x Genomics, Pleasanton, CA) following the manufacturer's protocol.

- 
- 8 Cells were diluted into the Chromium Single Cell A Chip to yield a recovery of 6,000 single-cell transcriptomes with <5% doublet rate.
  - 9 Libraries were sequenced on the NextSeq 500 (Illumina, San Diego, CA).