



MAR 01, 2024

OPEN ACCESS



Protocol Citation: Renuka Ravi Gupta, Nona Farbehi, Helaine Grazielle Santos Vieira, Helen E. King, henderson, Vikram Khurana, Gist Croft, Robert J Weatheritt, Lorenz Studer, Joseph Powell 2024. LENTIVIRAL PRODUCTION FOR PCRISPRi DUAL GUIDE mDA NEURON LIBRARY. [protocols.io](https://protocols.io/view/lentiviral-production-for-pcrispr-i-dual-guide-mda-c9uez6te) <https://protocols.io/view/lentiviral-production-for-pcrispr-i-dual-guide-mda-c9uez6te>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development
We are still developing and optimizing this protocol

Created: Feb 27, 2024

🌐 LENTIVIRAL PRODUCTION FOR PCRISPRi DUAL GUIDE mDA NEURON LIBRARY

Renuka Ravi Gupta^{1,2,3,4}, Nona Farbehi^{1,2,3,5}, Helaine Grazielle Santos Vieira¹, Helen E. King¹, henderson^{3,6}, Vikram Khurana^{3,7}, Gist Croft^{3,8}, Robert J Weatheritt¹, Lorenz Studer^{3,6}, Joseph Powell^{1,2,3,4}

¹Garvan Institute of Medical Research, Sydney, NSW 2010, Australia;

²Garvan Weizmann Center for Cellular Genomics, Garvan Institute of Medical Research, Sydney, NSW 2010, Australia;

³Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815, USA;

⁴School of Medical Science, University of New South Wales, Sydney, NSW, 2052, Australia;

⁵Graduate School of Biomedical Engineering, University of New South Wales, Sydney, NSW, 2052, Australia;

⁶The Centre for Stem Cell Biology, Developmental Biology Program, Sloan Kettering Institute for Cancer Research, New York, NY, USA;

⁷Ann Romney Center for Neurologic Diseases, Brigham and Women's Hospital and Harvard Medical School, Boston, MA 02115, USA;

⁸The New York Stem Cell Foundation Research Institute, New York, NY, USA



Renuka Ravi Gupta

ABSTRACT

This protocol outlines the production of lentiviral supernatant from a CRISPRi Library plasmid. The plasmid is transfected into HEK293T cells and used for CRISPRi perturbation in both human pluripotent stem cells (hPSCs) and hPSC-derived dopaminergic neural cells. To enhance the viral titer, the virus is concentrated using the Lenti-X concentrator from Takara-Bio following collection of the viral supernatant from HEK 293T cells.

ATTACHMENTS

[LENTIVIRAL PRODUCTION FOR pCRISPRi dual-guide mDA Neuron Library.docx](#)

Last Modified: Mar 01, 2024

MATERIALS

PROTOCOL integer ID: 95846

Keywords: ASAPCRN, Lentiviral Production, CRISPRi machinery, dopaminergic neuron differentiation, Perturb-Seq

Funders Acknowledgement:
ASAP
Grant ID: ASAP-000472

| A | B | C |
|---|--|-----------|
| MATERIAL | COMPANY | CATALOG |
| T25 Flask | Corning | 430639 |
| 15ml polypropylene centrifuge tubes | Falcon | 352096 |
| 5ml serological pipettes | Corning | 4487 |
| 10ml serological pipettes | Corning | 4488 |
| DNA Low-bind tubes 1.5ml | Eppendorf | 022431021 |
| DMEM High glucose | ThermoFisher Scientific | 11965092 |
| FBS | Bovogen | 2008A |
| Lipofectamine 3000 | ThermoFisher Scientific | L3000015 |
| psPAX2 | Addgene | 12260 |
| pCAG-VSV-G | Addgene | 35616 |
| Transfer Plasmid (Plasmid of interest) pCRISPRi dual-guide mDA Neuron Library | Prepared by Robert Weatherhitt Lab at Garvan | NA |
| Opti-MEM I Reduced Serum Medium | ThermoFisher Scientific | 31985062 |
| Lenti-X Concentrator | Takara-Bio | 631231 |

1. Reagent Composition for HEK293T culture media

| A | B |
|-----------|--------|
| REAGENT | AMOUNT |
| DMEM (1X) | 450ml |
| 10% FBS | 50 ml |

1. Aliquot the media into 50 ml tubes and store it at 4C.
2. Do not use Pen-Strep in your media for transfection protocol.

Day 0: Seeding HEK 293T

- 1 Seed HEK 293T cells into T25 flasks in 5ml of HEK culture media at a cell seeding density of 1.5×10^6 cells/flask.

Day 1: Transfection OF CRISPRi plasmid pool

- 2 After 24 hours, check whether the HEK293T cells are at least 60% confluent and proceed with transfection and bring the Opti-MEM to room temperature.

- 3 In TUBE A, mix the following reagents. Make this volume of reagent per T25 flask.

| A | B |
|--------------------|-------------------|
| TUBE A | AMOUNT IN μ l |
| Lipofectamine 3000 | 15 |
| Opti-MEM | 235 |
| Total | 250 |

- 4 In TUBE B, mix the following reagents. Make this volume of reagent per T25 flask.

| A | B |
|--|-------------|
| TUBE B | AMOUNT |
| psPAX2 | 4 μ g |
| VSV-G | 2 μ g |
| pCRISPRi dual-guide mDA Neuron Library | 4 μ g |
| P3000 | 20 μ l |
| Opti-MEM | X μ l |
| Total | 250 μ l |

- 5 Transfer tube B contents to tube A, mix gently by pipetting and incubate for 10 mins at room temperature.

- 6 Replace the media in the T25 flask with fresh 5 ml of HEK culture media during incubation.
- 7 Add the whole mix (~500ul) gently, drop by drop, making sure to cover the entire area of the T25 flask.
- 8 Gently swirl the media in the flask to ensure the proper distribution of the packaging mixture, and incubate the flask in the incubator for 24 hours.

Day 2: Media change

- 9 After 24 hours, aspirate the media and gently add 5 ml of fresh HEK culture media

Note

Caution: HEK293T cells tend to lift off easily, so be extremely careful while changing the media and add the media on the opposite wall of the flask where the cells are not attached.

- 10 Return the flask to the incubator.

Day 3: Harvesting the supernatant and concentrating the virus

- 11 Harvest the lentivirus-containing supernatant into 15 ml tubes.
- 12 Centrifuge the supernatant at 500 g for 10 mins.

13 Transfer the clarified supernatant (5 ml) into a sterile 15 ml falcon tube.

14 Depending on the volume of the supernatant, combine 1 volume of Lenti-X Concentrator with 3 volumes of clarified supernatant.

| A | B |
|---------------------------------|------------------------|
| VOLUME OF CLARIFIED SUPERNATANT | VOLUME OF CONCENTRATOR |
| x = 5 ml | x/3 = 1.6 ml |

15 Mix by gentle inversion.

16 Incubate mixture overnight at 4C.

Day 4: Centrifuge and Resuspension of pellet

17 Centrifuge the concentrated virus at 1500g for 45 mins at 4C.

Note

After the centrifugation, an off-white pellet will be visible.

18 Carefully remove the supernatant without disturbing the pellet. The residual supernatant can be removed with a pipette tip or by brief centrifugation at 1500g.

19 Gently resuspend the pellet in 1/10th of the original volume (500ul) with DMEM.

Note

The pellet will be sticky at first but will suspend quickly.

20 Detailed protocol for concentration for the virus can be found in this link: [Lenti-X™ Concentrator Protocol-at-a-Glance](#)