



Oct 17, 2022

Modified protocol to improve Bodo saltans yield in culture.

Zhu-Hong Li¹, Fatma Gomaa², Virginia Edgcomb³, Roberto Docampo⁴¹The University of Georgia; ²Harvard University; ³Woods Hole Oceanographic Institution;⁴University of Georgia

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.9vyh67w

Protist Research to Optimize Tools in Genetics (PROT-G)

Zhu-Hong Li
The University of Georgia

ABSTRACT

This protocol is a modified version to improve the yield of Bodo saltans cell density in culture. The original protocol is: <https://www.protocols.io/view/bodo-saltans-culture-protocol-sh6eb9e>

ATTACHMENTS

[Modified protocol to improve Bodo Saltans yield in culture.docx](#)

DOI

dx.doi.org/10.17504/protocols.io.9vyh67w

EXTERNAL LINK

<https://doi.org/10.1111/1462-2920.15918>

PROTOCOL CITATION

Zhu-Hong Li, Fatma Gomaa, Virginia Edgcomb, Roberto Docampo 2022.
Modified protocol to improve Bodo saltans yield in culture.. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.9vyh67w>



MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

Gomaa F, Li Z, Beaudoin DJ, Alzan H, Girguis PR, Docampo R, Edgcomb VP, Array. Environmental Microbiology 24(7). doi: [10.1111/1462-2920.15918](https://doi.org/10.1111/1462-2920.15918)

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Nov 28, 2019

LAST MODIFIED

Oct 17, 2022

PROTOCOL INTEGER ID

30360

- 1 Prepare the bacteria-Bodo saltans medium as described in the protocol (<https://www.protocols.io/view/bodo-saltans-culture-protocol-sh6eb9e>) .
- 2 Collect Bodo cells from a T25 tissue culture flask by centrifugation at 1200 x g for 6 minutes. After removing the supernatant, add deionized water to resuspend the pelleted cells. After two washes with water, resuspend the cells in culture medium. Usually the Bodo saltans/bacteria ratio is >1 after two washes, as determined by FACS analysis.
- 3 Inoculate a T25 tissue culture flask (50 ml) containing 10 to 15 ml of fresh medium with 10 mg/ml puromycin. Puromycin at this concentration has no inhibitory effect on the growth of Bodo saltans but it can slow down bacterial growth. Transfer 100 ml-0.5 ml cells from step 2 into the flask. Incubate horizontally at 18C with loosely adjusted cap.
- 4 After 2-4 days, when the density of Bodo cells reaches to $\sim 3 \times 10^6/\text{ml}$, collect and wash the cells for downstream applications.
- 5 The advantages of this protocol are improved yield and purity. 10 mg/ml puromycin can efficiently prevent bacterial population overgrowth in the flask. The medium is much cleaner and there are no aggregates formed by overgrown bacteria. Usually, we can only get $1 \times 10^6/\text{ml}$ of Bodo cells using conventional culture conditions, and need filtration and 4-5 washes to remove all bacteria from the culture for other applications, such as electroporation or DNA extraction. After these modifications Bodo cells can reach $3 \times 10^6/\text{ml}$ in the culture resulting in a purer population of cells after only 2 washes.
- 6

