



2 ▼

Dec 02, 2021

RNA Extraction Protocol for Shorea V.2

mpfsum¹¹University of Malaya

1



dx.doi.org/10.17504/protocols.io.b2ipqcdn

MP2Lab



mpfsum

RNA extraction protocol using CTAB method optimized for leaf and bud samples from *Shorea curtisii*.

Adapted from extraction protocol for *Shorea beccariana* (see attached publication).

DOI

dx.doi.org/10.17504/protocols.io.b2ipqcdn

mpfsum 2021. RNA Extraction Protocol for Shorea. **protocols.io**<https://dx.doi.org/10.17504/protocols.io.b2ipqcdn>

mpfsum



protocol

Kobayashi MJ, Takeuchi Y, Kenta T, Kume T, Diway B, Shimizu KK. Mass flowering of the tropical tree *Shorea beccariana* was preceded by expression changes in flowering and drought-responsive genes. *Mol Ecol*. 2013 Sep;22(18):4767-82. doi: 10.1111/mec.12344. Epub 2013 May 8. PMID: 23651119; PMCID: PMC3817532.

protocol ,

Dec 02, 2021

Dec 02, 2021

55599

Lysis

1h 10m

1 Weigh BUD **40 mg (25-40 mg)** or LEAF **60 mg (50-70 mg)** and transfer into a 1.5 mL

15m

tube.

- 2 Ground in CTAB buffer (3% CTAB, 1.4 M NaCl, 20 mM EDTA, 100 mM Tris-HCl pH 8.0, 0.2% ^{10m}β-mercaptoethanol) using tissuelyzer.

- 3 Incubate at **60 °C** for **00:45:00** . 45m

Precipitation 25m

- 4 Add 1 volume of chloroform. 5m
Mix by vortex but not too vigorous.
Ensure that the two phases of the mixture are homogenized well.

- 5 Centrifuge at **5000 rpm, Room temperature, 00:01:00** . 5m
Transfer the supernatant into a new 1.5 mL tube.

- 6 Add 2/3 volume of isopropanol. 5m
Mix gently (turning the tubes upside down **00:02:00**).

- 7 Centrifuge **5000 rpm, Room temperature, 00:05:00** 10m
Discard solution carefully.



Formation of white to translucent pellet.

Washing 25m

- 8 Add **700 µL** ethanol. 5m
Mix gently (flicking the tube) until the pellet no longer stick to the tube.

- 9 Centrifuge **5000 rpm, Room temperature, 00:05:00** . 10m
Discard solution carefully.

- 10 Air dry the pellet 🕒00:01:00 .
Repeat Washing if necessary.

5m

Resuspension 25m

- 11 Resuspend the pellet in 📏50 µL RNase-free water.
Flick the tube until the pellet is no longer visible.

If the pellet is difficult to dissolve, leave it for ~30 min before flicking again.
If it still does not dissolve, add more RNase-free water AND adjust the subsequent steps accordingly.

Genomic DNA Digestion 55m

- 12 Removed gDNA using
🔗Turbo DNA-free Kit Invitrogen - Thermo
Fisher Catalog #AM1907
Follow the manufacturer's protocol.

Purification 40m

- 13 🔗RNeasy Plant Mini
Purification is done using 📦Kit Qiagen Catalog #74904
adjustment to the protocol.
Add 📏350 µL RLT buffer.
Mix gently (flicking).

5m

with minor

- 14 Add 📏250 µL absolute ethanol.
Mix gently (flicking).

5m












- 15 Transfer the solution into RNeasy Mini Spin column.
Centrifuge 🕒10000 rpm, Room temperature, 00:00:15 .
Discard the flow through.

3m

- 16 Add 📏700 µL RW1 buffer to the column.
Centrifuge 🕒10000 rpm, Room temperature, 00:00:15 .

3m

Discard the flow through.

- 17 Add  **500 µL** RPE buffer to the column. 3m
Centrifuge  **10000 rpm, Room temperature, 00:00:15** .
Discard the flow through.
- 18 Add  **500 µL** RPE buffer to the column. 5m
Centrifuge  **10000 rpm, Room temperature, 00:02:00** .
Discard the flow through.
- 19 Replace the collection tube. 3m
Centrifuge  **14000 rpm, Room temperature, 00:01:00** to further dry the membrane.
- 20 Replace the collection tube with a new 1.5 mL tube. 8m
Add  **50 µL** RNase-free water directly onto the membrane.
Incubate at  **Room temperature** for  **00:05:00** .
- 21 Centrifuge  **10000 rpm, Room temperature, 00:01:00** . 5m
Store the extracted RNA in  **-20 °C** or  **-80 °C** for longer storage.

Quality Assessment

- 22 Evaluate the RNA quality using Agilent 2100 Bioanalyzer.