



Sep 08, 2020

Colony PCR -- CHEM 584

Ken Christensen¹¹Brigham Young University

In Development

dx.doi.org/10.17504/protocols.io.bk5xky7nKen Christensen
Brigham Young University

ABSTRACT

A quick way to screen transformants for an insert. This is faster and less expensive than doing a restriction digest.

DOI

dx.doi.org/10.17504/protocols.io.bk5xky7n

PROTOCOL CITATION

Ken Christensen 2020. Colony PCR -- CHEM 584. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bk5xky7n>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 08, 2020

LAST MODIFIED

Sep 08, 2020

PROTOCOL INTEGER ID

41879

After you have colonies from a transformation

- 1 Draw grid on clean agar plate (use the appropriate antibiotic). You will use this to inoculate an overnight culture for a positive clone.
- 2 Select colonies to pick. Streak portion of colony to numbered sector and place the remainder in a correspondingly numbered PCR tube with 50 μ L of autoclaved water.
- 3 Heat at 95 °C for 10 minutes. *** This can be done in PCR machine, a heat block, or a boiling water bath*
- 4 Spin solution on high for 10 minutes to pellet cellular debris, and remove 4 μ L for PCR.
- 5 Setup and run the PCR reaction, then assess the results using agarose gel electrophoresis. Be sure to include a positive control (e.g., plasmid) if available. Negative controls are also useful if one is available.

