

Oct 09, 2024

## Sinai SCENT TMC- Methylation Array (EPIC)

DOI

[dx.doi.org/10.17504/protocols.io.81wgbzqyygpk/v1](https://dx.doi.org/10.17504/protocols.io.81wgbzqyygpk/v1)

Sojin Kim<sup>1</sup>

<sup>1</sup>Icahn School of Medicine at Mount Sinai



Sojin Kim

Icahn School of Medicine at Mount Sinai

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.81wgbzqyygpk/v1](https://dx.doi.org/10.17504/protocols.io.81wgbzqyygpk/v1)

External link: [https://mtsinai-my.sharepoint.com/personal/sojin\\_kim\\_mssm\\_edu/\\_layouts/15/onedrive.aspx?id=%2Fpersonal%2Fsojin%5Fkim%5Fmssm%5Fedu%2FDocuments%2FExperiments%2FMethylation%20Array%5FBiorepository%2Ffinanium%2Dhd%2Dmethylation%2Dguide%2D15019519%2D01%2Epdf&parent=%2Fpersonal%2Fsojin%5Fkim%5Fmssm%5Fedu%2FDocuments%2FExperiments%2FMethylation%20Array%5FBiorepository](https://mtsinai-my.sharepoint.com/personal/sojin_kim_mssm_edu/_layouts/15/onedrive.aspx?id=%2Fpersonal%2Fsojin%5Fkim%5Fmssm%5Fedu%2FDocuments%2FExperiments%2FMethylation%20Array%5FBiorepository%2Ffinanium%2Dhd%2Dmethylation%2Dguide%2D15019519%2D01%2Epdf&parent=%2Fpersonal%2Fsojin%5Fkim%5Fmssm%5Fedu%2FDocuments%2FExperiments%2FMethylation%20Array%5FBiorepository)

**Protocol Citation:** Sojin Kim 2024. Sinai SCENT TMC- Methylation Array (EPIC). protocols.io  
<https://dx.doi.org/10.17504/protocols.io.81wgbzqyygpk/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** July 31, 2024

**Last Modified:** October 09, 2024

**Protocol Integer ID:** 104381



## Abstract

High throughput, proven procedure for bisulfite conversion of DNA  
96-well desulphonation and recovery of bisulfite-treated DNA



## High-Throughput Bisulfite Conversion of DNA

- 1 Add 5  $\mu$ l of M-Dilution Buffer to each DNA sample in a Conversion Plate and adjust the total volume to 50  $\mu$ l with water. Mix each sample by pipetting up and down.  
  
Example: For 14  $\mu$ l of a DNA sample add 5  $\mu$ l M-Dilution Buffer and 31  $\mu$ l water.
- 2 Incubate the Conversion Plate containing the samples at 37°C for 15 minutes.
- 3 After the above incubation, add 100  $\mu$ l of the prepared CT Conversion Reagent to each sample and mix.  
  
Note: The CT Conversion Reagent is light sensitive, so try to minimize the reaction's exposure to light whenever possible.
- 4 Incubate the Conversion Plate in the dark at 50°C for 12-16 hours (e.g., using a thermal cycler).  
  
Note: See Appendix (page 11) for alternative incubation conditions (e.g., when using the Illumina Infinium® Methylation Assay).
- 5 Incubate the sample at 0-4°C (e.g., on ice or using a thermal cycler) for 10 minutes. Samples may be kept at 4°C for up to 20 hours.
- 6 Add 400  $\mu$ l of M-Binding Buffer to each well of a Zymo-Spin™ I-96 Binding Plate on a Collection Plate.  
  
Note: The capacity of each well of the Binding Plate is 1.1 ml. The capacity of each well of the Collection Plate is 800  $\mu$ l. Empty the Collection Plate whenever necessary to prevent contamination of the Binding Plate contents by flow-through.
- 7 Load the samples (from Step 5) into the wells of the Zymo-Spin™ I-96 Binding Plate containing the M-Binding Buffer. Mix by pipetting up and down.
- 8 Centrifuge at  $\geq 3,000 \times g$  (5,000  $\times g$  max.) for 5 minutes. Discard the flow-through.
- 9 Add 400  $\mu$ l of M-Wash Buffer to each well and centrifuge at  $\geq 3,000 \times g$  for 5 minutes.
- 10 Add 200  $\mu$ l of M-Desulphonation Buffer to each well and let stand at room temperature (20-30°C) for 15-20 minutes. After this incubation, centrifuge at  $\geq 3,000 \times g$  for 5 minutes.



## Infinium HD Assay Methylation Protocol

- 11 Follow the Illumina Infinium HD Assay Methylation Protocol Guide. (see the reference)

### Protocol references

Zymo Research EZ-96 DNA Methylation Kit (Deep-Well Format)

Infinium® HD Assay Methylation Protocol Guide