



Jul 31, 2021

PROCEDURE TO ISOLATE AND CULTURE NEURONS FROM EMBRYONIC MOUSE CORTEX

In 1 collection

Odetta Antico¹, Miratul M. Muqit¹¹Medical Research Council Protein Phosphorylation and Ubiquitylation Unit, School of Life Science, University of Dundee, Dow Street, Dundee, DD1 5EH, UK

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.bsr6nd9e

ASAP2020

alessi

m.muqit

ABSTRACT

Mutations in PINK1 cause early-onset Parkinson's disease. PINK1 becomes stabilised and active upon mitochondrial depolarisation. This leads to phosphorylation of ubiquitin and Parkin via Serine 65 residues and a feed forward mechanism whereby PINK1 phosphorylates newly formed polyubiquitin chains, generating phospho-ubiquitin, which further promotes Parkin recruitment and activation. Once activated, Parkin ubiquitylates proteins at the outer face of the outer mitochondrial membrane (OMM) and then initiates a downstream pathway that eventually leads to mitophagy, a mitochondria-specific type of autophagy. Notably, much of previous investigation into PINK1/Parkin activity has been performed in non-neuronal human cancer cells where Parkin and/or PINK1 is over-expressed. Here we report a protocol for generation of mouse embryonic cortical neuronal cultures that produce high cell yields and can be used for studying endogenous PINK1 and Parkin signalling by biochemical methods and proteomics.

ATTACHMENTS

[Mouse PINK1 pathway protocol \(166 - 337\).pdf](#)

DOI

dx.doi.org/10.17504/protocols.io.bsr6nd9e

PROTOCOL CITATION

Odetta Antico, Miratul M. Muqit 2021. PROCEDURE TO ISOLATE AND CULTURE NEURONS FROM EMBRYONIC MOUSE CORTEX . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bsr6nd9e>

COLLECTIONS ①

**Cell-based analysis of PINK1-Parkin pathway activation in primary mouse cortical neurons**

KEYWORDS

Neurons, PINK1, Parkin, Mitochondrial stress, ubiquitin

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Feb 25, 2021

LAST MODIFIED

Jul 31, 2021

OWNERSHIP HISTORY

Feb 25, 2021  Jesintha Maniraja

Feb 25, 2021  Urmilas

Mar 01, 2021  Urmilas

May 05, 2021  m.muqit

PROTOCOL INTEGER ID

47646

PARENT PROTOCOLS

Part of collection

[Cell-based analysis of PINK1-Parkin pathway activation in primary mouse cortical neurons](#)

MATERIALS TEXT

For culture

1. E16.5 mouse embryos (8-10 embryos, either sex; we routinely use C57BL/6j mice or PINK1 wild-type and knockout mice)

CRITICAL All experiments must be conducted in accordance with the relevant institutional and governmental guidelines and regulations.

2. Dissection medium:

 [HBSS, calcium, magnesium, no phenol red Gibco - Thermo](#)

Fisher Catalog #14025050

3. Digestion medium:

 [Trypsin-EDTA Gibco - Thermo](#)

Fisher Catalog #25300054

 [DNase I Merck Millipore](#)

Sigma Catalog #11284932001

in HBSS

A	B
Trypsin-EDTA	0.025%
DNase I	0.125 mg/mL


4. Dissociation medium:

 [Neurobasal™ Medium Gibco - Thermo](#)

Fisher Catalog #21103049

 [Fetal Bovine Serum qualified heat inactivated Brazil Gibco - Thermo](#)

Fisher Catalog #10500064

 [B-27™ Supplement \(50X\) serum free Gibco - Thermo](#)

Fisher Catalog #17504044

[GlutaMAX™ Supplement Gibco - Thermo](#)

Fisher Catalog #35050061

A	B
Neurobasal medium	
Foetal Bovine Serum (FBS) heat-inactivated	10%
B27 supplement, serum free	1X
GlutaMAX supplement	1%

5. Culturing medium:

[Neurobasal™ Medium Gibco - Thermo](#)

Fisher Catalog #21103049

[B-27™ Supplement \(50X\) serum free Gibco - Thermo](#)

Fisher Catalog #17504044

[GlutaMAX™ Supplement Gibco - Thermo](#)

Fisher Catalog #35050061

A	B
Neurobasal medium	
B27 supplement, serum free	1X
GlutaMAX supplement	1%

[Poly-L-lysine](#)

6. [hydrobromide Merck Catalog #P2636](#)

[Water sterile-filtered BioReagent suitable for cell](#)

7. [culture Merck Catalog #W3500](#)

[Trypan Blue solution Sigma -](#)

8. **Aldrich Catalog #T8154**

For biochemistry

1. Mitochondrial depolarisation:

[Antimycin A from Streptomyces sp. Sigma -](#)

10 µM Aldrich Catalog #A8674

[Oligomycin A Sigma -](#)

1 µM Aldrich Catalog #75351

in DMSO

[☒ Dimethyl sulfoxide Sigma](#) –

Aldrich Catalog #D2650

2. Lysis Buffer:

[☒ cOmplete™ EDTA-free Protease Inhibitor](#)

Cocktail Roche Catalog #11873580001

[☒ 2-Chloroacetamide Sigma](#) –

Aldrich Catalog #C0267

A	B
Tris-HCl (pH 7.5)	50 mM
Sucrose	250 mM
EDTA	1 mM
EGTA	1 mM
Sodium orthovanadate	1 mM
Sodium β-glycerophosphate	10 mM
Sodium fluoride	50 mM
Sodium pyrophosphate	5 mM
TritonX	1% (w/v)
Protease inhibitor cocktail supplemented with 200 mM chloroacetamide	

[☒ DPBS no calcium no magnesium Gibco - Thermo](#)

3. Fischer Catalog #14190094

[☒ Coomassie Protein Assay Reagent Thermo](#)

4. Scientific Catalog #1856209

[☒ NuPAGE™ LDS Sample Buffer \(4X\) Invitrogen - Thermo](#)

5. Fisher Catalog #NP0008

[☒ 2-Mercaptoethanol Sigma](#) –

6. Aldrich Catalog #M6250

[☒ PageRuler™ Prestained Protein Ladder 10 to 180 kDa Thermo Fisher](#)

7. Scientific Catalog #26616

[☒ Immobilon-P PVDF](#)

8. Membrane Merck Catalog #IPVH00010

[☒ Amersham™ Protran® Western blotting membranes](#)

9. nitrocellulose Merck Catalog #GE10600041

10.

[☒ NuPAGE™ 4 to 12% Bis-Tris 1.0 mm Mini Protein Gel 10-well Invitrogen - Thermo](#)

Fisher Catalog #NP0321BOX

[NuPAGE™ 4 to 12% Bis-Tris 1.0 mm Midi Protein Gel 20-well](#) **Invitrogen - Thermo**

Fisher Catalog #WG1402BOX

[NuPAGE™ MOPS SDS Running Buffer \(20X\)](#) **Invitrogen - Thermo**

11. **Fisher Catalog #NP000102**

12. 1 X Towbin transfer buffer: [25 mM](#) Tris, [192 mM](#) Glycine, 20% methanol

13. 1X Tris Buffered-Saline (TBS): [500 mM](#) Tris, [150 mM](#) sodium chloride, pH 7.6, at [25 °C](#) .

14. 1X Tris-Buffered Saline, 0.1% Tween® 20 Detergent (TBST)

15. 5 % non-fat milk in TBST

[Bovine Serum Albumin Fraction V](#) **Sigma –**

16. 5% **Aldrich Catalog #10735094001**

17. **Primary antibodies:**

[Phospho-Ubiquitin \(Ser65\) \(E2J6T\) Rabbit mAb](#) **Cell Signaling**

Technology Catalog #62802

[Purified anti-Ubiquitin](#)

Antibody BioLegend Catalog #646302

Anti-Parkin phospho-Ser65 Rabbit mAb by Epitomics in collaboration with the Michael J Fox Foundation for Research,

[Parkin Antibody \(PRK8\)](#) **Santa Cruz**

Biotechnology Catalog #32282

[GAPDH Antibody \(6C5\)](#) **Santa Cruz**

Biotechnology Catalog #32233

18. **Secondary Antibodies:**

[Goat anti-Rabbit IgG \(H L\) Secondary Antibody HRP](#) **Invitrogen - Thermo**

Fisher Catalog #31460

[Rabbit anti-Mouse IgG \(H L\) Secondary Antibody HRP](#) **Invitrogen - Thermo**

Fisher Catalog #31450

19. ECL™ Western Blotting Reagents (Merck, Cytiva, #RPN2106)

[ECL™ Western Blotting](#)

Reagents Merck Catalog #RPN2106

20. Amersham Hyperfilm ECL (Merck, Cytiva, #28906837)

[Hyperfilm™](#)

ECL™ Merck Catalog #28906837

STOCK SOLUTION PREPARATION:

- **Poly-L-lysine solution:** Dissolve **10 mg/ml** of poly-L-lysine in sterile water; filter, aliquot and store at **-20 °C**. The solution is stable for 2–3 months.
- **DNaseI:** Dissolve **100 mg/ml** (wt/vol) DNase in sterile double-distilled water; filter, aliquot and store at **-20 °C**. The solution is stable for 2–3 months.
- **Antimycin A:** Prepare **50 mM** of Antimycin A in DMSO; aliquot and store at **-20 °C**.
- **Oligomycin A:** Prepare **10 mM** of Antimycin A in DMSO; aliquot and store at **-20 °C**.

EQUIPMENT

- Dumont #5 Forceps Biologie Inox (Fine Science Tool #11252-20)

Dumont #5 Forceps
Forceps
Fine Science Tools 11252-20 [↗](#)

- Dumont #5XL Forceps Standard Inox (Fine Science Tool #11253-10)

Dumont #5XL Forceps
Forceps
Fine Science Tools 11253-10 [↗](#)

- Dumont #7 Fine Forceps Biologie Inox (Fine Science Tool #11274-20)

Dumont #7 - Fine Forceps
Forceps
Fine Science Tools 11274-20 [↗](#)

- Dumont #5 45 Forceps Standard Dumoxel (Fine Science Tool #11251-35)

Dumont #5/45 Forceps
Forceps
Fine Science Tools 11251-35 [↗](#)

- Fine Scissors ToughCut Straight 9cm (Fine Science Tool #14058-09)

Fine Scissors - ToughCut®

Scissor

Fine Science Tools 14058-09 [↗](#)

- Dissecting microscope Nikon (P-PS32 Plain Stand) with fiber light illumination unit
- Cell Counter-DeNovix CellDrop™
- 37 °C water bath
- Laminar flow cell culture hood
- Cell culture incubator 5% CO₂, 95% humidity HERAcell®CO₂ incubator (150 L)
- Microcentrifuges, Micro Star 17R (VWR #521-1647)

Microcentrifuges, ventilated/refrigerated,

Micro Star 17 / 17R

Microcentrifuges

VWR 521-1647 [↗](#)

- Probe sonicator, Branson Digital Sonifier.
- XCell SureLock™ 4Midi-Cell running tank (Invitrogen™ #WR0100)

XCell4 SureLock™ Midi-Cell

Midi-Cell running tank

Invitrogen WR0100 [↗](#)

- XCell SureLock™ Mini-Cell running tank (Invitrogen™ #EI0001)

XCell SureLock Mini-Cell Electrophoresis

System

Electrophoresis System

Invitrogen EI0001

- Mini Trans-Blot®Cell transfer tank (BIORAD # 1703930)

Mini Trans-Blot Electrophoretic Transfer
Cell
Electrophoresis System
Mini Trans-Blot 1703930 [↗](#)

- Trans-Blot®Cell transfer tank (BIORAD # 1703939)


Trans-Blot Cell With Plate Electrodes and
Super Cooling Coil
Trans-Blot 1703939 [↗](#)

- ChemiDoc MP Imaging System (BIORAD)
- ECOMAX™ X-ray Processor

CONSUMABLES

- Cell culture multidishes, 6 well (Thermo Scientific™ #140675)

Nunc 6-Well Plate, Round
cell culture plate
Thermo Fisher Scientific 140675 [↗](#)
Nunc 6-Well Plate, Round



- Cell strainer 40µM (Greiner Bio-one #5420400)

EASYSTRAINER 40 µM, FOR 50 ML TUBES,
FOR TUBES 227XXX/210XXX, GREEN,
STERILE, SINGLE PACKED
Cell strainer
Greiner Bio-one 542040 [↗](#)

- Stericups 0.22µm, 250 mL and 500 mL (Merck #SCGPU02RE, #SCGPU05RE)

Stericup-GP Sterile Vacuum Filtration
System
Sterile Vacuum Filtration System
Merck SCGPU05RE [↗](#)

Stericup-GP Sterile Vacuum Filtration
System
Sterile Vacuum Filtration System
Merck SCGPU02RE [↗](#)

- 50mL Stripette® Serological Pipets (Corning #4490)

Stripette™ Serological Pipets
Serological Pipet
Stripette 4490 [↗](#)

- 25mL Stripette® Serological Pipets (Corning #4489)

25 mL Stripette™ Serological Pipets,
Polystyrene, Individually Paper/Plastic
Wrapped, Sterile, 25/Bag, 200/Case
Serological Pipet
Stripette 4489 [↗](#)

- 10mL Stripette® Serological Pipets (Corning #4488)

Stripette™ Serological Pipets
Serological Pipet
Stripette 4488 [↗](#)

- 5mL Stripette® Serological Pipets (Corning #4487)

5 mL Stripette™ Serological Pipets,
Polystyrene, Individually Paper/Plastic
Wrapped, Sterile, 50/Bag, 200/Case
Serological Pipet

Stripette 4487 [↗](#)

- 15ml CELLSTAR® tubes (Greiner bio-one. Catalog# 188271)

15 ml CELLSTAR® Polypropylene Tube
Polypropylene Tube

CELLSTAR 188271 [↗](#)

- 50ml CELLSTAR® tubes (Greiner bio-one. Catalog# 227261)

50 ml CELLSTAR® Polypropylene Tube
Polypropylene Tube

CELLSTAR 227261 [↗](#)

- Standard 1ml and 200µl Pipette tips (Greiner bio-one # 686271, #685261)

100 - 200 µl Pipette Tips
Pipette Tip

Greiner bio-one 685261 [↗](#)

1000 µl Pipette Tips
Pipette Tip

Greiner bio-one 686271 [↗](#)

- Syringe filter (0.22µm. Sartorius, Item # ST16541-Q)

Minisart® Syringe Filter, Polyethersulfone (PES), Pore Size 0.22 µm, Non-Sterile, Female Luer Lock, Male Luer Slip, Pack Size 500
Syringe Filter

Sartorius ST16541-Q [↗](#)

- Syringes (50ml) (Terumo™ # 8SS50L1)

Terumo™ 3-Part 50mL Luer Lock Syringes
Luer Lock Syringe

Terumo 8SS50L1 [↗](#)

- 1.5ml Eppendorf tubes (Eppendorf™ # 0030120086)

Eppendorf® Safe-Lock micro test tubes
Micro test tube

Eppendorf 0030120086 [↗](#)

- Disposable Cell Lifter (FisherBrand #08100240)

Fisherbrand™ Cell Lifters
Cell Lifter

Fisherbrand 08100240 [↗](#)

Coating of multiwell dishes with poly-L-lysine



1

⏏TIMING ⌚ 00:30:00 , 1 d before culture in this section.

In a sterile laminar hood, prepare poly-L-lysine [M]1 mg/ml in sterile from [M]10 mg/ml stock solution.

2 

Add enough poly-L-lysine solution to cover the bottom of the well, ensure that the volume of poly-L-lysine covers the bottom of the dish entirely, usually 1 mL for 1 single well.

3 

Wrap the dish in aluminum foil to prevent evaporation and leave it ⌚ Overnight at ⚡ Room temperature .

It is important to make sure that the poly-L-lysine does not dry out during incubation. If you are working in the laminar hood, be sure to turn the blower off during the incubation

Washing of dishes after coating

4

⏏TIMING ⌚ 01:00:00 , day of culture in this section.

Aspirate the poly-L-lysine carefully.

5 

Add 1 mL of sterile water into each well and aspirate; repeat this step three times.

6 Aspirate water and leave dishes under the hood until completely dry.

Dissection of cortex from E16.5 mouse embryos





7 

⚡TIMING ⌚01:00:00 , ⌚00:02:00 – ⌚00:05:00 for each embryo in this section.

Use sterilized instruments by autoclave or washing them with 70% (vol/vol) ethanol.

Dry thoroughly if ethanol is used.

8 

Prepare 60-mm dishes with dissection medium. If you are culturing from individual embryos, prepare  15 mL tubes with  1 mL of dissection medium. If you are culturing pooled embryos, prepare  15 mL tubes with  5 mL of dissection medium.

9 Euthanize the embryos pup by decapitation and separate the head from the body.

10 Place the head on a dish with dissection medium and hold down the sides with forceps.

11 Under a dissecting microscope, dissect the skin on the top of the head and hold down the skin on either side with the forceps.

12 Cut open the skull by making an incision at the base of the brain. Separate the two halves of the skull and remove carefully.

Take care to not cut through the brain tissue when removing the skull bone

13 

By using forceps, pinch off the brain from the base and peel off the meninges carefully, ideally as a single piece. Check for the remaining pieces of meninges and remove them completely.

The meninges can be sticky and difficult to remove, it is important to ensure that the meninges are completely removed, so that they do not contribute any non-neuronal cells to the culture.





14 Separate the two halves of the brain by making a sagittal cut along the midline.

15 

Orient the tissue so that the hippocampus is on the top. The hippocampus can be identified by its C-shaped structure and opacity, which differ from the neighboring cortical tissue. Cut away hippocampus and flat the cortex to expose the striatum. Using an angled forceps scoop out the striatum and cortex can also be dissected out and processed for cortical cultures.

It is important to ensure that the hippocampus and striatum are completely removed, so that they do not affect cortical distribution of neurons and contribute to an increase of glia cells.

16  

To culture neurons from individual embryo, put each brain in a separate  15 mL tube containing  1 mL of dissection medium. For pooled cultures, put cortices from two embryos into one  15 mL tube containing  5 mL of dissection medium.

Ensure that the cortices are submerged in the medium and do not let them dry out at any point. It is important that the dissection be done as quickly as possible to ensure cell viability and health.




17 For individual cultures, collect a piece of tail for genotyping.

Cell dissociation and plating



42m

18  

⬠TIMING ~  01:00:00 -  01:30:00 in this section.

Prepare digestion medium by adding  125 µl of DNase I (stock solution  10 mg/ml) to  10 mL of Trypsin 0,05%.

19  

Add 1:1 digestion medium to dissection medium containing the cortices, for individual mouse culture add  1 mL of digestion media, for pooled mouse culture add  5 mL of digestion media.

Note. It is common to use HBSS Ca²⁺ and Mg²⁺ free buffer. In our hands, we have noticed that using HBSS plus Ca²⁺ and Mg²⁺ ensure a milder but optimal concentration of trypsin that allows to reduce excessive number of dead cells and higher yield of cells. At the same time, it stimulates DNase I activity.

20



30m

Incubate at **37 °C** in a water bath for **00:30:00**.

It is important to ensure that this incubation does not proceed for longer than **00:30:00**

21



Inactivate trypsin digestion by adding dissociation medium, **2 mL** for individual mouse culture or **5 mL** for pooled mouse culture.

22




5m

Centrifuge at **1200 rpm** for **00:05:00**.

23

Resuspend cortices in **2 mL** (for individual embryos) or **4 mL** (for pooled embryos) of dissociation medium.

24



Carefully triturate the tissue 20 times to dissociate the cells gently and obtain a homogenous cell suspension, by using a p1000 pipette.

The trituration of cells should be done slowly and carefully to minimize damage to cells. It is best to avoid any bubbling during the procedure.

25

Remove any chunks of tissue by using a cell strainer and distribute drop by drop the cell suspension.

26



7m

Centrifuge at **700 rpm** for **00:07:00** to pellet down the cells and resuspend them in the culturing medium.

Note: It is important to remove FBS from the medium to reduce the proliferation of glial cells.

27 Resuspend the dissociated cells in **2 mL** (for individual embryos) or **4 mL** (for pooled embryos) of culturing medium.

28

Take a **15 µl** aliquot, add 1:1 ratio Trypan Blu and determine the density of cells and cell viability to the cell counter. Typical yields are $\sim 2 \times 10^6$ /mL and viability >80 %.

29 5.0×10^5 cells/well plates are plated out on 6-well multidishes, containing **2 mL** of pre-warmed culturing media.

Maintenance of neurons

30

◇TIMING **72:00:00** Up to **504:00:00** in this section.

Every 5 days aspirate 1/3 of the media from each well and replace it with fresh culturing medium warmed to **37 °C**. These neurons can be maintained in culture for up to DIV 28 and be used anytime during this period (DIV0-28). The culture can be used for examining PINK1 activation at 21DIV, since at this stage they express functional activity of Parkin and PINK1.

It is important to not replace the entire medium, neurons secrete factors that promote growth and survival.