



Sep 09, 2021

Longmire lysis buffer

Abigail Wells¹, Linda Park¹¹NWFSC-NOAA

1 Works for me



Share

dx.doi.org/10.17504/protocols.io.bx4xpqxn

Ana Ramón-Laca

ABSTRACT

Lysis buffer recipe (Longmire et al 1997):
To make 1 liter

DOI

dx.doi.org/10.17504/protocols.io.bx4xpqxn

PROTOCOL CITATION

Abigail Wells, Linda Park 2021. Longmire lysis buffer. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bx4xpqxn>

MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

Use of "lysis buffer" in DNA isolation and its implications for museum collections J. Longmire, M. Maltbie and R. J. Baker Occasional Papers Museum of Texas Tech University 1997 Vol. 163 DOI: 10.5962/bhl.title.143318

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 09, 2021

LAST MODIFIED

Sep 09, 2021

PROTOCOL INTEGER ID

53111

PARENT PROTOCOLS

In steps of

[Collection and preservation of eDNA from marine water samples](#)[eDNA extraction: phenol-chloroform-isoamyl alcohol DNA purification from filters stored in Longmire buffer](#)

MATERIALS TEXT

H₂O

[M]**1 Molarity (M)** Tris-HCL, **pH8**

[M]**0.5 Molarity (M)** EDTA, **pH8**

[M]**5 Molarity (M)** NaCl

[M]**20 Mass Percent** SDS (w/v)

To make 1 L of Longmire

- 1 975 ml double-distilled water
- 2 100 ml of 1 M Tris-HCL, pH 8.0
- 3 200 ml of 0.5 M EDTA, pH8.0
- 4 2 ml of 5 M NaCl
- 5 25 ml of 20% SDS (w/v) Filter the buffer with an autofill PES bottle top filtration device (sterile 500ml, 0.22 µm)