



MAR 06, 2023

# 🌐 Preparation and transformation of electrocompetent cells

Andreas Sagen<sup>1</sup>

<sup>1</sup>University of Oslo



Andreas Sagen

University of Oslo, The National Institute of Occupational H...

## ABSTRACT

A protocol for electroporation of E. coli. Other bacteria may work, with optimization of transformation buffer and settings.

## MATERIALS

Centrifuge  
LAF  
Vortex  
Incubator  
Electroporator

## OPEN ACCESS

**Protocol Citation:** Andreas Sagen 2023. Preparation and transformation of electrocompetent cells.

**protocols.io**

<https://protocols.io/view/preparation-and-transformation-of-electrocompetent-cp26vqhe>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** In development  
We are still developing and optimizing this protocol

**Created:** Feb 26, 2023




**Last Modified:** Mar 06, 2023

**PROTOCOL integer ID:**  
77630

**Keywords:** bacteria, electrocompetent, GYT medium

## Preparation of GYT medium

1 In a sterile flask, add  400 mL distilled water


- 2 Measure and add  0.6 g yeast extract,  1.25 g tryptone, and  50 mL glycerol

Materials:


 Select Yeast Extract **Merck MilliporeSigma (Sigma-Aldrich) Catalog #Y0875**

 Tryptone **Merck Millipore (EMD Millipore) Catalog #T9410**

 Glycerol **MP Biomedicals Catalog #194680**


- 3 Filter sterilize through a 0.22- $\mu$ m filter, and store in aliquots at  4 °C



## Preparation of cells

- 4 Inoculate 500 mL of prewarmed LB medium from 25 mL overnight E. coli culture. Incubate at  37 °C and 300 rpm. Measure OD every 20 minutes, until OD<sub>600</sub>=0.4

### Note



The density is usually archived after ~2.5 hours of incubation with DH5 $\alpha$



- 5 Transfer the culture to appropriate centrifugation containers, and cool on ice for  00:30:00 30m



- 6 Centrifuge culture at  1000 rcf, 4°C, 00:15:00, and resuspend pellet in  200 mL ice-cold water 15m

### Note

Combine the culture into a smaller number of centrifugation containers

- 7 Centrifuge culture at  1000 rcf, 4°C, 00:20:00, and resuspend pellet in  100 mL ice-cold 10% glycerol 20m

8 Centrifuge culture at  1000 rcf, 4°C, 00:20:00 , and resuspend pellet in  25 mL ice-cold 10% glycerol 20m

9 Centrifuge culture at  1000 rcf, 4°C, 00:20:00 , and resuspend pellet in  2 mL ice-cold GYT medium 20m

#### Note



This is best done by gentle swirling rather than pipetting or vortexing

10 Measure OD<sub>600</sub> of a 1:100 dilution of the cell suspension, and dilute the cells to between  $2-3 \times 10^{10}$  cells/mL

#### Note




1.0 OD<sub>600</sub> =  $\sim 2.5 \times 10^8$  cells/ml

11 Transfer 40 µl of the suspension to an ice-chilled electroporation cuvette and test whether arching occurs when an electrical discharge is applied. If so, wash the remainder of the cell suspension once more with ice-cold GYT medium to ensure that the conductivity of the bacterial suspension is sufficiently low

12 Dispense  40 µL aliquots of the cell suspension into sterile, ice-cold microfuge , drop into a bath of liquid nitrogen, and transfer to a  -80 °C freezer

## Transformation

1h 5m

13 Pre-chill cuvettes on ice for  00:05:00 , and pre-heat LB plates with an appropriate selection agent and SOC medium at  37 °C for  01:00:00 1h 5m

14 Mix an appropriate amount plasmid to an aliquot of electrocompetent cells and transfer to a pre-chilled cuvette

### Note

100 pg pUC19 plasmid is appropriate in most cases, but a specific amount plasmid depend on many different factors and have to be optimized

### Protocol



NAME




**S. O. C. medium**

CREATED BY

**Andreas Sagen**

**PREVIEW**

**15** Incubate plasmid-suspension mix for 1 minute, then perform electroporation with optimized settings

**16** Flush  1 mL pre-heated S. O. C. medium, then transfer to a culture tube and recover at  37 °C and 200 rpm shaking for  01:00:00

1 h

**17** Plate an appropriate amount culture on selection agar plates